Amies Transport Medium w/o Charcoal

Semi-solid medium for transport of clinical specimens and other microorganisms

INTENDED PURPOSE

Semi-solid medium for transport of clinical specimens and other microorganisms.

DESCRIPTION

Amies Transport Medium w/o Charcoal is used for collecting, transporting and preserving microbiological specimens. It is formulated to maintain the viability of microorganisms without significant increase in growth, being nonnutritive, phosphate buffered and semi-solid.

TYPICAL FORMULA

Composition	(g/l)
Sodium Chloride	3.0
Calcium Chloride	0.1
Potassium Chloride	0.2
Magnesium Chloride	0.1
Monopotassium Phosphate	0.2
Disodium Phosphate	1.15
Sodium Thioglycollate	1.0
Agar	4.0
Final pH = 7.2±0.2 at 25 °C	

METHOD PRINCIPLE

Amies Transport Medium w/o Charcoal provides a reduced environment due to the presence of sodium thioglycollate and the small amount of agar.

The Chloride salts supply essential electrolytes for transport and osmotic balance. Phosphates act as a buffer system. Sodium thioglycollate lowers the redox potential of the medium, thereby allowing a better conservation of anaerobic bacteria. The survival of bacteria in a transport medium depends on various factors such as bacteria type and concentration in the specimen, transport medium formulation, the temperature and duration of transport, and inoculation to appropriate culture media within 24 hours. Optimal growth and typical morphology can only be expected if direct inoculation and appropriate cultivation are followed.

PREPARATION

Suspend 9.8 g of powder in 1 liter of distilled or deionized water.

Heat to boiling, shake until completely dissolved and dispense into final containers as appropriate.

Sterilize in autoclave at 121°C for 15 minutes.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, culture media, swab, quality control organisms.

SPECIMENS

Collect specimens with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

After collecting the clinical specimen with a swab, insert one third of the swab into the center of the medium. Then cut the stick and screw down the cap to force the swab to the bottom of the medium.

Make sure the cap is screwed firmly on the bottle and keep cool during the transport period.

STORAGE

The powder is very hygroscopic: store the powder at 10-30°C, in a dry environment, tightly closed in its original container and use it before the expiry date on the label or until signs of deterioration or contamination are evident. Store prepared tubes at 2-8 °C.

SHELF LIFE

4 years

QUALITY CONTROL

Appearance of dehydrated medium: beige, free-flowing, homogeneous.

Appearance of prepared medium: clear, colorless to light amber.

Control strain		Inoculum	Incubation	Characteristic reactions
<i>Haemophilus influenzae</i>	ATCC® 49766	50-100 CFU	Microaerobic 18-24 h, 35 ± 2°C	Good growth
<i>Neisseria gonorrhoeae</i>	ATCC® 19424			Good growth
<i>Streptococcus pneumoniae</i>	ATCC® 49619			Good growth
<i>Streptococcus pyogenes</i>	ATCC® 19615			Good growth

PERFORMANCE CHARACTERISTICS

Performance testing of Amies Transport Media w/o Charcoal was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

- Some growth of contaminants may also occur during the long period of transport.
- After transportation, the specimen should be inoculated in the proper medium as soon as possible. For optimum results, the time lapse between sample collection and inoculum onto culture medium should be reduced to the minimum.
- It may not be suitable for the transport of fastidious organisms.
- Gonococci survive well in Amies Transport Medium for at least 6 to 12 hours provided they are not exposed to temperature extremes. By 24 hours, the numbers of gonococci decrease to an extent that may prevent their recovery if small numbers were present initially in the specimen.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Amies Transport Media w/ o Charcoal	500 g	610191
Amies Transport Media w/o Charcoal	100 g	620191

Revision History

Revision	Release Date	Change Summary
1	2024-09-16	Added notice to report any malfunction, defect or incident.
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.

*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

[liofilchem.com/ifu-sds](https://www.liofilchem.com/ifu-sds)

Amies Transport Medium w/o Charcoal

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