## Chocolate Agar

Highly nutrient medium for isolation and cultivation of fastidious microorganisms from clinical specimens.

### DESCRIPTION

Chocolate Agar is a non-selective medium used for the isolation and cultivation of fastidious organisms, especially *Neisseria* and *Haemophilus* spp, from clinical specimens and other materials of sanitary importance.

| TYPICAL FORMULA                    | (g/l)   |
|------------------------------------|---------|
| Enzymatic Digest of Casein         | 7.5     |
| Enzymatic Digest of Animal Tissues | 7.5     |
| Corn Starch                        | 1.0     |
| Dipotassium Phosphate              | 4.0     |
| Monopotassium Phosphate            | 1.0     |
| Sodium Chloride                    | 5.0     |
| Agar                               | 17.0    |
| Defibrinated Horse Blood           | 50.0 ml |
| Vitalex Growth Supplement          | 10.0 ml |
| Final pH 7.2 ± 0.2 at 25°C         |         |

### METHOD PRINCIPLE

Enzymatic digest of casein and enzymatic digest of animal tissues provide amino acids, carbon, nitrogen, vitamins and minerals for organisms growth. Corn starch absorbs any toxic metabolites produced. The phosphates are buffering agents. Sodium chloride maintains the osmotic balance of the medium. Agar is the solidifying agent. Blood, that is a good source of hemin and Vitalex Growth Supplement providing cofactors, vitamins and NAD, are incorporated into the medium to enhance recovery of *Haemophilus* and *Neisseria* spp.

### PREPARATION

|                          |                                                                                                                                                                                                                                                                                                                                                 |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <u>Medium in bottles</u> | Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes. |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### TEST PROCEDURE

Inoculate the medium by streaking the specimen\* over the agar surface. Slant tubes are primarily used for the cultivation and maintenance of pure cultures. Incubate at 35 ± 2°C for 18-48 h in an aerobic atmosphere supplemented with carbon dioxide.

\*Specimens must be sent to the laboratory by using appropriate transport media at 20-25°C (DO NOT REFRIGERATE) and inoculated as soon as possible after collection.

### INTERPRETING RESULTS

Examine for bacterial growth. To identify the isolated microorganisms perform adequate biochemical or immunological tests.

### APPEARANCE

Opaque, chocolate brown.

### STORAGE AND TRANSPORT CONDITIONS

Store bottles, tubes and prepared plates at 2-8°C away from light. However, our stability studies have shown that the transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, does not alter in any way the performance of the product. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

### SHELF LIFE

Medium in tubes/bottles: 1 year.  
Ready-to-use plates: 5 months.

**QUALITY CONTROL**

Plates are inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: 50-100 CFU.

Incubation conditions: 35 ± 2°C for 18-48 hours, in CO<sub>2</sub> enriched atmosphere.

**QC Table.**

| Microorganism                   |             | Growth |
|---------------------------------|-------------|--------|
| <i>Haemophilus influenzae</i>   | ATCC® 10211 | Good   |
| <i>Neisseria gonorrhoeae</i>    | ATCC® 43069 | Good   |
| <i>Neisseria meningitidis</i>   | ATCC® 13090 | Good   |
| <i>Streptococcus pneumoniae</i> | ATCC® 6305  | Good   |
| <i>Streptococcus agalactiae</i> | ATCC® 13813 | Good   |

**WARNING AND PRECAUTIONS**

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used only by properly trained operators.

**DISPOSAL OF WASTE**

Disposal of waste must be carried out according to national and local regulations in force.

**BIBLIOGRAPHY**

1. Clinical Microbiology Procedures Handbook, 3<sup>rd</sup> ed. (2010) ASM. Washington, D.C. 2,540 pages.
2. Murray, P. R., E. J. Baron, J. H. Jorgensen, M. L. Landry, and M. A. Pfaller (ed.) Manual of Clinical Microbiology, 9<sup>th</sup> ed. (2007) ASM Press, Washington, D. C.
3. Chapin, C.K., G.V. Doern (1983). J. Clin. Microbiol. 17: 1163-1165.
4. Thayer J.D., R.L. Jackson (1975) Public Health Rep. 82:361.
5. Martin, J.E., Armstrong J.H., Smith P.B. (1974). Appl. Microbiol. 27: 802-805.
6. Seth, A. (1970). Brit. J. Vener. Dis. 46: 201-202.

**PRESENTATION**

|                |                           | Contents           | Ref.   |
|----------------|---------------------------|--------------------|--------|
| Chocolate Agar | 90 mm ready-to-use plates | 20 plates          | 10023  |
| Chocolate Agar | 90 mm ready-to-use plates | 100 plates         | 10023* |
| Chocolate Agar | Slant tubes               | 10 x 9 ml tubes    | 30099  |
| Chocolate Agar | Slant tubes               | 20 x 9 ml tubes    | 31099  |
| Chocolate Agar | Bottles                   | 6 x 500 ml bottles | 470120 |
| Chocolate Agar | Bottles                   | 6 x 100 ml bottles | 442220 |

**TABLE OF SYMBOLS**

|                             |                                                                                                            |                                                                                                                       |                                                                                                                           |                                                                                                                 |
|-----------------------------|------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>LOT</b> Batch code       | <b>IVD</b> <i>In vitro</i> Diagnostic Medical Device                                                       |  Manufacturer                      |  Use by                               |  Fragile, handle with care |
| <b>REF</b> Catalogue number |  Temperature limitation |  Contains sufficient for <n> tests |  Caution, consult Instruction For Use |  Do not reuse              |

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