INTEGRAL SYSTEM STAFILOCOCCI REF. 71718

- ***Sistema colorimetrico per l'identificazione biochimica e l'antibiogramma degli stafilococchi***
- ***Colorimetric system for biochemical identification and susceptibility testing of staphylococci***
- ***Système colorimétrique pour l'identification biochimique et l'antibiogramme de staphylocoques***
- ***Sistema colorimétrico para la identificación bioquímica y el antibiograma de los estafilococos***
- ***Sistema colorimétrico para a identificação bioquímica e o antibiograma das estafilococos***
- ***Σύστημα για τον βιοχημικό προσδιορισμό και το αντιβιογράμμα των σταφυλοκόκκων***

Code F00713
Rev. 1 / 15.03.2004

INTEGRAL SYSTEM STAFILOCOCCI

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INTEGRAL SYSTEM STAFILOCOCCI

Colorimetric system for biochemical identification and susceptibility testing of staphylococci

DESCRIPTION

INTEGRAL SYSTEM STAFILOCOCCI is a 24-well system containing desiccated biochemical and antibiotic substrates for the identification and susceptibility testing of staphylococci. The system is inoculated with the bacterial suspension of the micro-organism under examination and incubated at $36\text{ }^{\circ}\text{C} \pm 1^{\circ}\text{C}$ for 18-24 hours. The tests for identification and for the susceptibility testing are interpreted by assessing the change in colour of the various wells.

KIT CONTENTS

The kit contains:

- 20 Systems INTEGRAL SYSTEM STAFILOCOCCI
- 20 Vials of Inoculum Broth (7.0 ml/vial)
- 20 Xylose Discs
- 1 Instruction sheet

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

- Vaseline oil for microbiology use (Vaselin oil 2 bottles of 50 ml)..... cod. 80278
- Physiological Solution for microbiology use (Physiological Solution 20 vials of 7.0 ml/vial)cod. 20095
- Sundry microbiology laboratory materials

CONFIGURATION

The configuration of the system is shown in Table n°1.

Table n°1

Well	BIOCHEMICAL IDENTIFICATION
1-ADC	Decarboxylation of arginine
2-UR	Hydrolysis of urea
3-MAL	Fermentation of maltose
4-TRE	Fermentation of trehalose
5-MAN	Fermentation of mannitol
6-XYL	Fermentation of xylose
7-SAC	Fermentation of sucrose
8-MNN	Fermentation of mannose
Well	SUSCEPTIBILITY TESTING (*)
9-CN	Gentamicin - 8 µg/ml
10-TOB	Tobramycin - 8 µg/ml
11-AUG	Amoxicillin + clavulanic acid – 8/4 µg/ml
12-AMS	Ampicillin + sulbactam – 32/16 µg/ml
13-CAZ	Ceftazidime - 32 µg/ml
14-CXM	Cefuroxime - 32 µg/ml
15-CIP	Ciprofloxacin - 4 µg/ml
16-LEV	Levofloxacin - 8 µg/ml
17-AZM	Azithromycin - 8 µg/ml
18-CLR	Clarithromycin - 8 µg/ml
19-ERY	Erythromycin - 8 µg/ml
20-ROX	Roxithromycin - 8 µg/ml
21-FOS	Fosfomicin - 200 µg/ml
22-RD	Rifampicin - 4 µg/ml
23-SXT	Co-Trimoxazole - 8 µg/ml
24-C	Growth control for the susceptibility testing

(*): The concentration of each antibiotic is according to the updated NCCLS-January 2004, Vol 24 N°1 ⁽¹⁾

PRINCIPLE OF THE METHOD

INTEGRAL SYSTEM STAFILOCOCCI makes possible simultaneous biochemical identification and assessment of the susceptibility testing of micro-organisms of the genus *Staphylococcus* isolated from clinical samples. The identification is based on biochemical tests performed on culture media containing specific substrates in the wells from **1-ADC** to **8-MNN**.

The susceptibility testing is assessed on the basis of the growth or inhibition of the micro-organisms in the media containing the antibiotic and a growth indicator in the wells from **9-CN** to **23-SXT**. Well **24-C** does not contain antibiotics but only the culture medium and the indicator; it serves as a microbe growth control for assessment of the susceptibility testing.

COMPOSITION

Table n°2

Well	Content
1-ADC	Culture medium for showing arginine decarboxylase
2-UR	Culture medium for showing hydrolysis of urea
3-MAL	Culture medium for showing fermentation of maltose
4-TRE	Culture medium for showing fermentation of trehalose
5-MAN	Culture medium for showing fermentation of mannitol
6-XYL	Culture medium for showing fermentation of xylose
7-SAC	Culture medium for showing fermentation of sucrose
8-MNN	Culture medium for showing fermentation of mannose
9-CN	Culture medium containing Gentamycin - 8 µg/ml
10-TOB	Culture medium containing Tobramycin - 8 µg/ml
11-AUG	Culture medium containing Amoxicillin + Clavulanic acid - 8/4 µg/ml
12-AMS	Culture medium containing Ampicillin + Sulbactam - 32/16 µg/ml
13-CAZ	Culture medium containing Ceftazidime - 32 µg/ml
14-CXM	Culture medium containing Cefuroxime - 32 µg/ml
15-CIP	Culture medium containing Ciprofloxacin - 4 µg/ml
16-LEV	Culture medium containing Levofloxacin - 8 µg/ml
17-AZM	Culture medium containing Azitromycin - 8 µg/ml
18-CLR	Culture medium containing Claritromycin - 8 µg/ml
19-ERY	Culture medium containing Erythromycin - 8 µg/ml
20-ROX	Culture medium containing Roxitromycin - 8 µg/ml
21-FOS	Culture medium containing Fosfomycin - 200 µg/ml
22-RD	Culture medium containing Rifampicin - 4 µg/ml
23-SXT	Culture medium containing Co-Trimoxazole - 8 µg/ml
24-C	Growth control for the susceptibility testing

Inoculum solution (g/l): Anhydrous glucose **30 g**, Tryptone **5 g**, Sodium chloride **9 g**, Yeast extract **5 g**, Agar **1 g**,
Distilled water **1000.0 ml**, pH **6.4 ± 0.2**

COLLECTION AND STORAGE OF THE SAMPLES

The colonies to be subjected to the biochemical identification test and susceptibility testing assessment with INTEGRAL SYSTEM STAFILOCOCCI must be taken from a culture medium, selective or non selective, used for the isolation of staphylococci.

TEST PROCEDURE

1. Take a system from the kit.
2. After checking that the colonies developed on the culture media presumably belong to the genus *Staphylococcus*, take one or more morphologically similar, well isolated, colonies from a solid culture medium and suspend them in 5 ml physiological solution for microbiological use in such a way as to obtain a turbidity equivalent to 0.5 Mc Farland (**Bacterial suspension**).
3. Trasfer:
 - 0.2 ml of **Bacterial suspension** into each of the first 8 wells, and to wells **1-ADC** and **2-UR** add 2 drops of oil of vaseline for microbiological use; insert a disc of **Xylose disc** in well **6-XYL**. (BIOCHEMICAL IDENTIFICATION).
 - 0.2 ml of the **Bacterial suspension** into the phial of **Inoculum Broth** contained into the kit and distribute 0.2 ml into each of the wells from **9-CN** to **24-C** (SUSCEPTIBILITY TESTING).
4. Cover the system with the lid provided and incubate at $36^{\circ}\text{C} \pm 1^{\circ}\text{C}$ for 18-24 hours.

INTERPRETATION OF THE RESULTS

BIOCHEMICAL IDENTIFICATION

Interpret the results of the first 8 wells using Table n°3 and form the 3-digit code following the instructions given to the paragraph **FORMATION OF THE NUMERICAL CODE**. Then use the code and Table n°5 to identify the bacteria.

SUSCEPTIBILITY TESTING

Observe the change in colour of the wells from **9-CN** to **24-C** and interpret the results using Table n°3. The control well (**24-C**) must be positive (yellow). If it is negative (red or orange), it is necessary to check the vitality of the inoculum and whether it has been correctly prepared; then to repeat the tests using a new system.

Table 3

Well	BIOCHEMICAL IDENTIFICATION	Well colour	
		Positive reaction	Negative reaction
1-ADC	Decarboxylation of arginine	Violet	Yellow
2-UR	Hydrolysis of urea	Red- Fucsia	Yellow-Orange
3-MAL	Fermentation of maltose	Yellow	Red
4-TRE	Fermentation of trehalose	Yellow	Red
5-MAN	Fermentation of mannitol	Yellow	Red
6-XYL	Fermentation of xylose	Yellow	Red
7-SAC	Fermentation of sucrose	Yellow	Red
8-MNN	Fermentation of mannose	Yellow	Red

SUSCEPTIBILITY TESTING		
WELL COLOUR	BACTERIAL GROWTH	INTERPRETATION
Red	Inhibited	S = Sensitive
Orange	Intermediate	I = Intermediate
Yellow	Good	R = Resistant

FORMATION OF THE NUMERICAL CODE

- 1) The 8 biochemical tests have been divided into two groups containing 3 tests and one group containing 2, each one is indicated with a positiveness value of 1,2,4.
 - Value 1: first test positive in each group (ADC, UR, MAL);
 - Value 2: second test positive in each group (TRE,MAN,XYL);
 - Value 4: third test positive in each group (SAC, MNN);
 - Value zero: negative reaction in every group.
- 2) Adding the number of positive reactions in each group, one obtains a three digit code which is used to identify the micro-organism under examination from the table of numerical codes, as in the example.

	Group I			Group II			Group III	
Well	ADC	UR	MAL	TRE	MAN	XYL	SAC	MNN
Positiveness code	1	2	4	1	2	4	1	2
Results	+	+	+	+	+	-	+	+
Codes addition	7			3			3	
NUMERICAL CODE: 733				IDENTIFICATION: Staphylococcus aureus				

TABLE OF REACTIONS

Table n°4

Micro-organism	Group I			Group II			Group III	
	ADC	UR	MAL	TRE	MAN	XYL	SAC	MNN
<i>Staphylococcus aureus</i>	V	+	+	+	+	-	+	+
<i>Staphylococcus epidermidis</i>	V	+	+	-	-	-	+	V
<i>Staphylococcus saprophyticus</i>	-	V	+	+	V	-	+	-
<i>Staphylococcus cohnii</i>	-	-	V	+	V	-	-	V
<i>Staphylococcus xylosus</i>	-	+	+	+	V	+	+	+
<i>Staphylococcus haemolyticus</i>	+	-	+	+	V	-	+	-
<i>Staphylococcus hominis</i>	V	+	+	V	-	-	V	-
<i>Staphylococcus warneri</i>	V	+	V	+	V	-	+	-
<i>Staphylococcus capitis</i>	V	-	-	-	+	-	V	+
<i>Staphylococcus simulans</i>	+	+	-	V	+	-	+	V
<i>Staphylococcus sciuri</i>	-	-	V	+	+	-	+	V

+ : Positive Reaction

- : Negative Reaction

V: Variable Reaction

TABLE OF NUMERICAL CODES

Table n°5

Code	Micro-organism	Code	Micro-organism	Code	Micro-organism
010	<i>S. cohnii</i>	333	<i>S. simulans</i>	611	<i>S. hominis</i>
012	<i>S. cohnii</i>	410	<i>S. cohnii</i>		<i>S. saprophyticus</i>
022	<i>S. capitis</i>	411	<i>S. saprophyticus</i>		<i>S. warneri</i>
023	<i>S. capitis</i>	412	<i>S. cohnii</i>	631	<i>S. warneri</i>
030	<i>S. cohnii</i>	430	<i>S. cohnii</i>		<i>S. saprophyticus</i>
031	<i>S. sciuri</i>	431	<i>S. saprophyticus</i>	633	<i>S. aureus</i>
032	<i>S. cohnii</i>		<i>S. sciuri</i>	653	<i>S. xylosus</i>
033	<i>S. sciuri</i>	432	<i>S. cohnii</i>	673	<i>S. xylosus</i>
122	<i>S. capitis</i>	433	<i>S. sciuri</i>	700	<i>S. hominis</i>
123	<i>S. capitis</i>	511	<i>S. haemolyticus</i>	701	<i>S. hominis</i>
211	<i>S. warneri</i>	531	<i>S. haemolyticus</i>		<i>S. epidermidis</i>
231	<i>S. warneri</i>	600	<i>S. hominis</i>	703	<i>S. epidermidis</i>
311	<i>S. warneri</i>	601	<i>S. epidermidis</i>	710	<i>S. hominis</i>
321	<i>S. simulans</i>		<i>S. hominis</i>	711	<i>S. hominis</i>
323	<i>S. simulans</i>	603	<i>S. epidermidis</i>		<i>S. warneri</i>
331	<i>S. warneri</i>	610	<i>S. hominis</i>	731	<i>S. warneri</i>
	<i>S. simulans</i>			733	<i>S. aureus</i>

QUALITY CONTROL

Each batch of INTEGRAL SYSTEM STAFILOCOCCI is subjected to quality control using the following ATCC reference strains .

Staphylococcus aureus..... ATCC 25923; *Staphylococcus xylosus*.....ATCC 35033
Staphylococcus epidermidisATCC 12228.

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardisation of the inoculum; application of the method to micro-organisms not in the genus *Staphylococcus*; mixed or contaminated cultures; use of expired systems or expired reagents; incorrect application of the technique.

LIMITS AND WARNINGS

In some cases it may be necessary to use supplementary confirmation tests for definitive identification of the micro-organisms.

PERFORMANCE

The bacterial identification results obtained with the INTEGRAL SYSTEM STAFILOCOCCI agree with those obtained using traditional test-tube methods⁽⁵⁾.

The antibiogram results obtained with the INTEGRAL SYSTEM STAFILOCOCCI agree with those obtained using the radial diffusion methods according to Bauer et al.⁽⁶⁾ recommended by the Food and Drug Administration (FDA) and by the National Committee for Clinical Laboratory Standards (NCCLS) U.S.A ^(1,5).

PRECAUTIONS

The product INTEGRAL SYSTEM STAFILOCOCCI cannot be classified as hazardous under current legislation nor does it contain harmful substances in concentrations $\geq 1\%$. It therefore does not require a Safety Data Sheet to be available. INTEGRAL SYSTEM STAFILOCOCCI is a disposable device to be used only for diagnostic use *in vitro*; it is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, INTEGRAL SYSTEM STAFILOCOCCI and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Code	Kit
INTEGRAL SYSTEM STAFILOCOCCI	71718	20 tests

TABLE OF SYMBOLS

<u>SYMBOL</u>	<u>MEANING</u>	<u>SYMBOL</u>	<u>MEANING</u>
	In Vitro Diagnostic Medical Device		Do not reuse
	Manufacturer		Contains sufficient for <n> tests
	Catalogue number		Fragile, handle with care
	Use by		Caution, consult accompanying documents
	Temperature limitation		Batch code

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1. NCCLS – Performance Standards for antimicrobial Susceptibility Testing, fourteenth Informational Supplement- January 2004, Vol 24 N°1.
2. Kloos W. E. and Jorgensen J. H. "Staphylococci". Manual of Clinical Microbiology. 4th edition, 15:143-153.
3. Baker C. et al. Inoculum Standardization in Antimicrobial Susceptibility Testing: Evaluation of Overnight Agar Cultures and the Rapid Inoculum Standardization System . 1983. *J. Clin. Microbiol.* **17**, 450-457.
4. Lampe M. F. et al. Relations of Early Readings of Minimal Inhibitory Concentrations to the Results of Overnight Tests. 1957. *Antimicrob. Agents Chemoth.* **8**, 429-433.
5. Dati di Archivio Liofilchem s.r.l. 'Integral System Stafilococchi' (Maggio 2003).
6. Bauer, A.W., W.M.M. Kirby, J.C. Sherris and M. Turk. 1966. Antimicrobic susceptibility testing by a standardized single disk method. *AM. J. Clin. Pathol.* 45:493-496

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