remel

RapID™ CB Plus System

REF R8311008..... 20

1. INTENDED USE

The RapID™ CB Plus is a qualitative micromethod using enzyme reactions to identify clinical isolates of *Corynebacterium* species and other irregular Gram-positive bacilli. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

A complete listing of the organisms addressed by the RapID CB Plus System is provided in the RapID CB Plus Differential Chart.

2. SUMMARY AND EXPLANATION

The RapID CB Plus System is comprised of (1) RapID CB Plus Panels and (2) RapID CB Plus Reagent. RapID CB Plus Panels are disposable plastic trays with 18 reaction cavities, which contain dehydrated reactants. The Panel allows the simultaneous inoculation of each cavity with a predetermined amount of inoculum. A suspension of the test organism in RapID Inoculation Fluid is used as the inoculum which rehydrates and initiates test reactions. After incubation of the panel, each test cavity is examined for reactivity by noting the development of a color. In some cases, reagents must be added to the test cavities to provide a color change. The resulting pattern of positive and negative test scores is used as the basis for identification of the test isolate by comparison with the probability values in the Differential Chart (Table 4), or by use of the RapID ERIC™ software.

3. PRINCIPLE

The tests used in the RapID CB Plus System are based on microbial degradation of specific substrates detected by various indicator systems. The reactions employed are a combination of conventional tests and single-substrate chromogenic tests, described in Table 1.

4. REAGENTS

RapID CB Plus Reagent (provided with kit)	(15 ml/Btl)
Reactive ingredient per liter:	
<i>p</i> -Dimethylaminocinnamaldehyde	0.06 g
RapID Inoculation Fluid (R8325106, supplied separately)	(2 ml/tube)
KCl	6.0 g
CaCl ₂	0.5 g
Deminerlized Water.....	1000.0 ml
RapID Nitrate A Reagent (R8309003, supplied separately)	(15 ml/Btl)
Sulfanilic Acid	8.0 g
Glacial Acetic Acid	280.0 ml
Deminerlized Water.....	720.0 ml
RapID Nitrate B Reagent (R8309004, supplied separately)	(15 ml/Btl)
N,N-Dimethyl-1-naphthylamine	6.0 g
Glacial Acetic Acid	280.0 ml
Deminerlized Water.....	720.0 ml

5. PRECAUTIONS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test panels after use. Directions should be read and followed carefully.

Non-disposable apparatus should be sterilised by any appropriate procedure after use, although the preferred method is to autoclave for 15 minutes at 121°C; disposables should be autoclaved or incinerated. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated area swabbed with a standard bacterial disinfectant or 70% alcohol. Do NOT use sodium hypochlorite. Materials used to clean spills, including gloves, should be disposed of as biohazardous waste.

Do not use reagents beyond the printed expiration dates.

Do not use if there is evidence of contamination or other signs of deterioration.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established. In the event of malfunction do not use device.

Caution!

- RapID CB Plus Reagent is toxic and may cause harm to the environment. Harmful by inhalation, contact with skin or eyes, or if swallowed. May impair fertility or cause harm to unborn child.
- RapID Nitrate A Reagent and RapID Nitrate B Reagent may cause irritation to skin, eyes, and respiratory system.

DANGER

US ONLY

US & EU

H315	Causes skin irritation
H319	Causes serious eye irritation
H335	May cause respiratory irritation
H336	May cause drowsiness or dizziness
H360	May damage fertility. May damage the unborn child
H373	May cause damage to organs through prolonged or repeated exposure
P201	Obtain special instructions before use
P202	Do not handle until all safety precautions have been read and understood
P281	Use personal protective equipment as required
P264	Wash face, hands and any exposed skin thoroughly after handling
P280	Wear eye/face protection
P260	Do not breathe dust/fume/gas/mist/ vapors/spray
P271	Use only outdoors or in a well-ventilated area
P308+313	IF exposed or concerned: Get medical attention/advice
P304+P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing
P302+P352	IF ON SKIN: Wash with plenty of soap and water
P332+P313	If skin irritation occurs: Get medical advice/attention
P362	Take off contaminated clothing and wash before reuse
P305+P351 +P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
P337+P313	If eye irritation persists: Get medical advice/attention
P405	Store locked up
P403+P233	Store in a well-ventilated place. Keep container tightly closed
P501	Dispose of contents/container to an approved waste disposal plant

- Refer to Safety Data Sheet, available on company website, and product labeling for information on potentially hazardous components, for detailed information on reagent chemicals.

Composition / information on ingredients

2-Methoxyethanol 109-86-4

Acetic acid 64-19-7

Hydrochloric acid 7647-01-0

WARNING! This product contains a chemical known in the State of California to cause birth defects or other reproductive harm.

Emergency Telephone Number
INFOTRAC - 24 Hour Number: 1-800-535-5053
Outside of the United States, call 24 Hour Number:
001-352-323-3500 (Call Collect)

6. STORAGE

2°C-8°C

RapID CB Plus System, RapID Nitrate A Reagent, and RapID Nitrate B Reagent should be stored in their original containers at 2-8°C until used. Allow products to equilibrate to room temperature before use. DO NOT interchange reagents among different RapID systems. Remove only the number of panels necessary for testing. Reseal the plastic pouch and promptly return to 2-8°C. Panels must be used the same day they are removed from storage. RapID Inoculation Fluid should be stored in its original container at room temperature (20-25°C) until used.

7. PRODUCT DETERIORATION

This product should not be used if (1) the expiration date has passed, (2) the plastic tray is broken or the lid is compromised, or (3) there are other signs of deterioration.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.⁹⁻¹¹

9. MATERIALS SUPPLIED

- 20 RapID CB Plus Panels
- 20 report forms
- RapID CB Plus Reagent (one plastic dropper-bottle containing reagent sufficient for 20 panels)
- 1 color guide
- 2 chipboard incubation trays
- Instructions for use (IFU).

10. CONTENTS SYMBOLS

CB Plus Panels	CB Plus Panels
Report Forms	RapID Report Forms
CB Plus Reagent	CB Plus Reagent
Incubation Trays	Incubation Trays

11. MATERIALS REQUIRED BUT NOT SUPPLIED

- Loop sterilization device
- Inoculating loop, swabs, collection containers
- Incubators, alternative environmental systems
- Supplemental media
- Quality control organisms
- Gram stain reagents
- Microscope slides
- Cotton swabs
- Catalase test reagent (3% hydrogen peroxide)
- RapID Inoculation Fluid - 2 ml (R8325106)
- McFarland #4 turbidity standard or equivalent (R20414)
- Pipettes
- RapID Nitrate A Reagent (R8309003)
- RapID Nitrate B Reagent (R8309004)
- ERIC (Electronic RapID Compendium, R8323600) (Optional)

12. PROCEDURE

Inoculum Preparation:

- Test organisms must be grown in pure culture and examined by Gram stain prior to use in the system.
Note: Spore-forming Gram-positive bacilli should NOT be tested on the RapID CB Plus System.
- Test organisms may be cultured on routine laboratory media commonly used for coryneforms. The following types of media are recommended: Non-selective media such as Trypticase Soy Agar with 5% sheep blood and Columbia blood agar.
Notes:
 - Plates used for inoculum preparation should be less than 72 hours old. Plates cultivated for 24 hours may be used if sufficient growth is available. Plates that are 72 hours old may be used for very slow-growing strains.
 - The use of media other than those recommended may compromise test performance.
- Perform a catalase test on the isolate. Record the result in the Catalase space on the panel lid.
- Observe the test isolate for the production of yellow pigment. Remove several colonies from the surface of the agar plate with a cotton or polyester-tipped swab and observe for pigment color on the swab. If the organism produces a yellow pigment, record a plus in the Yellow Pigment space on the panel lid.
Note: Only record yellow pigmentation. The production of any other color should be recorded as negative in the pigment space.
- Using a cotton swab or inoculating loop, suspend sufficient growth from the agar plate culture in RapID Inoculation Fluid (2 ml) to achieve a visual turbidity equal to a #4 McFarland turbidity standard or equivalent.

Notes:

- Suspensions significantly less turbid than a #4 McFarland may result in aberrant reactions.
- Bacterial suspensions that are slightly more turbid than a #4 McFarland standard will not affect test performance and are recommended for stock cultures and quality control strains. However, suspensions prepared with a turbidity far greater than a #4 McFarland standard may compromise test performance.
- Suspensions should be mixed thoroughly and vortexed if required.
- Suspensions should be used within 15 minutes of preparation.

- An agar plate may be inoculated for purity and any additional testing that may be required using a loopful of the test suspension from the inoculation fluid tube. Incubate the plate for at least 18-24 hours at 35-37°C.

Inoculation of RapID CB PLUS Panels:

- Peel back the lid of the panel over the inoculation port by pulling the tab marked "Peel to Inoculate" up and to the left.
- Using a pipette, gently transfer the entire contents of the Inoculation Fluid tube into the upper right-hand corner of the panel. Reseal the inoculation port of the panel by pressing the peel-back tab back in place.
- After adding the test suspension, and while keeping the panel on a level surface, tilt the panel back away from the test cavities at approximately 45° (see below).



Table 1. Principles and Components of the RapID CB Plus System


Cavity #	Test Code	Reactive Ingredient	Quantity	Principle	Bibliography #
1	GLU	Glucose	2%	Utilization of the carbohydrate substrate produces acidic products which lower the pH and change the indicator.	1-5
2	SUC	Sucrose	2%		
3	RIB	Ribose	2%		
4	MAL	Maltose	2%		
5	αGLU	<i>p</i> -Nitrophenyl-α,D-glucoside	0.1%	Enzymatic hydrolysis of the colorless aryl-substituted glycoside or phosphoester, releases yellow <i>o</i> - or <i>p</i> -nitrophenol which is detected by the formation of a yellow color.	1-3, 6, 7
6	βGLU	<i>p</i> -Nitrophenyl-β,D-glucoside	0.1%		
7	NAG	<i>p</i> -Nitrophenyl-n-acetyl-β,D-glucosaminide	0.1%		
8	GLY1	<i>p</i> -Nitrophenyl-glycoside	0.1%		
9	ONPG	<i>o</i> -Nitrophenyl-β,D-galactoside	0.1%		
10	PHS	<i>p</i> -Nitrophenyl phosphate	0.2%	Hydrolysis of the fatty acid ester releases acidic products which lower the pH and change the indicator.	5, 7
11	EST	Fatty acid ester	2%		
12	PRO	Proline-β-naphthylamide	0.05%	Enzymatic hydrolysis of the aryl-amide substrate releases free β-naphthylamine which is detected with the RapID CB Plus Reagent.	3, 6-8
13	TRY	Tryptophane-β-naphthylamide	0.05%		
14	PYR	Pyrrrolidine-β-naphthylamine	0.05%		
15	LGLY	Leucyl-glycine-β-naphthylamide	0.05%		
16	LEU	Leucine-β-naphthylamide	0.05%	Hydrolysis of urea produces basic products which change the pH indicator.	1-3, 5
17	URE	Urea	1.2%		
18	NIT	Potassium nitrate	0.6%	Utilization of nitrate results in the formation of nitrite which is detected by the addition of RapID Nitrate Reagents.	1, 3, 5, 9

- While tilted back, gently rock the panel from side to side to evenly distribute the inoculum along the rear baffles as illustrated below.

- While maintaining a level, horizontal position (best achieved by using the bench top against the reaction cavity bottoms), slowly tilt the panel forward toward the reaction cavities until the inoculum flows along the baffles into the reaction cavities (see below). This should evacuate all of the inoculum from the rear portion of the panel.

Note: If the panel is tilted too quickly, air may be trapped at the test cavity junction, restricting fluid movement.

- Return the panel to a level position. If necessary, gently tap the panel on the bench top to remove any air trapped in the cavities.

Notes:

- Examine the test cavities. Test cavities should appear bubble-free and uniformly filled. Slight irregularities in test cavity fills are acceptable and will not affect test performance. If the panel is grossly misfilled, a new panel should be inoculated and the misfilled panel discarded.
- Complete the inoculation of each panel receiving inoculation fluid before inoculating additional panels.
- Do not allow the inoculum to rest in the back portion of the panel for prolonged periods without completing the procedure.

Incubation of RapID CB Plus Panels:

Incubate inoculated panels at 35-37°C in a non-CO₂ incubator for at least 4 hours but not more than 6 hours. For ease of handling, panels may be incubated in the chipboard incubation trays provided with the kit.

Scoring of RapID CB Plus Panels:

CAT (catalase) and PIG (yellow pigment) are not performed in the RapID CB Plus panel. However, they are recorded as part of the microcode score.

- While firmly holding the RapID CB Plus panel on the bench top, peel off the label lid over the reaction cavities by pulling the

RapID CB Plus Panel Test Location

Cavity #	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
Test Code	GLU	SUC	RIB	MAL	αGLU	βGLU	NAG	GLY1	ONPG	PHS	EST	PRO	TRY	PYR	LGLY	LEU	URE	NIT
												RapID CB Plus Reagent						NIT A NIT B

Table 2. Interpretation of RapID CB Plus System Tests*

Cavity #	Test Code	Reagent	Reaction		Comments
			Positive	Negative	
Tests contained in the RapID CB Plus Panel:					
1	GLU	None	Yellow, gold, or light yellow-orange	Red or dark orange	Only the development of a distinct yellow, gold, or yellow-orange should be scored as positive. Any other shade is scored as negative.
2	SUC				
3	RIB				
4	MAL				
5	αGLU				
6	βGLU	None	Yellow	Clear, tan	Any development of a yellow color should be scored as positive. Shades of yellow color may be different from cavity to cavity. Do not compare the color of one cavity to another.
7	NAG				
8	GLY1				
9	ONPG				
10	PHS				
11	EST	None	Yellow or yellow-orange	Red or orange	Only the development of a distinct yellow or yellow-orange should be scored as positive.
12	PRO	RapID CB Plus Reagent	Purple, red, or dark pink	Yellow or orange	Only significant color development should be scored as positive. Very pale colors should be scored as negative.
13	TRY				
14	PYR				
15	LGLY				
16	LEU				
17	URE	None	Red or dark red-orange	Yellow or orange	Only the development of a distinct red or deep red-orange should be scored as positive.
18	NIT	RapID Nitrate A and B	Red or pink	Clear or very pale pink	Any red to medium pink color should be scored as positive.
Tests not contained in the RapID CB Plus Panel that are part of the microcode:					
19	CAT	3% Hydrogen Peroxide	Bubbles	No bubbles	Any significant bubbling should be scored as positive.
20	PIG	None	Yellow pigment	No pigment or any other color	Only yellow-pigmented isolates are scored as positive.

*NOTE: Panels should be read by looking down through the reaction cavities against a white background.

Table 3. Quality Control Chart for RapID CB Plus Panels

Organism	GLU	SUC	RIB	MAL	αGLU	βGLU	NAG	GLY1	ONPG	PHS	EST	PRO	TRY	PYR	LGLY	LEU	URE	NIT
<i>Paenibacillus polymyxa</i> ATCC® 842	+	+	+	+	+	+	+	V	+	+	+	–	–	–	–	–	–	+
<i>Corynebacterium pseudodiphtheriticum</i> ATCC® 10701	–	–	–	–	–	–	–	–	–	V	–	+	+	–	+	+	+	+
<i>Arcanobacterium pyogenes</i> ^a ATCC® 19411 (also known as <i>Trueperella pyogenes</i>)	+	–	V	V	+	–	V	+	+	–	V	+	V	+	V	+	–	–

+, positive; –, negative; V, variable

^a **Key indicator strains** demonstrate acceptable performance of the most labile substrate in the system and reactivity in a significant number of wells, according to Clinical and Laboratory Standards Institute recommendations for streamlined quality control.¹⁶

to batch. As a result, culture media may influence constitutive enzymatic activity of designated quality control strains. If quality control strain results differ from the patterns indicated, a subculture onto medium from a different batch or from another manufacturer will often resolve quality control discrepancies.

15. LIMITATIONS

- The use of RapID CB PLUS System and the interpretation of results requires the knowledge of a competent microbiologist, familiar with laboratory procedures, who is trained in general microbiological methods and who judiciously makes use of training, experience, specimen information, and other pertinent procedures before reporting the identification obtained using this system.
- The RapID CB Plus System must be used with pure cultures of test organisms. The use of mixed microbial populations or direct testing of clinical material without culture will result in aberrant results.
- The RapID CB Plus System is designed for use with the taxa listed in the RapID CB Plus Differential Chart. Gram-positive bacilli that form spores or Gram-negative bacilli should not be tested.
- Expected values listed for RapID CB Plus System tests may differ from conventional test results or previously reported information.
- The accuracy of the RapID CB Plus System is based on statistical use of a multiplicity of specially designed tests and an exclusive, proprietary database. The use of any single test found in the RapID CB Plus System to establish the identification of a test isolate is subject to the error inherent in that test alone.

16. PERFORMANCE CHARACTERISTICS

The RapID CB Plus System performance characteristics have been established by laboratory testing of clinical, reference, and stock cultures.¹²⁻¹⁵ Among 435 strains tested, 422 (97%) RapID CB Plus results agreed with the reference results. Among the 422 strains that agreed, 395 (91%) agreed without additional testing. Twenty-seven (6%) probability overlap microcodes were encountered which required additional testing for resolution of choices and complete identification. Seven (1.6%) strains provided questionable microcodes and six (1.4%) RapID CB Plus results did not agree with reference identifications.

17. BIBLIOGRAPHY

- Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey and Scott’s Diagnostic Microbiology. 9th ed. Mosby, St. Louis, MO.
- Coyle, M.B. and B.A. Lipsky. 1990. Clin. Microbiol. Rev. 3:224-

Table 4 - RapID CB Plus Differential Chart

Organism	GLU	SUC	RIB	MAL	αGLU	βGLU	NAG	GLY1	ONPG	PHS	EST	PRO	TRY	PYR	LGLY	LEU	URE	NIT	CAT	PIG
<i>Corynebacterium accolens</i>	95	11	29	0	0	0	0	0	0	0	0	99	1	1	0	0	7	51	99	0
<i>Corynebacterium afermentans subsp. afermentans</i>	0	0	0	0	0	0	0	0	0	88	38	86	26	0	3	47	1	1	98	0
<i>Corynebacterium afermentans subsp. lipophilum</i>	0	0	0	0	0	0	0	0	0	0	9	98	0	0	0	0	0	0	98	0
<i>Corynebacterium amycolatum</i>	97	41	70	80	4	0	0	0	0	97	88	99	93	0	9	11	18	76	98	0
<i>Corynebacterium argentoratense</i>	98	0	5	0	0	0	0	0	0	14	69	99	90	0	7	82	0	0	98	0
<i>Corynebacterium auris</i>	0	0	0	0	0	90	0	0	0	96	95	99	71	1	8	76	1	1	98	0
<i>Corynebacterium bovis</i>	28	8	3	13	0	0	0	2	99	19	9	90	9	0	9	71	8	0	98	0
<i>Corynebacterium cystitidis</i>	90	0	0	0	0	0	0	90	0	95	5	33	45	9	28	88	99	0	98	0
<i>Corynebacterium diphtheriae</i>	98	2	81	90	99	0	0	0	0	1	24	99	21	1	30	41	0	28	98	0
<i>Corynebacterium glucuronolyticum</i>	95	86	91	24	0	7	0	99	0	6	11	99	86	6	48	78	71	19	98	0
<i>Corynebacterium jeikeium</i> (JK)	88	2	78	9	1	1	0	0	0	85	92	0	58	1	41	69	1	2	98	0
<i>Corynebacterium kutscheri</i>	96	93	82	96	93	99	0	0	0	0	58	96	90	99	98	95	99	81	98	0
<i>Corynebacterium matruchotii</i>	70	62	74	65	2	80	2	2	20	8	15	90	82	93	91	0	0	98	98	0
<i>Corynebacterium minutissimum</i>	98	73	76	99	3	9	0	0	4	83	27	99	78	1	74	88	0	0	98	0
<i>Corynebacterium propinquum</i> (ANF3)	0	0	0	0	0	0	0	0	0	91	63	6	79	33	70	57	2	99	98	0
<i>Corynebacterium pseudodiphtheriticum</i>	0	0	0	0	1	0	0	0	0	22	9	96	66	31	72	82	99	86	98	0
<i>Corynebacterium pseudotuberculosis</i>	84	4	85	30	18	0	0	0	0	4	7	99	78	0	69	74	99	0	98	13
<i>Corynebacterium renale</i>	88	2	11	0	0	0	12	99	7	5	2	98	14	1	0	0	99	0	90	58
<i>Corynebacterium striatum</i>	97	96	49	0	0	2	0	0	1	81	37	99	90	6	70	82	1	98	98	0
<i>Corynebacterium ulcerans</i>	97	0	95	88	97	0	8	2	1	88	40	99	99	1	33	62	98	0	98	0
<i>Corynebacterium urealyticum</i> (D-2)	1	0	3	1	0	0	0	0	0	0	38	1	66	1	7	66	99	0	98	0
<i>Corynebacterium xerosis</i>	98	96	90	90	96	0	4	0	0	84	22	96	76	0	51	88	1	61	99	0
CDC Group F-1	98	1	96	88	0	0	0	0	0	1	1	98	79	0	2	39	99	24	96	0
CDC Group G (G/LD)	96	82	88	0	0	0	0	0	0	77	21	99	31	71	68	82	0	50	96	0
CDC Group 11	90	0	8	2	0	0	0	0	0	98	32	98	76	0	0	6	1	91	97	0
<i>Listeria grayi/murrayi</i>	95	4	95	90	98	99	99	56	0	1	92	1	47	3	8	12	1	58	98	0
<i>Listeria innocua</i>	96	2	0	92	98	99	99	31	0	1	96	0	92	0	12	8	0	0	98	0
<i>Listeria ivanovii</i>	96	12	33	95	76	95	99	0	0	2	88	1	37	0	8	21	1	0	98	0
<i>Listeria monocytogenes</i>	96	2	7	96	99	99	99	14	0	0	96	0	80	0	1	9	1	0	96	0
<i>Listeria seeligeri</i>	96	2	0	95	98	98	98	0	0	4	96	1	42	0	10	22	1	0	95	0
<i>Listeria welshimeri</i>	96	5	1	96	69	99	99	73	0	2	96	1	81	0	15	33	2	0	95	0
<i>Actinomyces israelii</i>	98	95	92	13	99	97	2	4	87	2	8	98	64	19	86	90	2	31	0	0
<i>Actinomyces naeslundii</i>	90	82	13	90	86	96	0	2	91	11	2	98	65	1	71	87	94	77	0	0
<i>Actinomyces viscosus</i>	96	98	62	96	92	90	0	7	91	5	6	92	1	0	76	88	82	96	98	0
<i>Arcanobacterium haemolyticum</i>	95	9	88	83	98	26	92	99	84	88	86	81	66	72	78	90	2	3	0	0
<i>Arcanobacterium pyogenes</i> (also known as <i>Trueperella pyogenes</i>)	92	24	82	77	98	0	69	89	90	53	18	90	3	2	83	90	1	0	0	0
<i>Bifidobacterium</i> spp. (Gp E)	91	90	39	90	95	74	15	0	43	2	44	83	16	5	11	3	1	0	0	0
<i>Brevibacterium casei</i>	0	0	5	1	93	0	35	0	0	49	16	99	88	28	78	86	1	7	99	9
<i>Brevibacterium</i> spp. (Gp B)	1	0	9	0	4	3	29	1	1	88	19	96	86	11	76	90	5	9	95	16
<i>Cellulomonas</i> spp. (Gp A3/A4)	98	98	66	98	98	99	90	7	88	96	81	89	80	96	75	80	1	84	96	76
<i>Corynebacterium otitidis</i> ^c	0	0	0	0	0	0	0	0	0	95	92	99	99	0	96	57	0	0	98	0
<i>Dermabacter hominis</i>	95	90	77	88	99	91	98	92	95	71	63	29	92	98	88	77	0	0	98	0
<i>Erysipelothrix rhusiopathiae</i>	88	0	9	0	0	12	99	2	93	80	2	85	73	93	30	19	0	0	3	0
<i>Exiguobacterium acetylicum</i>	98	98	78	98	99	98	1	0	91	20	1	0	2	0	0	0	1	0	99	92
<i>Leifsonia aquatica</i> ^b	28	0	8	9	99	94	96	88	16	4	2	98	88	2	91	87	1	8	98	95
<i>Microbacterium</i> spp.	98	98	87	92	99	90	98	1	96	66	36	8	77	5	80	29	1	5	98	86
<i>Oerskovia</i> spp.	98	97	41	90	96	98	90	33	81	11	7	98	92	28	79	90	1	45	98	94
<i>Rhodococcus equi</i>	18	2	12	2	80	82	0	2	18	94	96	91	88	4	35	87	61	55	95	0
<i>Rothia</i> spp.	92	90	0	77	99	55	0	0	0	3	1	98	90	88	91	96	0	81	98	0
<i>Schaalia odontolyticus</i> ^c	92	90	12	65	81	66	3	2	86	14	38	98	13	0	92	82	1	91	0	0
<i>Trueperella bernardiae</i> ^a	60	0	40	45	96	0	99	0	0	5	0	90	15	4	63	74	1	0	0	0
<i>Winkia neuii</i> ^d	74	59	62	54	99	2	2	74	92	1	57	99	13	1	40	98	1	18	96	0

^aPreviously designated as *Arcanobacterium bernardiae* and prior to that *Actinomyces bernardiae* (group 2)

^bPreviously designated as *Corynebacterium aquaticum*

^cPreviously designated as *Actinomyces odontolyticus*

^dPreviously designated as *Actinomyces neuii* (Gp 1)

^ePreviously designated as *Turicella otitidis*

Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA

www.thermofisher.com/microbiology
Tel: (800) 255-6730 • International: (913) 888-0939

www.oxoid.com/IFU
Europe +800 135 79 135 • US 1 855 2360 190
CA 1 855 805 8539 • ROW +31 20 794 7071

Version	Date of modifications introduced
IFU8311008	August 2023 Updated to meet IVDR requirements

Printed in the UK

18. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
	Unique Device Identifier
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment
	Manufacturer

RapID™ is a trademark of Thermo Fisher Scientific and its subsidiaries.
ERIC™ is a trademark of Thermo Fisher Scientific and its subsidiaries.
ATCC® is a registered trademark of American Type Culture Collection.

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Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, SUA