Egg Yolk Emulsion

REF SR0047C

EN

*See additional IFUs for Blood Agar Base CM0055B, *Bacillus cereus* Selective Agar Base CM0617B

Intended Use

Egg Yolk Emulsion (SR0047C) is an enrichment supplement intended for use in culture media. It is intended to be added directly to nutrient media to aid the identification of organisms such as *Clostridium*, *Bacillus* and *Staphylococcus* species, by their lipase activity. Egg Yolk Emulsion (SR0047C), when added to *Bacillus cereus* Selective Agar Base (CM0617B), is intended to be used for the isolation of *Bacillus cereus* from faecal samples. Egg Yolk Emulsion (SR0047C), when added to Blood Agar Base (CM0055B and CM0055T), is intended to be used for the isolation of *Staphylococcus* and *Clostridium* species from clinical samples, including wound swabs.

The device is for professional use only, is not automated and nor is it a companion diagnostic.

Summary and Explanation

Over 260 species and several subspecies comprise the genus *Bacillus*, including *Bacillus subtilis* and *Bacillus cereus*. Although they are often found in the environment and as laboratory contaminants, there are a few which are known to cause infections in humans, such as *B. cereus* which cause foodborne illness, which in turn causes severe nausea, vomiting and diarrhoea^{1,2}. Moreover, *B. cereus* shows resistance to most β -lactam antibiotics, with resistance to penicillin being an identifying feature of *B. cereus*¹.

Staphylococci are Gram-positive cocci, widespread in nature and generally non-pathogenic³. However, some *staphylococci* can cause disease in humans³. Infection with *S. aureus* can cause infections ranging from superficial skin infections, such as cellulitis and abscesses, to life-threatening illnesses such as pneumonia, toxic shock syndrome and sepsis³.

Clostridium species are Gram-positive spore-forming anaerobic bacilli which are linked to serious human infections⁴. *Clostridium perfringens* is the most commonly isolated *Clostridium* species⁴. *C. perfringens* is a possible cause of acute gastrointestinal infections, which can range in symptoms from diarrhoea to gas gangrene and necrotizing enterocolitis⁴.

Egg Yolk Emulsion devices can be added to culture media to help detect lecithinase, lipase, and proteolytic activity of microorganisms such as *Bacillus cereus*, *Staphylococcus* species and *Clostridium* species. Lecithinase positive strains produce zones of opacity around the colonies, whilst lecithinase negative strains do not produce zones of opacity.

Principle of Method

Egg Yolk Emulsion (SR0047C) can be added to culture media to help detect lecithinase, lipase, and proteolytic activity of microorganisms such as *Bacillus cereus*, *Staphylococcus aureus* and *Clostridium perfringens*.

Typical Formula

	Milligrams per litre
10% v/v Egg Yolk Emulsion	100.0ml

Materials Provided

SR0047C: 1 x 100ml of Egg Yolk Emulsion

Materials Required but Not Supplied

Media such as Blood Agar Base CM0055B, *Bacillus cereus* Selective Agar Base CM0617B

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Storage

- Store product in its original packaging between 2°C and 8°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or vials.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 7, ID 8, ID 9 and S7.

Procedure

- Add the required volume of Egg Yolk Emulsion (SR0047C) to sterile molten medium that has been cooled to 50°C. For example, 50mL should be added to 500ml of sterile Blood Agar Base (CM0055T) or *Bacillus cereus* Selective Agar Base (CM0617B), prepared as directed. and cooled to 50°C.
- Mix well and aseptically dispense into sterile containers.

Refer to the instructions for use for CM0055T, CM0617B and other media being used. IFUs can be found at www.thermofisher.com


Bibliography

1. Public Health England. 2018. 'Identification of *Bacillus* species.' UK Standards for Microbiology Investigations ID 9 Issue 3.1. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-9-identification-of-bacillus-species>.
2. Public Health England. 2020b. 'Gastroenteritis.' UK Standards for Microbiology Investigations S 7 Issue 2. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>.
3. Public Health England. 2020a. 'Identification of *Staphylococcus* species, *Micrococcus* species and

Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.

4. Public Health England. 2016. 'Identification of Clostridium species.' UK Standards for Microbiology Investigations ID 8 Issue 4.1. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

ATCC Licensed
Derivative

©2022 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road,
Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2023-12-04

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Bibliografija

1. Public Health England. 2018. 'Identification of Bacillus species.' UK Standards for Microbiology Investigations ID 9 Issue 3.1. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-9-identification-of-bacillus-species>.
2. Public Health England. 2020b. 'Gastroenteritis.' UK Standards for Microbiology Investigations S 7 Issue 2. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>.
3. Public Health England. 2020a. 'Identification of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2016. 'Identification of Clostridium species.' UK Standards for Microbiology Investigations ID 8 Issue 4.1. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>.

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Oxoid Limited, Wade Road, Basingstoke,
RG24 8PW, UK

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Litteratur

1. Public Health England. 2018. 'Identification of Bacillus species.' UK Standards for Microbiology Investigations ID 9 Issue 3.1. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-9-identification-of-bacillus-species>.
2. Public Health England. 2020b. 'Gastroenteritis.' UK Standards for Microbiology Investigations S 7 Issue

2. Accessed 21 July 2022.
<https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>.
3. Public Health England. 2020a. 'Identification of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. Accessed 21 July 2022.
<https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. Public Health England. 2016. 'Identification of Clostridium species.' UK Standards for Microbiology Investigations ID 8 Issue 4.1. Accessed 21 July 2022.
<https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>.

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
www.thermofisher.com

Thermo
SCIENTIFIC

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10 mahuprotsendiline Egg Yolk Emulsion

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Bibliográfia

- Public Health England. 2018. 'Identification of Bacillus species.' UK Standards for Microbiology Investigations ID 9 Issue 3.1. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-9-identification-of-bacillus-species>.
- Public Health England. 2020b. 'Gastroenteritis.' UK Standards for Microbiology Investigations S 7 Issue 2. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>.
- Public Health England. 2020a. 'Identification of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
- Public Health England. 2016. 'Identification of Clostridium species.' UK Standards for Microbiology Investigations ID 8 Issue 4.1. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>.

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Literatūra

1. Public Health England. 2018. 'Identification of *Bacillus* species.' UK Standards for Microbiology Investigations ID 9 Issue 3.1. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-9-identification-of-bacillus-species>.
2. Public Health England. 2020b. 'Gastroenteritis.' UK Standards for Microbiology Investigations S 7 Issue 2. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>.

- Public Health England. 2020a. 'Identification of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
- Public Health England. 2016. 'Identification of Clostridium species.' UK Standards for Microbiology Investigations ID 8 Issue 4.1. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>.

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Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, JK

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Zdroje

1. Public Health England. 2018. 'Identification of *Bacillus* species.' UK Standards for Microbiology Investigations ID 9 Issue 3.1. Accessed 21 July 2022.

- <https://www.gov.uk/government/publications/smi-id-9-identification-of-bacillus-species>.
- Public Health England. 2020b. 'Gastroenteritis.' UK Standards for Microbiology Investigations S 7 Issue 2. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>.
 - Public Health England. 2020a. 'Identification of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
 - Public Health England. 2016. 'Identification of Clostridium species.' UK Standards for Microbiology Investigations ID 8 Issue 4.1. Accessed 21 July 2022. <https://www.gov.uk/government/publications/smi-id-8-identification-of-clostridium-species>.

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Oxoid Limited, Wade Road,
Basingstoke, Hampshire, RG24
8PW, UK

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK