Yersinia Selective Supplement

REF SR0109E

EN

This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Oxoid Yersinia Selective Agar Base (product codes: CM0653B and CM0653T).

Intended Use

Yersinia Selective Supplement (SR0109E) is a selective supplement that is intended to be used in the preparation of Yersinia Selective Agar Base (Dehydrated) (CM0653B and CM0653T) for the isolation of *Yersinia enterocolitica* from faecal samples.

Yersinia selective Supplement, when added to Yersinia Selective Agar Base (Dehydrated), is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having yersiniosis.

The device is for professional use only, is not automated and nor is it a companion diagnostic.

Summary and Explanation

Yersinia enterocolitica are Gram-negative, facultatively anaerobic, coccoid shaped bacteria and are the main causative agent of the foodborne zoonotic disease yersiniosis.^{1,2} In humans, *Y. enterocolitica* have primarily been isolated from faecal samples, but they have also been found in blood and lymph node tissues.¹

Y. enterocolitica infections are recognized as a main cause of yersiniosis.² In 2019, yersiniosis was the fourth most reported foodborne zoonotic disease in Europe, with more than 7,000 confirmed cases, with 99% of these infections caused by *Y. enterocolitica*.³ Infections present with a range of manifestations, including gastroenteritis, acute diarrhoea, mesenteric adenitis, terminal ileitis, 'pseudo-appendicitis', and in rare cases sepsis.^{1,2} Consequently, it is important to isolate *Yersinia enterocolitica* in clinical samples as it is a considerable human pathogen associated with yersiniosis.^{1,2}

Isolation and presumptive identification of *Y. enterocolitica* is conventionally carried out using a combination of microscopic appearance, cultural appearance on *Y. enterocolitica* selective agar, growth requirements and serotyping.

Principle of Method

Yersinia Selective Supplement (SR0109E) when used with Yersinia Selective Agar Base (Dehydrated) (CM0653B and CM0653T) is intended for the isolation of *Yersinia enterocolitica*, through the inclusion of a combination of selective agents. Novobiocin helps inhibit the growth of most Gram-positive microorganisms, as well as some strains of *Proteus*. Cefsulodin is a bactericidal antimicrobial agent with activity against *Pseudomonas aeruginosa*, whilst Irgasan is a broad-spectrum antimicrobial agent.

Typical Formula

Cefsulodin	7.5mg
Irgasan	2.0mg
Novobiocin	1.25mg

Materials Provided

SR0109E: 10x freeze dried vials, each to prepare 500ml of medium.

Materials Required but Not Supplied

- CM0653B: 500g dehydrated Yersinia Selective Agar Base.
- CM0653T: 5kg dehydrated Yersinia Selective Agar Base.
- Inoculating loops.
- Swabs.
- Collection containers.
- Petri dishes.

Storage

- Store product in its original packaging at 2–8°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

Signal Word: Danger

Hazard Statements

- H317 - May cause an allergic skin reaction
- H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled
- H319 - Causes serious eye irritation
- H400 - Very toxic to aquatic life
- H410 - Very toxic to aquatic life with long lasting effects

Precautionary Statements

P280 - Wear protective gloves/protective clothing
 P302 + P352 - IF ON SKIN: Wash with plenty of soap and water
 P333 + P313 - If skin irritation or rash occurs: Get medical advice/attention
 P285 - In case of inadequate ventilation wear respiratory protection
 P342 + P311 - If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician
 P273 - Avoid release to the environment
 P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing

Each vial is for single use. Do not re-use.

For in vitro diagnostic use only.

For professional use only.

Inspect the product packaging before first use.

Do not use the product if there is any visible damage to the packaging or vials.

Do not use the product beyond the stated expiry date.

Do not use the device if signs of contamination are present. Do not use the device if the colour has changed or there are other signs of deterioration.

It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local the recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 21 and S 7.

Procedure


- Allow supplement to equilibrate to room temperature before use.
- Aseptically add 2ml of 1:1 ethanol:sterile distilled water to 1 vial and mix gently to dissolve.
- Aseptically add the vial contents to 500ml of sterile Yersinia Selective Agar Base (CM0653) prepared as directed and cooled to 50°C.
- Mix well and aseptically dispense into sterile containers.


Refer to the IFU for CM0653B and CM0653T. IFUs can be found at www.thermofisher.com.

Bibliography

1. Public Health England. 2015. 'Identification of *Yersinia* species'. UK Standards for Microbiology Investigations. ID 21 Issue 3. Accessed 11 February 2022.
2. Public Health England. 2020. 'Gastroenteritis'. UK Standards for Microbiology Investigations. S 7 Issue 2. Accessed 11 February 2022.
3. European Centre for Disease Prevention and Control. 2021. 'Yersiniosis - Annual epidemiological report for 2019'. ECDC. Accessed 11 February 2022.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use

	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

 Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Modifications introduced
3.0	2024-07-16

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
1. Public Health England. 2015. 'Identification of *Yersinia* species'. UK Standards for Microbiology Investigations. ID 21 Issue 3. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-id-21-identification-of-yersinia-species-from-faeces>
2. Public Health England. 2020. 'Gastroenteritis'. UK Standards for Microbiology Investigations. S 7 Issue 2. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>
3. European Centre for Disease Prevention and Control. 2021. 'Yersiniosis - Annual epidemiological report for 2019'. ECDC. Accessed 11 February 2022. <https://www.ecdc.europa.eu/en/publications-data/yersiniosis-annual-epidemiological-report-2019>

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, H.B.

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Bibliografia

1. Public Health England. 2015. 'Identification of Yersinia species'. UK Standards for Microbiology Investigations. ID 21 Issue 3. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-id-21-identification-of-yersinia-species-from-faeces>
2. Public Health England. 2020. 'Gastroenteritis'. UK Standards for Microbiology Investigations. S 7 Issue 2. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>
3. European Centre for Disease Prevention and Control. 2021. 'Yersiniosis - Annual epidemiological report for 2019'. ECDC. Accessed 11 February 2022. <https://www.ecdc.europa.eu/en/publications-data/yersiniosis-annual-epidemiological-report-2019>

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.