Nitrocefin Solution

REF SR0112C

EN

Intended Use

Nitrocefin Solution (SR0112C) is a reagent for the rapid detection of beta-lactamase production in microorganisms. Nitrocefin Solution (SR0112C) is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of being infected with microorganisms which possess beta-lactamase activity. The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Nitrocefin degradation should be used to give a rapid indication of beta-lactam inactivating systems and the result so obtained will, in most cases, predict the outcome of susceptibility tests with beta-lactam anti-microbials. However, it should not entirely replace conventional susceptibility testing as other factors also influence the results of such tests, and on occasion intrinsic resistance to beta-lactam antimicrobials has not been correlated with production of beta-lactamase¹. Nitrocefin degradation has been found to be highly efficient in detecting beta-lactamase producing isolates of *Neisseria gonorrhoeae*^{2,3}, *Haemophilus influenzae*^{2,4,5,6} and *staphylococci*^{5,6}.

Excellent results have also been obtained with certain anaerobic bacteria, notably with *Bacteroides* species^{8,9,10}. It should be emphasized that the efficacy of the Nitrocefin tests in predicting the beta-lactam susceptibilities of other microorganisms is at present unproven. Another chromogenic cephalosporin, PADAC (Hoechst-Roussel) was not as effective as Nitrocefin in detecting *staphylococcal* beta-lactamase⁷.

Principle of Method

Nitrocefin is the chromogenic cephalosporin. This compound exhibits a rapid distinctive colour change from yellow (max at pH 7.0 = 390nm) to red (max at pH 7.0 = 486nm) as the amide bond in the beta-lactam ring is hydrolysed by a beta-lactamase, it is sensitive to hydrolysis by all known lactamases produced by Gram-positive and Gram-negative bacteria. This characteristic reaction forms the basis of several methods suitable for diagnostic use.

Typical Formula

SR112A (Rehydration Fluid)	
Phosphate buffer	1.9 ml
Dimethylsulphoxide	0.1 ml
SR0112B (Nitrocefin)	
Nitrocefin	1.0 mg/l

Materials Provided

5 vials lyophilized nitrocefin (1mg) (SR0112B) and 5 vials rehydrating fluid (2mL) (SR112A).

Materials Required but Not Supplied

- Inoculated media
- Glass slides
- Incubator
- Inoculating loop
- Sonicator
- Whatman No. 1 filter paper disc
- Petri dish

Storage

- Store product in its original packaging between -20°C and 8°C.
- When prepared, the solution may be stored at -20°C for up to two weeks.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Warnings and Precautions

- IF ON SKIN: Wash with plenty of soap and water.
- If skin irritation or rash occurs: Get medical advice/attention.

- IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- In case of inadequate ventilation wear respiratory protection.
- IF INHALED: Remove person to fresh air and keep at rest in a position comfortable for breathing.
- If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician
- Avoid release to the environment.
- DO NOT USE SUPPLEMENTS IF YOU ARE PREGNANT OR CONTEMPLATING PREGNANCY.
- Each vial is for single use. Do not re-use.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or vials.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.

It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage There are no specimen collection and/or transport materials provided with this device. Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) Q 5.

Procedure

Reconstitution

- Allow supplement to equilibrate to room temperature before use.
- Reconstitute the contents of one vial of lyophilised Nitrocefin SR0112B by adding the entire contents (2ml) of one vial of rehydration fluid SR112A. This yields a working Nitrocefin solution of 500mg/ml, (approx 10-3 M) suitable for most applications

Direct Plate Method

- Add one drop of the Nitrocefin solution on to the surface of the colony. If the isolate is a high beta-lactamase producer then the colony and the surrounding area will quickly turn red. To detect a weak beta-lactamase producer the plate should then be incubated for 30 minutes before being reported as negative.

Slide Method

- Add one drop of the Nitrocefin solution on to a clean glass slide.
- Using a sterile loop, pick one colony from the plate and emulsify into the Nitrocefin drop. Report as positive if the colour changes from yellow to red within 30 minutes (protect the slide from desiccation during the waiting period).

Broth Method

- Add four drops of Nitrocefin solution to 1ml of the grown culture. Report as positive if the colour changes to red within 30 minutes.

Broken Cell Method

- Sonicate 1ml of the culture in order to break open the cells. Add 4 drops of Nitrocefin solution. Report as positive if the colour changes to red within 30 minutes.

Paper Disc Spot Test


- A Whatman No.1 filter paper disc (diameter 7cm) is placed in a petri dish and impregnated with Nitrocefin solution (0-5ml). This impregnated paper is generally usable for one day but should be kept away from light to avoid spontaneous degradation.
- An isolated colony is applied to the impregnated paper with a loop; a pink to red reaction developing within 15 minutes indicates beta-lactamase presence.
- Detection of beta-lactamase activity on gels
- beta-lactamase activities from bacterial cells and the technique for analytical iso-electric focusing have been described by Matthew et al.¹¹ The developed gels are stained by applying Whatman No.54 paper impregnated with the Nitrocefin solution¹¹. Focused bands in the gel with beta-lactamase activity appear pink on a yellow background.
- Determination of beta-lactamase activity by spectrophotometric assay. The working solution of Nitrocefin (500µg/ml) is diluted tenfold in buffer (0.1M phosphate; 1mM EDTA, pH 7.0). Spectrophotometric assays for

beta-lactamase are carried out measuring changes in wavelength at 486nm. The molar extinction coefficient of Nitrocefin at this wavelength is 20,500. Test samples of the finished product for performance with control cultures.

Bibliography

1. Public Health England. 2018. "Identification of Campylobacter Species." UK Standards for Microbiology Investigations.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Made in the United Kingdom
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union

©2025 Thermo Fisher Scientific Inc. All rights reserved.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
3.0	09-June-2025 Minor formatting changes

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©2025 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.