<https://www.thermofisher.com>

EN

Oxoid™ Amoxycillin/Clavulanic Acid Discs 30 µg/3 µg (AMC30/AMC3)

REF CT0223B / CT0538B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Amoxycillin/Clavulanic Acid Antimicrobial Susceptibility Test (AST) Discs are 6mm paper discs that contain specific amounts of amoxycillin and clavulanic acid. The discs are clearly marked on both sides. The discs are labelled on both sides with details of the antimicrobial (AMC) and amount present (3 µg or 30 µg).

Discs are supplied in cartridges containing 50 discs each, and there are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Amoxycillin/Clavulanic Acid Antimicrobial Susceptibility Test (AST) Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which amoxycillin and clavulanic acid have been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA³, CLSI M100^{1ac}, or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Negative

- *Enterobacteriaceae*
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*

Gram-Positive

- N/A

Species with published breakpoints according to current CLSI literature:

Gram-Negative

- *Enterobacterales*
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*

Gram-Positive

- N/A

Species with published breakpoints according to current EUCAST literature:

Gram-Negative

- *Enterobacterales*
- *Enterobacterales* (uncomplicated UTI only)
- *Haemophilus influenzae*
- *Moraxella catarrhalis*
- *Pasteurella multocida*
- *Burkholderia pseudomallei*

Gram-Positive

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Amoxycillin/Clavulanic AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Amoxycillin/Clavulanic Acid AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Amoxycillin/Clavulanic Acid AST Discs consist of 6mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Amoxycillin/Clavulanic AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to table 2.

Table 1. Materials Provided with CT0223B/CT0538B		Table 2. Description of Amoxycillin/Clavulanic Acid AST Disc Reagents	
Component description	Material Description	Reagent	Description of Function
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.	Amoxycillin	White or almost white crystalline powder (milled). Antimicrobial agent. It is effective against susceptible bacteria during the stage of active multiplication. It acts through the inhibition of cell wall biosynthesis that leads to the death of the bacteria. It involves inhibition of bacterial protein synthesis by binding to 30S ribosomes.
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.		
Foil	Foil individually sealing each cartridge with its desiccant.	Potassium Clavulanate	White to pale brown powder. Clavulanic acid is a β -lactam drug that functions as a mechanism-based β -lactamase inhibitor. While not effective by itself as an antibiotic, when combined with penicillin-group antibiotics, it can overcome antibiotic resistance in bacteria that secrete β -lactamase, which otherwise inactivates most penicillins.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.		

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Amoxycillin/Clavulanic Acid AST Discs can be dispensed using a Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Amoxycillin/Clavulanic Acid AST discs have a shelf-life of 12 (CT0223B) and 36 (CT0538B) months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 1. Raw data of zone sizes, in accordance with CLSI methodology, taken from 3 batches of Amoxycillin/Clavulanic Acid (30 µg) AMC30 Antimicrobial Susceptibility Discs (CT0223B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3230872	<i>Escherichia coli</i>	35218	MHA	17	22	19.5	19	19	19	-0.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	28	36	32	33	33	33	1 ± 0.58 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	15	23	19	20	21	21	1.67 ± 0 (CV=0.35)
	<i>Escherichia coli</i>	25922	MHA	18	24	21	20	20	20	-1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3225171	<i>Escherichia coli</i>	35218	MHA	17	22	19.5	18	18	18	-1.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	28	36	32	33	33	33	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	15	23	19	20	20	20	1 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	18	24	21	18	18	18	-3 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3207304	<i>Escherichia coli</i>	35218	MHA	17	22	19.5	18	18	18	-1.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	28	36	32	34	33	33	1.34 ± 0.58 (CV=0.43)
	<i>Haemophilus influenzae</i>	49247	HTM	15	23	19	18	19	19	-0.34 ± 0.58 (CV=1.73)
	<i>Escherichia coli</i>	25922	MHA	18	24	21	18	18	18	-3 ± 0 (CV=0)

Table 2. Raw data of zone sizes, in accordance with EUCAST methodology, taken from 3 batches of Amoxicillin/Clavulanic Acid (30 µg) AMC30 Antimicrobial Susceptibility Disks Discs (CT0223B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3001661	<i>Escherichia coli</i>	25922	MHA	18	24	21	18	18	18	-1 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218	MHA	17	22	19.5	18	18	18	-0.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3001662	<i>Escherichia coli</i>	25922	MHA	18	24	21	18	18	18	-3 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218	MHA	17	22	19.5	18	18	18	-1.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3001695	<i>Escherichia coli</i>	25922	MHA	18	24	21	18	18	18	-3 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218	MHA	17	22	19.5	18	18	18	-1.5 ± 0 (CV=0)

Table 3. Raw data of zone sizes, in accordance with EUCAST methodology, taken from 3 batches of Amoxicillin/Clavulanic Acid (3 µg) AMC3 Antimicrobial Susceptibility Disks Discs (CT0538B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3221260	<i>Haemophilus influenzae</i>	49766	HTM	17	23	20	21	21	21	1 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	29213	MHA	19	25	22	21	21	21	-1 ± 0 (CV=0)

Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3190104	<i>Haemophilus influenzae</i>	49766	HTM	17	23	20	17	17	17	-3 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	29213	MHA	19	25	22	20	20	20	-2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3164310	<i>Haemophilus influenzae</i>	49766	HTM	17	23	20	19	19	19	-1 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	29213	MHA	19	25	22	20	20	20	-2 ± 0 (CV=0)

Clinical Performance Characteristics

A study⁴ evaluated quality of Amoxycillin/Clavulanic Acid AST on two occasions (2014 and 2017). Amoxycillin/Clavulanic Acid AST discs produced mean zone size values within 1mm of the expected target value.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established. In the event of malfunction do not use device.

References

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ahman J, Matuschek E, Kahlmeter G. The quality of antimicrobial discs from nine manufacturers—EUCAST evaluations in 2014 and 2017. Clinical Microbiology and Infection. 2019 Mar;25(3):346–52.


© 2025 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number

	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period-after-opening – the number of days the product can be used after opening
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged
	Authorized representative in the European Community/ European Union
	Unique device identifier
	European Conformity Mark
	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

Revision Information

Version	Date of issue and modifications introduced
3.0	2025-05-01 Typographical error corrected

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- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Åhman J, Matuschek E, Kahlmeter G. The quality of antimicrobial discs from nine manufacturers—EUCAST evaluations in 2014 and 2017. *Clinical Microbiology and Infection*. 2019 Mar;25(3):346–52.

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Viited

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ähman J, Matuschek E, Kahlmeter G. The quality of antimicrobial discs from nine manufacturers—EUCAST evaluations in 2014 and 2017. Clinical Microbiology and Infection. 2019 Mar;25(3):346–52.

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