	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Ampicillin Discs 2 µg/10 µg (AMP2/AMP10)

REF CT0002B / CT0003B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Thermo Scientific™ Oxoid™ Ampicillin Antimicrobial Susceptibility Test (AST) Discs are 6mm discs prepared by impregnating high quality absorbent paper with specified concentrations of ampicillin. The discs are clearly marked on both sides with the corresponding concentration code – AMP2 or AMP10. The code designates the antimicrobial agent (ampicillin) and drug content (2 µg or 10 µg).

Discs are supplied in cartridges containing 50 discs each, and there are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack.

Oxoid Ampicillin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disk should only be used once.

Intended Use

Oxoid Ampicillin AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which ampicillin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

Published clinical breakpoints in the current version of the FDA, CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Streptococcus* spp. (β- Hemolytic Group)
- *Enterococcus* spp.

Gram-Negative

- Enterobacteriaceae
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*

Species with published breakpoints according to current CLSI M100 literature (all associated with **10 µg**):

Gram-Positive

- *Enterococcus* spp.
- *Streptococcus* spp. β- Haemolytic Group
- *Streptococcus* spp. Viridans Group

Gram-Negative

- Enterobacterales
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Staphylococcus* spp. (*Staphylococcus saprophyticus*)
- *Enterococcus* spp.
- *Streptococcus pneumoniae* (indications other than meningitis)
- Viridans group *Streptococci*
- *Listeria monocytogenes* (intravenous)
- *Aerococcus sanguinicola* and *urinae*

Gram-Negative

- Enterobacterales

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

- *Haemophilus influenzae* (indications other than meningitis)

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Oxoid Ampicillin AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. A suitable therapeutic agent for in vivo use can be determined using filter paper discs impregnated with specified concentrations of antimicrobial agents placed on the surface of a suitable test medium. Pure cultures of clinical isolates are inoculated onto the test medium and the antimicrobial susceptibility test (AST) disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognized zone diameter ranges for the specific antimicrobial agent(s)/organism(s) combination under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Ampicillin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests. Calibrators or control materials are not sold with the device. Users are instructed to run Positive Control tests with each run to ensure correct performance.

Materials Provided

Oxoid Ampicillin AST Discs consist of 6mm diameter paper discs impregnated with a specific antimicrobial concentration. The discs are marked on both sides with an alpha-numeric code identifying the antimicrobial agent and concentration. Oxoid Ampicillin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0002B/CT0003B.	
Component Description	Material Description
Cartridge	K-Resin plastic cartridge for discs.
Spring/Cap/Plunger	Assembly components (cartridge).
Desiccant	Molecular sieve tablets 500 mg. Light beige to brown small lozenge shaped FBE tablets.
Foil	Hard temper foil. Unprinted.

Table 2. Description of Ampicillin AST Disc Reagent	
Reagent	Description of Function
Ampicillin	White to off white solid. Antimicrobial agent. Ampicillin is a beta-lactam penicillin that binds penicillin-binding proteins inside the bacterial cell wall, subsequently inhibiting cell wall synthesis. This ultimately causes cell lysis.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Ampicillin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Ampicillin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2 to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this document.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Ampicillin (10 µg) AMP10 Antimicrobial Susceptibility Test Discs (CT0003B).

CLSI Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3282779	<i>Staphylococcus aureus</i>	25923	MHA	27	35	31	34	34	34	3 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	15	22	18.5	21	21	21	1.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218	MHA	6	6	6	6	6	6	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	13	21	17	19	19	19	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5% SB	30	36	33	36	36	36	3 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
3304656	<i>Staphylococcus aureus</i>	25923	MHA	27	35	31	33	33	33	2 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	15	22	18.5	20	20	20	1.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218	MHA	6	6	6	6	6	6	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	13	21	17	19	19	19	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5% SB	30	36	33	35	35	35	2 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
3285642	<i>Staphylococcus aureus</i>	25923	MHA	27	35	31	30	30	30	-1 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	15	22	18.5	19	19	19	0.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218	MHA	6	6	6	6	6	6	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	13	21	17	18	18	18	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5% SB	30	36	33	36	36	36	3 ± 0 (CV=0)

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the 3 most recent batches of Ampicillin (2 µg) AMP2 Antimicrobial Susceptibility Test Discs (CT0002B).

EUCAST Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3326457	<i>Staphylococcus aureus</i>	29213	MHA	15	21	18	16	16	16	-2 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	29212	MHA	15	21	18	19	19	19	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	MHA-F	6	12	9	6	6	6	-3 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49766	MHA-F	19	25	22	24	24	24	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	25	31	28	28	28	28	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3326455	<i>Staphylococcus aureus</i>	29213	MHA	15	21	18	16	16	16	-2 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	29212	MHA	15	21	18	19	19	19	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	MHA-F	6	12	9	6	6	6	-3 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49766	MHA-F	19	25	22	23	23	23	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	25	31	28	28	28	28	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3326456	<i>Staphylococcus aureus</i>	29213	MHA	15	21	18	15	15	15	-3 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	29212	MHA	15	21	18	18	19	19	0.67 ± 0.58 (CV=0.87)
	<i>Haemophilus influenzae</i>	49247	MHA-F	6	12	9	6	6	6	-3 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49766	MHA-F	19	25	22	23	23	23	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	25	31	28	29	29	29	1 ± 0 (CV=0)

Clinical Performances Characteristics

A study collected blood samples from HIV positive patients at the Lagos University Teaching Hospital in Nigeria, from October 2005 to September 2006.⁴ The presence of HIV antibodies in the samples was confirmed and isolates were identified through standard biochemical procedures and confirmed using automated methods. Antimicrobial susceptibility testing was performed through the disc diffusion method for a variety of antibiotics, including the Oxoid Ampicillin (10 µg) AMP10 AST Discs. The susceptibility testing was undertaken according to the National Committee for Clinical Laboratory Standards (NCCLS), now the CLSI, guidelines. The isolates were categorised as resistant, intermediate, and susceptible, and the following reference strains were used as quality control: *E. coli* NTCC® 10148 and *Staphylococcus aureus* ATCC® 12600.

The study collected blood samples from 201 HIV patients, of which 26 were positive for a bacterial pathogen. Coagulase negative *Staphylococci* (CoNS) was identified in 15 (58.0%) of cases and non-typhoidal *Salmonellae* identified in six (12.9%) of cases. Other pathogens isolates included *E. coli* and *Pseudomonas aeruginosa*. Susceptibility results found most *Staphylococci* isolates to be resistant to ampicillin (73.3%), high resistance rates were also seen to tetracycline and penicillin G. In contrast, lower resistance was seen among the *Staphylococci* isolates to vancomycin and oxacillin. Similar results were found for the Gram-negative isolates, of which 63.0% were resistant to ampicillin; the same resistance rate was seen for tetracycline, chloramphenicol and co-trimoxazole.

The study found BSIs to be frequent among HIV-infected patients. The study demonstrated that CoNS and Gram-negative rod causing BSIs were frequently susceptible to many antibiotics, however, were often resistant to ampicillin. This finding was in keeping with previous findings. This study did not report any inconsistent or out of specification results for the Oxoid Ampicillin (10 µg) AMP10 AST Discs. Therefore, the study has demonstrated that the Oxoid Ampicillin AST Discs are an effective method for the in vitro analysis of Gram-negative and Gram-positive pathogens causing BSIs and for detecting emerging resistance.

Another study investigated clinical isolates of *Staphylococcus aureus* to determine its nasal colonisation and subsequent antimicrobial susceptibility using the disc diffusion method of various antibiotics.⁵

The study collected paediatric nasal swabs from children aged one to six years in Debre Markos town, Ethiopia, from April to June 2015. The isolates were identified using standard microbiological procedures. Antimicrobial susceptibility testing was performed through disc diffusion for a variety of antibiotics, including the Oxoid Ampicillin (10 µg) AMP10 AST Discs (CT0003B). Susceptibility testing results were interpreted according to the CLSI guidelines and the isolates were categorised as sensitive, intermediate, and resistant. The study used the reference strain, *Staphylococcus aureus* (ATCC-25923), to ensure quality control.

Of the 400 collected nasal swabs, 52 (13%) clinical isolates of *Staphylococcus aureus* were found to be prevalent among the children. Susceptibility testing found just 7.7% of isolates to be susceptible to ampicillin, this was also seen for penicillin. Full susceptibility was seen to oxacillin, cefoxitin, gentamycin and ceftriaxone. Many isolates demonstrated multiple drug resistance to ampicillin and other antibiotics.

The study found that the frequency of *Staphylococcus aureus* in the children was relatively low compared to previous studies, at just 13% of the 400 subjects being prevalent for the bacteria. It was also determined that the clinical isolates were highly resistant to older beta-lactam antibiotics such as penicillin and ampicillin. This study did not report any inconsistent or out of specification results for the Oxoid Ampicillin AST Discs. This study, therefore, demonstrates that the Oxoid Ampicillin AST Discs are an effective method for the in vitro analysis of *S. aureus*.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. A. Adeyemi, A. Sulaiman, D. Solomon, O. Chinedu and I. Victor. (2010). Bacterial bloodstream infections in HIV-infected adults attending a Lagos teaching hospital. *Journal of Health, Population, and Nutrition* 28 (4): 318-326. DOI: 10.3329/JHPN.V28I4.6037.
5. A. Reta, M. Wubie and G. Mekuria. (2017). Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus aureus* Among Pre-school Children in Ethiopia. *BMC Research Notes*, 10 (1) DOI: 10.1186/s13104-017-3079-6.


Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision information

Version	Date of issue and modifications introduced
1.0	2023-01-12.

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. A. Adeyemi, A. Sulaiman, D. Solomon, O. Chinedu and I. Victor. (2010). Bacterial bloodstream infections in HIV-infected adults attending a Lagos teaching hospital. Journal of Health, Population, and Nutrition 28 (4): 318-326. DOI: 10.3329/JHPN.V28I4.6037.
5. A. Reta, M. Wubie and G. Mekuria. (2017). Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus aureus* Among Pre-school Children in Ethiopia. BMC Research Notes, 10 (1) DOI: 10.1186/s13104-017-3079-6.

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Referencie

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).

- a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
- b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. A. Adeyemi, A. Sulaiman, D. Solomon, O. Chinedu and I. Victor. (2010). Bacterial bloodstream infections in HIV-infected adults attending a Lagos teaching hospital. Journal of Health, Population, and Nutrition 28 (4): 318-326. DOI: 10.3329/JHPN.V28I4.6037.
5. A. Reta, M. Wubie and G. Mekuria. (2017). Nasal Colonization and Antimicrobial Susceptibility Pattern of Staphylococcus aureus Among Pre-school Children in Ethiopia. BMC Research Notes, 10 (1) DOI: 10.1186/s13104-017-3079-6.

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© 2023 Thermo Fisher Scientific Inc. Всички права запазени. ATCC® е търговска марка на ATCC. Всички други търговски марки са собственост на Thermo Fisher Scientific Inc. и неговите дъщерни дружества. Тази информация няма за цел да насърчава използването на тези продукти по начин, който може да наруши правата на интелектуална собственост на други.