<https://www.thermofisher.com>

EN

Oxoid™ Ampicillin/Sulbactam Discs 20 µg (SAM20)

REF CT0520B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Ampicillin/Sulbactam Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent ampicillin/sulbactam. The discs are labelled on both sides with details of the antimicrobial (SAM) and amount present (µg): SAM20 (20µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Ampicillin/Sulbactam AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which ampicillin/sulbactam has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- N/A

Gram-Negative

- Enterobacteriaceae
- *Acinetobacter* spp.
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*

Species with published breakpoints according to current CLSI literature:

Gram-Positive

- N/A

Gram-Negative

- Enterobacterales
- *Acinetobacter* spp.
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- N/A

Gram-Negative

- Enterobacterales

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Ampicillin/Sulbactam AST Discs may be used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology, refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Ampicillin/Sulbactam AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Ampicillin/Sulbactam AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Ampicillin/Sulbactam AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT05020B	
Component Description	Material Description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Ampicillin/Sulbactam AST Disc Reagents	
Reagent	Description of Function
Ampicillin/Sulbactam	White crystalline powder. Antimicrobial agent. Ampicillin acts by binding to bacterial penicillin-binding proteins (PBPs) and thus prevent cell wall synthesis. Sulbactam is an irreversible inhibitor of β -lactamase enzyme produced by bacterial cells and therefore both drugs together inhibit bacterial growth.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Ampicillin/Sulbactam AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Ampicillin/Sulbactam AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Results summary for tested QC isolates in accordance with CLSI methodology

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3331953	<i>Staphylococcus aureus</i>	25923™	MHA	29	37	33	34	33	33	0.3 ± 0.6 (CV=1.7)
	<i>Escherichia coli</i>	25922™	MHA	19	24	21.5	20	20	20	-1.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218™	MHA	13	19	16	13	13	13	-3.0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247™	HTM	14	22	18	18	18	18	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3304693	<i>Staphylococcus aureus</i>	25923™	MHA	29	37	33	32	33	34	0 ± 1 (CV=0)
	<i>Escherichia coli</i>	25922™	MHA	19	24	21.5	20	20	20	-1.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218™	MHA	13	19	16	13	13	13	-3.0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247™	HTM	14	22	18	20	20	21	2.3 ± 0.6 (CV=0.3)
	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3278311	<i>Staphylococcus aureus</i>	25923™	MHA	29	37	33	37	37	37	4.0 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922™	MHA	19	24	21.5	19	19	19	-2.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218™	MHA	13	19	16	13	13	13	-3.0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247™	HTM	14	22	18	18	18	18	0 ± 0 (CV=0)

Table 4. Results summary for tested QC isolates in accordance with EUCAST methodology

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3331953	<i>Escherichia coli</i>	25922™	MHA	19	24	21.5	20	20	20	-1.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218™	MHA	13	19	16	13	13	13	-3.0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3304693	<i>Escherichia coli</i>	25922™	MHA	19	24	21.5	20	20	20	-1.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218™	MHA	13	19	16	13	13	13	-3.0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3278311	<i>Escherichia coli</i>	25922™	MHA	19	24	21.5	19	19	19	-2.5 ± 0 (CV=0)
	<i>Escherichia coli</i>	35218™	MHA	13	19	16	13	13	13	-3.0 ± 0 (CV=0)

Clinical Performance Characteristics

A study was conducted to document the occurrence of resistant *Staphylococcus aureus* strains in paediatric patients, as well as the antimicrobial susceptibility profiles when tested against penicillins or against penicillins in combination with other antibiotics.⁴ 564 isolates were collected. Antimicrobial susceptibility testing using 9 antibiotics was performed using the disc diffusion method. The susceptibility profiles of the isolates were analyzed to identify patterns related to gender. (355 isolates of *S. aureus* were collected from males; 209 from females.) 362 of the isolates were methicillin susceptible *S. aureus* (MSSA), 202 were methicillin resistant (MRSA). Within these two groups, the vast majority of *S. aureus* isolates were susceptible to piperacillin-tazobactam, ampicillin-sulbactam and co-amoxiclav. None of the MRSA isolates showed sensitivity to any penicillin when used alone. For MSSA isolates, sensitivity against penicillin alone was significantly lower than if combined with another antibiotic. This study demonstrates that Oxoid Ampicillin/Sulbactam AST Discs were effectively utilised and the most sensitive disc after piperacillin/tazobactam for the clinical analysis of MSSA and MRSA from clinical samples.

In 2018, a study was conducted to investigate the antibiotic susceptibility patterns of *Proteus mirabilis* isolates collected from specimens taken from different body sites of hospital patients.⁵ Antimicrobial susceptibility testing was performed via Kirby-Bauer disc diffusion method using discs containing 17 different antibiotics with results interpreted through CLSI guidelines. Out of 400 isolates, 47 were confirmed positive for *P. mirabilis*. The results reflected a uniform distribution of the bacteria throughout the body. High resistance was seen for 5 of the antibiotics including ampicillin-sulbactam. The resistance against all the other antibiotics tested was classified as intermediate. The patterns of susceptibility and the uniformity of the body distribution within patients corresponds with findings of previous studies. Oxoid Ampicillin/Sulbactam AST Discs proved to be a valuable tool for the analysis of *P. mirabilis*'s susceptibility and distribution over the body.

This study was designed to assess how antibiotic pressure impacts *Escherichia coli* and *Klebsiella* spp. uropathogens isolated from hospital patients.⁶ A total of 148 urine samples were collected from patients with symptomatic urinary tract infections (UTIs). Bacterial isolation and identification were performed through standard microbiological procedures. Antimicrobial susceptibility testing was performed through disc diffusion against several antibiotics, including Oxoid Ampicillin/Sulbactam SAM20 AST discs. Despite the significant presence of Gram-positive isolates, only isolates for *E. coli* and *Klebsiella* spp. were included for testing, being the species of bacteria more commonly responsible for antibiotics resistance in UTI cases. Both *E. coli* and *Klebsiella* spp. showed high resistance against the antibiotics belonging to the penicillin group. The resistance against ampicillin/sulbactam was also high, with 85.7% of *E. coli* isolates and 75% of *Klebsiella* spp. isolates being resistant to ampicillin/sulbactam. The high prevalence of Gram-negative species within the bacterial isolates was consistent with the findings of previous studies, as well as the majority of the strains represented by *E. coli* and *Klebsiella* spp.: this is due to the high pathogenicity of these bacteria, as well as their high resistance that also make them difficult to treat. Oxoid Ampicillin/Sulbactam SAM20 AST Discs proved to be a valuable tool in this study for monitoring antibiotic susceptibility against strains of *Klebsiella* spp. and *E. coli* in an environment where there is excessive use of broad-spectrum antibiotics.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Nosheen, S., H. Ejaz, A. Zafar, and H. Ikram. 2017. Antibacterial Activity of Penicillins Alone and in Combination with Different Agents Against *Staphylococcus Aureus*. *Pakistan Journal of Pharmaceutical Sciences* 30 (2): 393-397.
- Serry F.M., Abdel-Latif H.K., Gomaa S.E., Abbas H.A. 2018. Antimicrobial Resistance of Clinical *Proteus mirabilis* Isolated from Different Sources. *Zagazig Journal Pharmaceutical Science*. June, 2018. ISSN: 1110-5089 / (on-line) 2356-9786
- Innocent, A. A., M. Koffi et al. 2019. Assessment of Drugs Pressure on *Escherichia coli* and *Klebsiella* spp. Uropathogens in Patients Attending Abobo-Avocatier Hospital, North of Abidjan (Côte d'Ivoire) *African Journal of Microbiology Research* 13(29):658-666.

© 2025 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use

	Contains sufficient for <n> tests
	Use-by date
	Period-after-opening – the number of days the product can be used after opening
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

Revision Information

Version	Date of issue and modifications introduced
2.0	2025-04-16. Typographical error corrected.

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Nosheen, S., H. Ejaz, A. Zafar, and H. Ikram. 2017. Antibacterial Activity of Penicillins Alone and in Combination with Different Agents Against *Staphylococcus Aureus*. *Pakistan Journal of Pharmaceutical Sciences* 30 (2): 393-397.
- Serry F.M., Abdel-Latif H.K., Gomaa S.E., Abbas H.A. 2018. Antimicrobial Resistance of Clinical *Proteus mirabilis* Isolated from Different Sources. *Zagazig Journal Pharmaceutical Science*. June, 2018. ISSN: 1110-5089 / (on-line) 2356-9786
- Innocent, A. A., M. Koffi et al. 2019. Assessment of Drugs Pressure on *Escherichia coli* and *Klebsiella* spp. Uropathogens in Patients Attending Aboho-Avocater Hospital, North of Abidjan (Côte d'Ivoire) *African Journal of Microbiology Research* 13(29):658-666.

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Referencije

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Nosheen, S., H. Ejaz, A. Zafar, and H. Ikram. 2017. Antibacterial Activity of Penicillins Alone and in Combination with Different Agents Against *Staphylococcus Aureus*. *Pakistan Journal of Pharmaceutical Sciences* 30 (2): 393-397.
- Serry F.M., Abdel-Latif H.K., Gomaa S.E., Abbas H.A. 2018. Antimicrobial Resistance of Clinical *Proteus mirabilis* Isolated from Different Sources. *Zagazig Journal Pharmaceutical Science*. June, 2018. ISSN: 1110-5089 / (on-line) 2356-9786
- Innocent, A. A., M. Koffi et al. 2019. Assessment of Drugs Pressure on *Escherichia coli* and *Klebsiella* spp. Uropathogens in Patients Attending Abobo-Avocater Hospital, North of Abidjan (Côte d'Ivoire) *African Journal of Microbiology Research* 13(29):658-666.

© 2025. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Litteratur

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Nosheen, S., H. Ejaz, A. Zafar, and H. Ikram. 2017. Antibacterial Activity of Penicillins Alone and in Combination with Different Agents Against *Staphylococcus Aureus*. *Pakistan Journal of Pharmaceutical Sciences* 30 (2): 393-397.
5. Serry F.M., Abdel-Latif H.K., Gomaa S.E., Abbas H.A. 2018. Antimicrobial Resistance of Clinical *Proteus mirabilis* Isolated from Different Sources. *Zagazig Journal Pharmaceutical Science*. June, 2018. ISSN: 1110-5089 / (on-line) 2356-9786
6. Innocent, A. A., M. Koffi et al. 2019. Assessment of Drugs Pressure on *Escherichia coli* and *Klebsiella* spp. Uropathogens in Patients Attending Aboho-Avocater Hospital, North of Abidjan (Côte d'Ivoire) *African Journal of Microbiology Research* 13(29):658-666.

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Viited

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Nosheen, S., H. Ejaz, A. Zafar, and H. Ikram. 2017. Antibacterial Activity of Penicillins Alone and in Combination with Different Agents Against *Staphylococcus Aureus*. *Pakistan Journal of Pharmaceutical Sciences* 30 (2): 393-397.
5. Serry F.M., Abdel-Latif H.K., Gomaa S.E., Abbas H.A. 2018. Antimicrobial Resistance of Clinical *Proteus mirabilis* Isolated from Different Sources. *Zagazig Journal Pharmaceutical Science*. June, 2018. ISSN: 1110-5089 / (on-line) 2356-9786
6. Innocent, A. A., M. Koffi et al. 2019. Assessment of Drugs Pressure on *Escherichia coli* and *Klebsiella* spp. Uropathogens in Patients Attending Abobo-Avocatier Hospital, North of Abidjan (Côte d'Ivoire) *African Journal of Microbiology Research* 13(29):658-666.

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Referências

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Nosheen, S., H. Ejaz, A. Zafar, and H. Ikram. 2017. Antibacterial Activity of Penicillins Alone and in Combination with Different Agents Against *Staphylococcus Aureus*. *Pakistan Journal of Pharmaceutical Sciences* 30 (2): 393-397.
5. Serry F.M., Abdel-Latif H.K., Gomaa S.E., Abbas H.A. 2018. Antimicrobial Resistance of Clinical *Proteus mirabilis* Isolated from Different Sources. *Zagazig Journal Pharmaceutical Science*. June, 2018. ISSN: 1110-5089 / (on-line) 2356-9786
6. Innocent, A. A., M. Koffi et al. 2019. Assessment of Drugs Pressure on *Escherichia coli* and *Klebsiella* spp. Uropathogens in Patients Attending Abovo-Avocater Hospital, North of Abidjan (Côte d'Ivoire) *African Journal of Microbiology Research* 13(29):658-666.

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