	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Cefotaxime Discs 30 µg/5 µg (CTX30/CTX5)

REF CT0166B / CT0407B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Thermo Scientific™ Oxoid™ Cefotaxime Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent cefotaxime. The discs are labelled on both sides with details of antimicrobial (CTX) and amount present (µg): CTX30 (30µg) or CTX5 (5µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack.

Oxoid Cefotaxime AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disk should only be used once.

Intended Use

Oxoid Cefotaxime AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which cefotaxime has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate, or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA, CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA tables:

Gram-Positive

- *Streptococcus* spp. β- Hemolytic Group
- *Streptococcus* spp Viridans Group

Gram-Negative

- *Enterobacteriaceae*
- *Neisseria gonorrhoeae*
- *Neisseria meningitidis*

Species with published breakpoints according to current CLSI M100 tables (all associated with 30 µg cefotaxime):

Gram-Positive

- *Streptococcus* spp. β-haemolytic Group
- *Streptococcus* spp. Viridans Group

Gram-Negative

- *Enterobacterales*
- *Acinetobacter* spp.
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*
- *Neisseria gonorrhoeae*
- *Neisseria meningitidis*

Species with published breakpoints according to current EUCAST tables (all associated with 5 µg cefotaxime):

Gram-Negative

- *Enterobacterales* (indications other than meningitis) and (meningitis)
- *Moraxella catarrhalis*
- *Pasteurella multocida*
- *Kingella kingae*

Gram-Positive

- *Streptococcus* spp. Viridans Group

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principles of Method

Oxoid Cefotaxime AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. A suitable therapeutic agent for in vivo use can be determined using filter paper discs impregnated with specified concentrations of antimicrobial agents placed on the surface of a suitable test medium. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology, refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the Cefotaxime AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Cefotaxime AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Oxoid Cefotaxime AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. The discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, please refer to Table 2.

Table 1. Materials Provided with CT0166B & CT0407B	
Component Description	Material Description
Cartridge	K-Resin plastic cartridge for discs.
Spring/Cap/Plunger	Assembly components (cartridge).
Carton Susceptibility Disc	Absorbent paper individual discs. 6mm.
Desiccant	Molecular sieve tablets 500mg. Light beige to brown small lozenge shaped FBE tablets.
Foil	Hard temper foil. Unprinted.

Table2. Description of Cefotaxime Sodium AST Disc Reagent	
Reagent	Description of Function
Cefotaxime Sodium	Antimicrobial agent. White to slightly yellow-white crystalline powder. The bactericidal activity of cefotaxime results from the inhibition of cell wall synthesis via affinity for penicillin-binding proteins (PBPs). Cefotaxime shows high affinity for penicillin-binding proteins in the cell wall including PBP 1b and PBP III.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Cefotaxime AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Cefotaxime AST Discs have a shelf-life of 24 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2 to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Results Summary for Tested QC Isolates for Oxoid Cefotaxime CTX5 5µg AST Discs in Accordance with EUCAST methodology.

EUCAST Methodology									
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Reading (mm)			Average (mm)+ SD(%)
						1	2	3	
3123682	<i>Escherichia coli</i>	ATCC® 25922	MHA	25	31	28	28	28	28 ±0%
	<i>Klebsiella pneumoniae</i>	ATCC® 700603	MHA	12	18	13	13	13	13 ±0%
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	28	34	34	33	33	33.3 ±0.57%
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	29	37	32	32	32	32 ±0%
3111870	<i>Escherichia coli</i>	ATCC® 25922	MHA	25	31	25	25	25	25 ±0%
	<i>Klebsiella pneumoniae</i>	ATCC® 700603	MHA	12	18	13	13	13	13 ±0%
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	28	34	30	30	30	30 ±0%
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	29	37	33	33	33	33 ±0%
3001689	<i>Escherichia coli</i>	ATCC® 25922	MHA	25	31	26	26	26	26 ±0%
	<i>Klebsiella pneumoniae</i>	ATCC® 700603	MHA	12	18	12	12	12	12 ±0%
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	28	34	29	29	29	29 ±0%
	<i>Haemophilus influenzae</i>	ATCC® 49766	MHA-F	29	37	31	31	31	31 ±0%

Table 4. Results Summary for Tested QC Isolates for Oxoid Cefotaxime CTX30 30µg AST Discs (CT0166B) in Accordance with CLSI Methodology.

CLSI Methodology									
Product Batch / (Date)	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Reading (mm)			Average (mm)+ SD(%)
						1	2	3	
3123682	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	25	31	28	28	28	28 ±0%
	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	35	33	33	33	33 ±0%
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	18	22	20	20	20	20 ±0%
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	31	39	33	33	32	32.6±0.57%
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5%SB	31	39	34	34	34	34 ±0%

	<i>Neisseria gonorrhoea</i>	ATCC® 49226	GCV	38	48	42	43	42	42.3 ±0.57%
	<i>Klebsiella pneumoniae</i>	ATCC® 700603	MHA	17	25	18	18	18	18 ±0%
311870	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	25	31	31	31	31	31±0%
	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	35	33	33	33	33 ±0%
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	18	22	21	21	21	21±0%
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	31	39	31	31	31	31 ±0%
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5%SB	31	39	34	34	34	34 ±0%
	<i>Neisseria gonorrhoea</i>	ATCC® 49226	GCV	38	48	39	39	39	39 ±0%
	<i>Klebsiella pneumoniae</i>	ATCC® 700603	MHA	17	25	19	19	19	19 ±0%
3001689	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	25	31	28	28	28	28 ±0%
	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	35	33	33	33	33 ±0%
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	18	22	20	20	20	20 ±0%
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	31	39	31	31	31	31±0%
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5%SB	31	39	36	36	36	36±0%
	<i>Neisseria gonorrhoea</i>	ATCC® 49226	GCV	38	48	39	39	39	39±0%
	<i>Klebsiella pneumoniae</i>	ATCC® 700603	MHA	17	25	20	20	20	20 ±0%

Clinical Performance Characteristics

A study was conducted to identify the antimicrobial susceptibility pattern of *Haemophilus influenzae* and *Streptococcus pneumoniae* pathogens causing lower respiratory tract infections (LRTIs).⁴ This study collected clinical samples from LRTI patients across the UK from September 2002 to March 2003. The isolates were identified using standard microbiological procedures and antimicrobial susceptibility testing through disc diffusion was performed. The AST was performed using a variety of antibiotics, including the Oxoid Cefotaxime (5 µg) CTX5 AST Discs according to the British Society of Antimicrobial Chemotherapy (BSAC) (guidelines now harmonised with EUCAST). The following reference strains were utilised to ensure quality control: *Haemophilus influenzae* (ATCC-49247), *Escherichia coli* (ATCC-25922), *Streptococcus pneumoniae* (ATCC-49619) and *Staphylococcus aureus* (ATCC-29213).

A total of 1,369 clinical isolates were collected. 100.0% (581) of *H. influenzae* isolates and 95.4% *S. pneumoniae* were susceptible to cefotaxime. Cefotaxime was found to be the most active cephalosporin antibiotic towards the isolates. This study exhibited very high rates of susceptibility of cefotaxime towards *H. influenzae* and *S. pneumoniae*, which indicates that cefotaxime could be a successful therapy for LRTIs in the UK. The study demonstrated that the Oxoid Cefotaxime (5 µg) CTX5 AST Discs can be used successfully for the in vitro clinical analysis of *H. influenzae* and *S. pneumoniae* causing LRTIs. Additionally, the study did not report any inconsistent or out of specification results for the Oxoid Cefotaxime AST Discs, therefore, further confirming clinical performance.

Numerous studies and research papers have used Oxoid Cefotaxime Antimicrobial Susceptibility Test Discs. This reinforces that the device is safe and effective for use described in the intended use.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. I. Morrissey, M. Robbins, L. Viljoen and D. F. J. Brown. Antimicrobial susceptibility of community-acquired respiratory tract pathogens in the UK during 2002/3 determined locally and centrally by BSAC methods. (Morrissey et al. 2005).

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use


	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-13

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- I. Morrissey, M. Robbins, L. Viljoen and D. F. J. Brown. Antimicrobial susceptibility of community-acquired respiratory tract pathogens in the UK during 2002/3 determined locally and centrally by BSAC methods. (Morrissey et al. 2005).

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