<https://www.thermofisher.com>

EN

Oxoid™ Cefoxitin Discs 30 µg (FOX30)

REF CT0119B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Thermo Scientific™ Oxoid™ Cefoxitin (30 µg) FOX30 Antimicrobial Susceptibility Test Discs (CT0119B) are 6mm discs prepared by impregnating high quality absorbent paper with specified concentrations of cefoxitin. The discs are clearly marked on both sides with the corresponding concentration code – FOX30. The code designates the antimicrobial agent (CEFOXITIN) and drug content (30 µg).

The Antimicrobial Susceptibility Test (AST) Discs are supplied in cartridges containing 50 discs each and there are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Oxoid Cefoxitin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disk should only be used once.

Intended Use

Oxoid Cefoxitin AST Discs are used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which cefoxitin has been shown to be active both clinically and *in vitro*. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate, or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Species with published breakpoints according to current FDA literature:

Gram-Negative

- *Neisseria gonorrhoeae*

Species with published breakpoints according to current CLSI literature:

Gram-Negative

- Enterobacterales
- *Neisseria gonorrhoeae*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Staphylococcus spp.* (*Staphylococcus aureus* and *Staphylococcus epidermidis*)

Gram-Negative

- Enterobacterales

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principles of Method

Oxoid Cefoxitin AST Discs are used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Cefoxitin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests.

Materials Provided

Oxoid Cefoxitin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Oxoid™ Cefoxitin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with (CT0119B).	
Component Description	Material Description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Cefoxitin AST Disc Reagents	
Reagent	Description of Function
Cefoxitin sodium salt	Antimicrobial agent. Off-white, light yellow crystalline powder. Mechanism of action: inhibition of cell wall synthesis.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Cefoxitin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Cefoxitin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided (with an unsaturated (orange) desiccant), or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2 - 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table3: Results summary for tested QC isolates, in accordance with EUCAST methodology
EUCAST Methodology

Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3141072	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	23	29	26	25	25	25	1 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	23	29	26	28	28	28	-2 ± 0 (CV=0)
	<i>Staphylococcus aureus (MRSA)</i>	NCTC® 12493	MHA	14	20	17	40	40	40	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3108143	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	23	29	26,	27	27	27	1 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	23	29	26	24	24	24	0 ± 0 (CV=0)
	<i>Staphylococcus aureus (MRSA)</i>	NTCC® 12493*	MHA	14	20	17	18	18	18	3 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3107411	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	23	29	26	29	29	29	3 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	23	29	26	25	25	25	-1 ± 0 (CV=0)
	<i>Staphylococcus aureus (MRSA)</i>	NCTC® 12493	MHA	14	20	17	18	18	18	1 ± 0 (CV=0)

Table4: Results summary for tested QC isolates, in accordance with CLSI methodology

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3141072	<i>Escherichi a coli</i>	ATCC® 25923	MHA	23	29	26	29	29	29	-2 ± 0 (CV=0)
	<i>Staphyloc occus aureus</i>	ATCC® 25922	MHA	23	29	26	24	24	24	3 ± 0 (CV=0)
	<i>Neisseria gonorrhoe ae</i>	ATCC® 49226	GCV	33	41	37	38	38	38	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3108143	<i>Escherichi a coli</i>	ATCC® 25923	MHA	23	29	26	29	29	29	0 ± 0 (CV=0)
	<i>Staphyloc occus aureus</i>	ATCC® 25922	MHA	23	29	26	26	26	26	3 ± 0 (CV=0)
	<i>Neisseria gonorrhoe ae</i>	ATCC® 49226	GCV	33	41	37	41	41	41	4 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3107411	<i>Escherichi a coli</i>	ATCC® 25923	MHA	23	29	26	28	28	28	-1 ± 0 (CV=0)
	<i>Staphyloc occus aureus</i>	ATCC® 25922	MHA	23	29	26	25	25	25	2 ± 0 (CV=0)
	<i>Neisseria gonorrhoe ae</i>	ATCC® 49226	GCV	33	41	37	40	40	40	3 ± 0 (CV=0)

Clinical Performances Characteristics

In a study by Oli et al. 2017, aimed to determine the type and antibiotic susceptibility pattern of bacterial uropathogens isolated from female patients attending Chukwuemeka Odumegwu Ojukwu University Teaching Hospital (COOUTH), Awka, Nigeria⁴. A total of 240 clean catch mid-stream urine samples were collected from females and identified using cultural, gram staining and biochemical characteristics. Out of the 240 samples, 214 had significant bacteriuria with *E. coli* (28.5%) and *Staphylococcus aureus* (28.0%) being the most common isolates in the cases. Antimicrobial susceptibility testing was carried out using the disk diffusion method and interpreted in accordance with CLSI guidelines. Each isolate was tested against a panel of several antibiotic discs including, Thermo Scientific™ Oxoid™ Cefoxitin (30µg) FOX30 Antimicrobial Susceptibility Disks among other disks from Thermo Scientific™ Oxoid™. Results showed a 13.3% susceptibility of *E. coli* to cefoxitin and a mean susceptibility of 22.05% across all organisms (Oli et al. 2017).

Serious Incidents


Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Oli, Angus N., Vivian B. Akabueze, Chijioke E. Ezeudu, George U. Eleje, Obiora S. Ejiofor, Ifeanyichukwu U. Ezebialu, Charlotte B. Oguejiofor, Ifeoma M. Ekejindu, George O. Emechebe, and Kenneth N. Okeke. 2017. "Bacteriology and Antibigram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number
	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period-after-opening – the number of days the product can be used after opening
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged
	Authorized representative in the European Community/ European Union
	Unique device identifier
	European Conformity Mark

	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

Revision information

Version	Date of issue and modifications introduced
2.0	2024-08-20

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- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oli, Angus N., Vivian B. Akabueze, Chijioke E. Ezeudu, George U. Eleje, Obiora S. Ejiofor, Ifeanyichukwu U. Ezebialu, Charlotte B. Oguejiofor, Ifeoma M. Ekejindu, George O. Emechebe, and Kenneth N. Okeke. 2017. "Bacteriology and Antibigram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." The Open Microbiology Journal 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>.

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Reference

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oli, Angus N., Vivian B. Akabueze, Chijioke E. Ezeudu, George U. Eleje, Obiora S. Ejiofor, Ifeanyichukwu U. Ezebialu, Charlotte B. Oguejiofor, Ifeoma M. Ekejindu, George O. Emechebe, and Kenneth N. Okeke. 2017. "Bacteriology and Antibigram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>.

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Litteratur

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oli, Angus N., Vivian B. Akabueze, Chijioke E. Ezeudu, George U. Eleje, Obiora S. Ejiofor, Ifeanyi Chukwu U. Ezebialu, Charlotte B. Oguejiofor, Ifeoma M. Ekejiindu, George O. Emechebe, and Kenneth N. Okeke. 2017. "Bacteriology and Antibigram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>.

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- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oli, Angus N., Vivian B. Akabueze, Chijioke E. Ezeudu, George U. Eleje, Obiora S. Ejiofor, Ifeanyichukwu U. Ezebialu, Charlotte B. Ogueji, Ifeoma M. Ekeji, George O. Emechebe, and Kenneth N. Okeke. 2017. "Bacteriology and Antibiotogram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>.

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Referencias

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oli, Angus N., Vivian B. Akabueze, Chijioke E. Ezeudu, George U. Eleje, Obiora S. Ejiofor, Ifeanyichukwu U. Ezebialu, Charlotte B. Oguejiofor, Ifeoma M. Ekejindu, George O. Emechebe, and Kenneth N. Okeke. 2017. "Bacteriology and AntibioGram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>.

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Referanser

- Clinical and Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Standards for Antimicrobial Susceptibility Testing - M100 - Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method - Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing - Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oli, Angus N., Vivian B. Akabueze, Chijioke E. Ezeudu, George U. Eleje, Obiora S. Ejiofor, Ifeanyichukwu U. Ezebialu, Charlotte B. Ogueji, Ifeoma M. Ekejindu, George O. Emechebe, and Kenneth N. Okeke. 2017. "Bacteriology and AntibioGram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>

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Referencie

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Oli, Angus N., Vivian B. Akabueze, Chijioke E. Ezeudu, George U. Eleje, Obiora S. Ejiofor, Ifeanyichukwu U. Ezebialu, Charlotte B. Oguejiofor, Ifeoma M. Ekejindu, George O. Emechebe, and Kenneth N. Okeke. 2017. "Bacteriology and Antibiogram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>.

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© 2023 Thermo Fisher Scientific Inc. Всички права запазени. ATCC® е търговска марка на ATCC. Всички други търговски марки са собственост на Thermo Fisher Scientific Inc. и неговите дъщерни дружества. Тази информация няма за цел да насърчава използването на тези продукти по начин, който може да наруши правата на интелектуална собственост на други.