	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Ceftazidime Discs 10 µg/30 µg (CAZ10/CAZ30)

REF CT1629B / CT0412B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Thermo Scientific™ Oxoid™ Ceftazidime (10 µg) CAZ10 Antimicrobial Susceptibility Test Discs (AST) (CT1629B) and Ceftazidime (30 µg) CAZ30 Antimicrobial Susceptibility Test Discs (CT0412B) are 6mm paper discs which contain specific amounts of the antimicrobial agent ceftazidime. The discs are labelled on both sides with details of the antimicrobial (CAZ) and amount present (µg): CAZ10 (10µg) or CAZ30 (30µg).

The discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Oxoid ceftazidime AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Oxoid AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which ceftazidime has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Species with published breakpoints according to current FDA literature:

Gram-Negative:

- *Enterobacteriaceae*
- *Pseudomonas aeruginosa*

Species with published breakpoints according to current EUCAST breakpoint tables (all associated with 10 µg ceftazidime):

Gram-Negative:

- *Enterobacterales*
- *Pseudomonas* spp.
- *Aeromonas* spp.
- *Burkholderia pseudomallei*

Species with published breakpoints according to current CLSI M100 (all associated with 30 µg ceftazidime):

Gram-Negative:

- *Enterobacterales*
- *Pseudomonas aeruginosa*
- *Acinetobacter* spp.
- *Burkholderia cepacia*
- *Haemophilus influenzae* and *parainfluenzae*

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principles of Method

Ceftazidime AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Ceftazidime AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI.

Materials Provided

Oxoid Ceftazidime AST Discs consist of 6mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Oxoid Ceftazidime AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, please refer to Table 2.

Table 1. Materials Provided with (CT1629B/CT0412B)	
Component description	Material description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.

Table2. Description of Ceftazidime AST Disc Reagents	
Reagent	Description of Function
Ceftazidime	Ceftazidime is a semisynthetic, broad-spectrum, beta-lactam antibiotic. It is bactericidal in action exerting its effect by inhibition of enzymes responsible for cell-wall synthesis.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Ceftazidime AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life & Storage Conditions

Unopened cartridges of Ceftazidime AST discs have a shelf-life of 12 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided (with an unsaturated (orange) desiccant), or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2 - 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Results summary for tested QC isolates, in accordance with CLSI methodology

CLSI Methodology										
Batch #	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3247844	<i>Staphylococcus aureus</i>	25923	MHA	16	20	18	20	20	20	2 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	25	32	28.5	29	30	30	1.17 ± 0.58 (CV=0.49)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	22	29	25.5	28	28	28	2.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	27	35	31	31	31	31	0 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	49226	GCV	35	43	39	42	42	42	3 ± 0 (CV=0)
Batch #	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3230873	<i>Staphylococcus aureus</i>	25923	MHA	16	20	18	19	19	19	1 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	25	32	28.5	25	25	25	-3.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	22	29	25.5	26	26	26	0.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	27	35	31	29	29	29	-2 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	49226	GCV	35	43	39	42	42	42	3 ± 0 (CV=0)
Batch #	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3207306	<i>Staphylococcus aureus</i>	25923	MHA	16	20	18	20	20	20	2 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	25	32	28.5	28	29	30	0.5 ± 1 (CV=2)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	22	29	25.5	27	27	27	1.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	27	35	31	28	28	28	-3 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	49226	GCV	35	43	39	42	40	41	2 ± 1 (CV=0.5)

Table 4. Results summary for tested QC isolates, in accordance with EUCAST methodology

EUCAST Methodology										
Batch #	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient variance) (%)
							1	2	3	
3247844	<i>Pseudomonas aeruginosa</i>	278853	MHA	21	27	24	24	24	24	0 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	23	29	26	23	23	23	-3 ± 0 (CV=0)
	<i>Klebsiella pneumoniae</i>	700603	MHA	6	12	9	7	8	8	-1.33 ± 0.58 (CV=-0.43)
Batch #	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3214574	<i>Pseudomonas aeruginosa</i>	278853	MHA	21	27	24	25	25	25	1 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	23	29	26	27	27	27	1 ± 0 (CV=0)
	<i>Klebsiella pneumoniae</i>	700603	MHA	6	12	9	10	10	10	1 ± 0 (CV=0)
Batch #	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3154410	<i>Pseudomonas aeruginosa</i>	278853	MHA	21	27	24	23	23	23	2 ± 0 (CV=0)
	<i>Escherichia coli</i>	25922	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Klebsiella pneumoniae</i>	700603	MHA	6	12	9	11	11	11	2 ± 0 (CV=0)

Clinical Performances Characteristics

A study by Ravi S Giriyaapur et al., (2011) compared ceftazidime (30 µg) susceptibility discs to cefotaxime and ceftriaxone discs from Hi-Media Laboratories when detecting the presence of Extended-Spectrum Beta Lactamase (ESBL) producing *Enterobacteriaceae*⁶. A total of 313 consecutive non-repeat culture isolates of *Enterobacteriaceae* were obtained from 280 different specimens from different specialty wards. Among the 313 *Enterobacteriaceae* isolates, *Klebsiella* sp. was the most common pathogen isolated from the tested samples, constituting 115/313 (36.74%) of the total isolates, followed by *Escherichia coli* 82/313 (26.17%), *Providencia* sp. 25/313 (7.98%), *Proteus vulgaris* 24/313 (7.66%) and *Proteus mirabilis* 18/313 (5.75%). Antimicrobial susceptibility testing was determined by the Kirby-Bauer disc diffusion method as per the NCCLS recommendations and screened for ESBLs by NCCLS phenotypic confirmatory combination disc diffusion test. Among the third generation cephalosporins used for screening, ceftazidime discs had a better sensitivity and specificity at 89.49% and 95.74% respectively, followed by ceftriaxone at 86.72% sensitivity and 95.4% specificity and finally cefotaxime at 79.89% sensitivity and 92.30% specificity.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Clinical Laboratory Standards Institute (CLSI).

- a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
- b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
- c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Giriapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing Enterobacteriaceae. Journal of Laboratory Physicians, vol. 3, no. 01, Jan. 2011, pp. 033–036.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Version	Date of issue and modifications introduced
---------	--

1.0	2023-01-13
-----	------------

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.

2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Giriyaapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing Enterobacteriaceae. Journal of Laboratory Physicians, vol. 3, no. 01, Jan. 2011, pp. 033–036.

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Referencije


- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.

2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Giriapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing Enterobacteriaceae. Journal of Laboratory Physicians, vol. 3, no. 01, Jan. 2011, pp. 033–036.

© 2023. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Reference

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Giriapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing Enterobacteriaceae. Journal of Laboratory Physicians, vol. 3, no. 01, Jan. 2011, pp. 033–036.

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The UKCA logo consists of the letters "UK" stacked above "CA" in a bold, sans-serif font.	Značka shody UK
---	-----------------

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Litteratur

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.

3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Giriyaapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing Enterobacteriaceae. Journal of Laboratory Physicians, vol. 3, no. 01, Jan. 2011, pp. 033–036.

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Viitteet

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.

- c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Giriyaapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing Enterobacteriaceae. Journal of Laboratory Physicians, vol. 3, no. 01, Jan. 2011, pp. 033–036.

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Références

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Giriyaapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing Enterobacteriaceae. Journal of Laboratory Physicians, vol. 3, no. 01, Jan. 2011, pp. 033–036.

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- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Giriapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing *Enterobacteriaceae*. Journal of Laboratory Physicians, vol. 3, no. 01, Jan. 2011, pp. 033–036.

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Bibliografia

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Giriapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing *Enterobacteriaceae*. *Journal of Laboratory Physicians*, vol. 3, no. 01, Jan. 2011, pp. 033–036.

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Referências

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.

- b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Giriyaapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing Enterobacteriaceae. Journal of Laboratory Physicians, vol. 3, no. 01, Jan. 2011, pp. 033–036.

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Referințe

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Giriyapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing Enterobacteriaceae. Journal of Laboratory Physicians, vol. 3, no. 01, Jan. 2011, pp. 033–036.

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Referencie

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.

- b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
- c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Giriyaapur, Ravi S, et al. (2011). Comparison of Disc Diffusion Methods for the Detection of Extended-Spectrum Beta Lactamase-Producing Enterobacteriaceae. Journal of Laboratory Physicians, vol. 3, no. 01, Jan. 2011, pp. 033–036.

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