EN

Oxoid™ Cefuroxime Sodium Discs 30 µg (CXM30)

REF CT0127B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Cefuroxime Sodium Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent cefuroxime. The discs are labelled on both sides with details of antimicrobial (CXM) and amount present (30 µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Cefuroxime Sodium AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disk should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which Cefuroxime Sodium has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- N/A

Gram-Negative

- *Enterobacteriaceae*
- *Haemophilus influenzae* and *parainfluenzae*

Species with published breakpoints according to current CLSI literature:

Gram-Positive

- N/A

Gram-Negative

- *Enterobacterales*
- *Haemophilus influenzae* and *parainfluenzae*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Streptococcus* spp. (Viridans group)

Gram-Negative

- *Enterobacterales* — *E. coli*, *Klebsiella* spp. (except *K. aerogenes*), *Raoultella* spp. and *P. mirabilis*
- *Haemophilus influenzae*
- *Moraxella catarrhalis*
- *Kingella kingae*

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Cefuroxime Sodium AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Cefuroxime Sodium AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Cefuroxime Sodium AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Cefuroxime Sodium AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to table 2.

Table 1. Materials Provided with CT0127B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Cefuroxime AST Disc Reagents	
Reagent	Description of Function
Cefuroxime	Off white to light yellow powder. Mechanism of action: Binds to specific penicillin-binding proteins (PBPs) located inside the bacterial cell wall that causes inhibition of the third and last stage of bacterial cell wall synthesis.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Cefuroxime Sodium AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Cefuroxime Sodium AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from 3 batches of Cefuroxime Sodium (30 µg) CXM30 Antimicrobial Susceptibility Test Discs (CT0127B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3108140	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	27	35	31	32	32	32	1 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC 25922	MHA	20	26	23	21	21	21	-2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC 49766	HTM	28	36	32	30	30	30	-2 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC 49226	GCV	33	41	37	37	37	37	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3104594	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	27	35	31	31	31	31	0 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC 25922	MHA	20	26	23	22	22	22	-1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC 49766	HTM	28	36	32	31	31	31	-1 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC 49226	GCV	33	41	37	39	38	38	1.33 ± 0.58 (CV=0.43)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3001685	<i>Staphylococcus aureus</i>	ATCC 25923	MHA	27	35	31	30	30	30	-1 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC 25922	MHA	20	26	23	22	22	22	-1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC 49766	HTM	28	36	32	30	30	30	-2 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC 49226	GCV	33	41	37	35	35	35	-2 ± 0 (CV=0)

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, for 3 batches of Cefuroxime Sodium (30 µg) CXM30 Antimicrobial Susceptibility Test Discs (CT0127B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3158891	<i>Escherichia coli</i>	ATCC 25922	MHA	20	26	23	21	21	21	-2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA-F	28	34	31	31	31	31	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC 49766	MHA-F	26	34	30	30	30	30	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3138686	<i>Escherichia coli</i>	ATCC 25922	MHA	20	26	23	22	22	22	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA-F	28	34	31	33	33	33	2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC 49766	MHA-F	26	34	30	30	30	30	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3108144	<i>Escherichia coli</i>	ATCC 25922	MHA	20	26	23	22	22	22	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA-F	28	34	31	32	32	32	1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC 49766	MHA-F	26	34	30	31	31	31	1 ± 0 (CV=0)

Clinical Performance Characteristics

A total of 145 bacteria strains were isolated from 460 samples collected from urine, wound, throat and anal swabs from 220 healthy volunteers and 240 patients in four hospitals in Enugu metropolis (Nigeria)⁴. For antibiotic susceptibility testing, antibiotic discs used were obtained from Oxoid (England) including cefuroxime (30 µg). Antibiotic sensitivity testing methods were described in the article and inhibition zone diameters (IZDs) were measured and recorded. Isolates were characterised using National Committee for Clinical Laboratory Standard (NCCLS) reference guidelines.

The study demonstrated the presence of ESBL producing bacteria strains from community and hospital subjects within Enugu metropolis (Nigeria). The authors felt their study provided a baseline for future projections and effective management of infectious diseases caused by the ESBL producing bacteria strains. This article shows that the Oxoid Cefuroxime Sodium discs are effective in the clinical analysis of Enterobacterales.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.


References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Afunwa, Ruth A., Odimegwu, Damian, C., Iroha, Romanus I. and Esimone, Charles O. (2011). *Bosnian Journal of Basic Medical Sciences*, Volume 11, Issue 2: 91–96.

© 2025 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number
	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period-after-opening – the number of days the product can be used after opening
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged

	Authorized representative in the European Community/ European Union
	Unique device identifier
	European Conformity Mark
	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

Revision Information

Version	Date of issue and modifications introduced
3.0	2025-02-04 Test media in table 4 updated

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Références

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Afunwa, Ruth A., Odimegwu, Damian, C., Iroha, Romanus I. and Esimone, Charles O. (2011). *Bosnian Journal of Basic Medical Sciences*, Volume 11, Issue 2: 91–96.

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Referanser

1. Clinical and Laboratory Standards Institute (CLSI). Se siste utgave:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Se siste versjon:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Afunwa, Ruth A., Odimegwu, Damian, C., Iroha, Romanus I. and Esimone, Charles O. (2011). *Bosnian Journal of Basic Medical Sciences*, Volume 11, Issue 2: 91–96.

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