<https://www.thermofisher.com>

Oxoid™ Ciprofloxacin Discs 5 µg (CIP5)

REF CT0425B

Antimicrobial Susceptibility Test Discs

EN

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Thermo Scientific™ Oxoid™ Ciprofloxacin (5 µg) CIP5 Antimicrobial Susceptibility Test Discs (CT0425B) are 6 mm paper discs which contain specific amounts of the antimicrobial agent ciprofloxacin. The discs are labelled on both sides with details of the antimicrobial (CIP) and amount present (5 µg).

Discs are supplied in cartridges containing 50 discs each, and there are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Ciprofloxacin Antimicrobial Susceptibility Test Discs (AST) can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Ciprofloxacin Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which ciprofloxacin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram negative

- Enterobacteriaceae
- *Salmonella* spp.
- *Pseudomonas aeruginosa*
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*
- *Neisseria gonorrhoeae*

Gram positive

- *Staphylococcus* spp.
- *Enterococcus* spp.
- *Streptococcus pneumoniae*
- *Streptococcus* spp. β-Hemolytic Group

Species with published breakpoints according to current CLSI M100 literature:

Gram negative

- Enterobacterales
- *Pseudomonas aeruginosa*
- *Acinetobacter* spp.
- *Haemophilus influenzae* and *Haemophilus parainfluenzae*
- *Neisseria gonorrhoeae*
- *Neisseria meningitidis*

Gram positive

- *Staphylococcus* spp.
- *Enterococcus* spp.

Species with published breakpoints according to current EUCAST literature:

Gram negative

- Enterobacterales
- *Pseudomonas* spp.

- *Acinetobacter* spp.
- *Haemophilus influenzae*
- *Moraxella catarrhalis*
- *Pasteurella multocida*
- *Campylobacter jejuni* and *Campylobacter coli*
- *Kingella kingae*
- *Aeromonas* spp.

Gram positive

- *Staphylococcus* spp. (*S. aureus* and *Coagulase-negative staphylococci*)
- *Enterococcus* spp. (uncomplicated UTI only)
- *Corynebacterium* spp.
- *Aerococcus sanguinicola* and *Aerococcus urinae* (uncomplicated UTI only)
- *Bacillus* spp.

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Ciprofloxacin AST Discs may be used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Ciprofloxacin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Ciprofloxacin AST Discs consist of 6mm diameter paper discs impregnated with a specific antimicrobial concentration. The discs are marked on both sides to indicate the agent and amount. Ciprofloxacin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to table 2

Table 1. Materials Provided with (CT0425B).	
Component Description	Material Description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Ciprofloxacin AST Disc Reagent	
Reagent	Description of Function
Ciprofloxacin	White to pale yellow, crystalline powder. Antimicrobial agent with broad-spectrum activity against Gram-positive and Gram-negative bacteria.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Ciprofloxacin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Ciprofloxacin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Clinical Performances Characteristics

A study by Åhman et al., 2019, evaluated the performance and quality of antimicrobial discs from nine manufacturers in 2014 and 2017⁴. Antimicrobial discs of 16 agents, including Thermo Scientific Oxoid Ciprofloxacin AST, from nine manufacturers were evaluated using EUCAST criteria. Discs were tested in triplicate on MHA against EUCAST quality control strains. Mean values were compared with targets and ranges in the EUCAST quality control tables. They concluded that Oxoid Ciprofloxacin AST discs produced mean values within ± 1 mm of the target value in both 2014 and in 2017.

A study by O'Halloran et al., 2017 aimed to determine the comparability of CLSI, EUCAST and Stokes AST methods for determining susceptibility of uropathogenic *E. coli* to various antibiotic discs, including Oxoid ciprofloxacin AST Discs⁵. A total of 100 *E. coli* isolates from boric acid urine samples from patients attending GP surgeries, in Ireland, were collected. A comparison of EUCAST and CLSI without Stokes showed a 94% agreement between methods for Ciprofloxacin.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.


References

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Åhman J, Matuschek E, Kahlmeter G. The quality of antimicrobial discs from nine manufacturers-EUCAST evaluations in 2014 and 2017. *Clin Microbiol Infect.* 2019 Mar;25(3):346–52.
- O'Halloran C, Walsh N, O'Grady MC, Barry L, Hooton C, Corcoran GD, et al. Assessment of the comparability of CLSI, EUCAST and Stokes antimicrobial susceptibility profiles for *Escherichia coli* uropathogenic isolates. *Br J Biomed Sci.* 2018 Jan;75(1):24–9.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests

	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark
	Period-after-opening – the number of days the product can be used after opening
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom


Revision Information

Version	Date of issue and modifications introduced
2.0	2023-11-29

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1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Åhman J, Matuschek E, Kahlmeter G. The quality of antimicrobial discs from nine manufacturers-EUCAST evaluations in 2014 and 2017. *Clin Microbiol Infect.* 2019 Mar;25(3):346–52.
5. O'Halloran C, Walsh N, O'Grady MC, Barry L, Hooton C, Corcoran GD, et al. Assessment of the comparability of CLSI, EUCAST and Stokes antimicrobial susceptibility profiles for *Escherichia coli* uropathogenic isolates. *Br J Biomed Sci.* 2018 Jan;75(1):24–9.

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Viited

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Ähman J, Matuschek E, Kahlmeter G. The quality of antimicrobial discs from nine manufacturers-EUCAST evaluations in 2014 and 2017. *Clin Microbiol Infect*. 2019 Mar;25(3):346–52.
5. O'Halloran C, Walsh N, O'Grady MC, Barry L, Hooton C, Corcoran GD, et al. Assessment of the comparability of CLSI, EUCAST and Stokes antimicrobial susceptibility profiles for *Escherichia coli* uropathogenic isolates. *Br J Biomed Sci*. 2018 Jan;75(1):24–9.

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Viitteet

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ahman J, Matuschek E, Kahlmeter G. The quality of antimicrobial discs from nine manufacturers-EUCAST evaluations in 2014 and 2017. *Clin Microbiol Infect.* 2019 Mar;25(3):346–52.
- O'Halloran C, Walsh N, O'Grady MC, Barry L, Hooton C, Corcoran GD, et al. Assessment of the comparability of CLSI, EUCAST and Stokes antimicrobial susceptibility profiles for Escherichia coli uropathogenic isolates. *Br J Biomed Sci.* 2018 Jan;75(1):24–9.

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Références

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ahman J, Matuschek E, Kahlmeter G. The quality of antimicrobial discs from nine manufacturers-EUCAST evaluations in 2014 and 2017. *Clin Microbiol Infect.* 2019 Mar;25(3):346–52.
- O'Halloran C, Walsh N, O'Grady MC, Barry L, Hooton C, Corcoran GD, et al. Assessment of the comparability of CLSI, EUCAST and Stokes antimicrobial susceptibility profiles for *Escherichia coli* uropathogenic isolates. *Br J Biomed Sci.* 2018 Jan;75(1):24–9.

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Referenzen

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Ahman J, Matuschek E, Kahlmeter G. The quality of antimicrobial discs from nine manufacturers-EUCAST evaluations in 2014 and 2017. *Clin Microbiol Infect.* 2019 Mar;25(3):346–52.
5. O'Halloran C, Walsh N, O'Grady MC, Barry L, Hooton C, Corcoran GD, et al. Assessment of the comparability of CLSI, EUCAST and Stokes antimicrobial susceptibility profiles for *Escherichia coli* uropathogenic isolates. *Br J Biomed Sci.* 2018 Jan;75(1):24–9.

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Nuorodos

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Åhman J, Matuschek E, Kahlmeter G. The quality of antimicrobial discs from nine manufacturers-EUCAST evaluations in 2014 and 2017. *Clin Microbiol Infect.* 2019 Mar;25(3):346–52.
- O'Halloran C, Walsh N, O'Grady MC, Barry L, Hooton C, Corcoran GD, et al. Assessment of the comparability of CLSI, EUCAST and Stokes antimicrobial susceptibility profiles for Escherichia coli uropathogenic isolates. *Br J Biomed Sci.* 2018 Jan;75(1):24–9.

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