
EN

Oxoid™ Clindamycin Discs 2 µg (DA2)

REF CT0064B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Thermo Scientific™ Oxoid™ Clindamycin Antimicrobial Susceptibility Test (AST) discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent Clindamycin (DA). The discs are labelled on both sides with details of antimicrobial (DA) and amount present (µg): DA (2 µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Clindamycin AST discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disk should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which Clindamycin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA tables:

Gram-Negative

- N/A

Gram-Positive

- *Staphylococcus* spp.
- *Streptococcus pneumoniae*
- *Streptococcus* spp. (β- Hemolytic Group)
- *Streptococcus* spp. (Viridans Group)

Species with published breakpoints according to current CLSI M100 tables:

Gram-Negative

- N/A

Gram-Positive

- *Staphylococcus* spp.
- *Streptococcus pneumoniae*
- *Streptococcus* spp. (β- Hemolytic Group)
- *Streptococcus* spp. (Viridans Group)

Species with published breakpoints according to current EUCAST tables:

Gram-Negative

- N/A

Gram-Positive

- *Staphylococcus* spp.
- *Streptococcus pneumoniae*
- *Streptococcus* spp. (Groups A, B, C and G)
- *Streptococcus* spp. (Viridans Group)
- *Corynebacterium* spp.
- *Bacillus* spp.

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Clindamycin AST discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Clindamycin AST discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Oxoid Clindamycin AST discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and concentration. Clindamycin AST discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0064B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Clindamycin AST Disc Reagent	
Reagent	Description of Function
Clindamycin	Antimicrobial agent. White crystalline powder. Reversibly binds to 50S ribosomal subunits preventing peptide bond formation thus inhibiting bacterial protein synthesis; bacteriostatic or bactericidal depending on drug concentration, infection site, and organism.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Clindamycin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Clindamycin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw Data of Zone Sizes in Accordance with CLSI Methodology Taken from 3 Batches of Oxoid Clindamycin (2 µg) DA2 AST Discs.

Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3108141	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	24	30	27	24	24	24	-3 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5 %SB	19	25	22	24	24	24	2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3001686	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	24	30	27	27	27	27	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5 %SB	19	25	22	23	23	23	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
2987565	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	24	30	27	25	25	25	-2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA+5 %SB	19	25	22	23	23	23	1 ± 0 (CV=0)

Table 4: Raw Data of Zone Sizes in Accordance with EUCAST Methodology Taken from 3 Batches of Oxoid Clindamycin (2 µg) DA2 AST Discs.

Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3108141	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	23	29	26	26	26	26	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	22	28	25	23	23	23	-2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3001686	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	23	29	26	27	27	27	1 ± 0 (CV=0)

	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	22	28	25	26	26	26	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
2987565	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	23	29	26	25	25	25	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA-F	22	28	25	23	23	23	-2 ± 0 (CV=0)

Clinical Performance Characteristics

A study was conducted to determine the nasal colonization, associated factors and antimicrobial susceptibility patterns of *S. aureus* isolates.⁴ A total of 400 nasal swab samples were taken and 52 isolates of *S. aureus* identified. Isolates were identified by colony morphology, Gram staining, catalase testing, coagulase testing and latex agglutination assay. Antimicrobial susceptibility testing was conducted using the Kirby-Bauer disk diffusion technique and zone interpretation was carried out in accordance with CLSI guidelines. A control strain of *S. aureus* ATCC® 25923 was used. Each isolate was tested for sensitivity against a panel of antibiotic disks, all supplied by Thermo Scientific Oxoid including Oxoid Clindamycin DA2 (2µg) AST disks. The results displayed no MRSA isolates and 100% susceptibility to clindamycin in all isolates. The results agreed with previous studies referenced.

Literature reviews provide further evidence that Oxoid Clindamycin AST discs are suitable for the assessment of antimicrobial susceptibility in a clinical setting as part of a diagnostic workflow to aid clinicians in treatment options.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References


- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre-School Children in Ethiopia." *BMC Research Notes* 10 (1): 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

© 2021 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number


	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period-after-opening – the number of days the product can be used after opening
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged
	Authorized representative in the European Community/ European Union
	Unique device identifier
	European Conformity Mark
	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

Revision Information

Version	Date of issue and modifications introduced
3.0	2024-10-23 Typographical error corrected

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Reference

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre-School Children in Ethiopia." *BMC Research Notes* 10 (1): 1–7.
<https://doi.org/10.1186/s13104-017-3079-6>

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Referenzen

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of Staphylococcus Aureus among Pre-School Children in Ethiopia." *BMC Research Notes* 10 (1): 1–7.
<https://doi.org/10.1186/s13104-017-3079-6>

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Viitteet

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre-School Children in Ethiopia." *BMC Research Notes* 10 (1): 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

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Nuorodos

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. “Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre-School Children in Ethiopia.” *BMC Research Notes* 10 (1): 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

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Atsauces

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of Staphylococcus Aureus among Pre-School Children in Ethiopia." *BMC Research Notes* 10 (1): 1–7. <https://doi.org/10.1186/s13104-017-3079-6>

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Referanser

- Clinical and Laboratory Standards Institute (CLSI). Se siste utgave:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Se siste versjon:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Reta, Alemayehu, Moges Wubie, and Getnet Mekuria. 2017. "Nasal Colonization and Antimicrobial Susceptibility Pattern of *Staphylococcus Aureus* among Pre-School Children in Ethiopia." *BMC Research Notes* 10 (1): 1–7.
<https://doi.org/10.1186/s13104-017-3079-6>

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