<https://www.thermofisher.com>

EN

Oxoid™ Erythromycin Discs 15 µg (E15)

REF CT0020B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Thermo Scientific™ Oxoid™ Erythromycin Antimicrobial Susceptibility Test Discs (CT0020B) are 6mm discs which contain specific amounts of the antimicrobial agent Erythromycin. The discs are labelled on both sides with details of the antimicrobial (E) and amount present (µg): E15 (15µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Erythromycin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Oxoid Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which erythromycin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA tables, CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- Staphylococcus spp.
- Enterococcus spp.
- Streptococcus pneumoniae
- Streptococcus spp. β-haemolytic group
- Streptococcus spp. Viridans

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- Staphylococcus spp.
- Streptococcus groups A, B, C and G
- Streptococcus pneumoniae
- Listeria monocytogenes
- Bacillus spp.

Gram-Negative

- Moraxella catarrhalis
- Campylobacter jejuni and coli
- Kingella kingae

Species with published breakpoints according to current CLSI literature:

Gram-Positive

- Staphylococcus spp.

- Streptococcus pneumoniae
- Enterococcus spp.
- Streptococcus spp. β -haemolytic Group
- Streptococcus spp. Viridans Group

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principles of Method

Oxoid Erythromycin AST Discs are used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Erythromycin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests.

Materials Provided

Oxoid Erythromycin AST Discs consist of 6mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Erythromycin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, please refer to table 2.

Table 1. Materials Provided with CT0020B	
Component Description	Material Description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table2. Description of Erythromycin AST Disc Reagents	
Critical Ingredient	Description of Function
Erythromycin	Antimicrobial agent. Erythromycin acts by inhibition of protein synthesis.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Erythromycin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life & Storage Conditions

Unopened cartridges of Erythromycin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2 - 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table3: Results summary for tested QC isolates in accordance with CLSI/CLSI VET/EUCAST methodology.

CLSI Methodology									
Product Batch / (Date)	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Reading (mm)			Average (mm)+ SD(%)
						1	2	3	
3108137 (10.08.2020)	<i>Staphylococcus aureus</i>	25923	MHA	22	30	27	27	27	27±0%
	<i>Streptococcus pneumoniae</i>	4961	MHA+5%SB	25	30	29	29	29	29±0%
3001681 (16.07.2020)	<i>Staphylococcus aureus</i>	25923	MHA	22	30	28	28	28	28±0%
	<i>Streptococcus pneumoniae</i>	4961	MHA+5%SB	25	30	29	29	29	29±0%
2987560 (18.06.2020)	<i>Staphylococcus aureus</i>	25923	MHA	22	30	27	27	27	27±0%
	<i>Streptococcus pneumoniae</i>	4961	MHA+5%SB	25	30	29	29	29	29±0%
CLSI VET Susceptibility Methodology									
3108137 (10.08.2020)	<i>Campylobacter jejuni</i>	33560	MHA+5%SB	26	38	31	31	31	31±0%
3001681 (16.07.2020)	<i>Campylobacter jejuni</i>	33560	MHA+5%SB	26	38	30	30	30	30±0%
2987560 (18.06.2020)	<i>Campylobacter jejuni</i>	33560	MHA+5%SB	26	38	28	28	28	28±0%
EUCAST Susceptibility Methodology									
3108137 (10.08.2020)	<i>Staphylococcus aureus</i>	29213	MHA	23	29	26	26	26	26±0%
	<i>Streptococcus pneumoniae</i>	4961	MHA-F	26	32	28	28	28	28±0%
3001681 (16.07.2020)	<i>Staphylococcus aureus</i>	29213	MHA	23	29	25	25	25	25±0%
	<i>Streptococcus pneumoniae</i>	4961	MHA-F	26	32	27	27	27	27±0%
2987560 (18.06.2020)	<i>Staphylococcus aureus</i>	29213	MHA	23	29	25	25	25	25±0%
	<i>Streptococcus pneumoniae</i>	4961	MHA-F	26	32	27	27	28	27.3±0.57%

Clinical Performances Characteristics

A retrospective study analysed various clinical samples, including blood, pus, urine, stool and nasal samples, collected from September 2011 to December 2014 at the Debre Markos Referral Hospital⁴. A total of 575 clinical samples were collected, of which, 100 (35.4%) were identified as *S. aureus* isolates and 16 (5.7%) as *S. pneumoniae*. 66 Gram positive isolates were tested against erythromycin, and 60 (90.9%) were found to be resistant, with 60 out of 62 (96.8) *S. aureus* being resistant. The study did not find any inconsistent or out of specification results for the antimicrobial susceptibility testing and demonstrated

that Oxoid Erythromycin AST Discs (CT0020B) are effective in the clinical analysis of Gram-positive bacteria isolated from a variety of infections.

References


1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. W. Ji, L. Zhang, Z. Guo, S. Xie, W. Yang, J. Chen, J. Wang, Z. Cheng, Xin. Wang, X. Zhu, J. Wang, H. Wang, J. Huang, N. Liang and D. McIver. 2017. Colonization prevalence and antibiotic susceptibility of Group B Streptococcus in pregnant women over a 6-year period in Dongguan, China. PLoS One 12(8). 10.1371/journal.pone.0183083.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number
	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period-after-opening – the number of days the product can be used after opening
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged

	Authorized representative in the European Community/ European Union
	Unique device identifier
	European Conformity Mark
	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

Revision information

Version	Date of issue and modifications introduced
2.0	2024-07-19

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- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- W. Ji, L. Zhang, Z. Guo, S. Xie, W. Yang, J. Chen, J. Wang, Z. Cheng, Xin. Wang, X. Zhu, J. Wang, H. Wang, J. Huang, N. Liang and D. McIver. 2017. Colonization prevalence and antibiotic susceptibility of Group B Streptococcus in pregnant women over a 6-year period in Dongguan, China. PLoS One 12(8). 10.1371/journal.pone.0183083.

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Referanser

- Clinical and Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Standards for Antimicrobial Susceptibility Testing - M100 - Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria. Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method - Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing - Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- W. Ji, L. Zhang, Z. Guo, S. Xie, W. Yang, J. Chen, J. Wang, Z. Cheng, Xin. Wang, X. Zhu, J. Wang, H. Wang, J. Huang, N. Liang and D. McIver. (2017). Colonization prevalence and antibiotic susceptibility of Group B Streptococcus in pregnant women over a 6-year period in Dongguan, China. PLoS One 12(8). 10.1371/journal.pone.0183083.

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Referencie

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- W. Ji, L. Zhang, Z. Guo, S. Xie, W. Yang, J. Chen, J. Wang, Z. Cheng, Xin. Wang, X. Zhu, J. Wang, H. Wang, J. Huang, N. Liang and D. McIver. 2017. Colonization prevalence and antibiotic susceptibility of Group B *Streptococcus* in pregnant women over a 6-year period in Dongguan, China. PLoS One 12(8). 10.1371/journal.pone.0183083.

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© 2023 Thermo Fisher Scientific Inc. Всички права запазени. ATCC® е търговска марка на ATCC. Всички други търговски марки са собственост на Thermo Fisher Scientific Inc. и неговите дъщерни дружества. Тази информация няма за цел да насърчава използването на тези продукти по начин, който може да наруши правата на интелектуална собственост на други.