	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Gentamicin Discs 10 µg/30 µg (CN10/CN30)

REF CT0024B/CT0072B

Antimicrobial Susceptibility Test Discs

Note: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Gentamicin (10 µg/30 µg) CN10 and CN30 Antimicrobial Susceptibility Test Discs (CT0024B and CT0072B) are 6mm discs which contain specific amounts of the antimicrobial agent gentamicin. The discs are labelled on both sides with details of the antimicrobial (CN) and amount present (µg); CN10 (10µg), CN30 (30µg).

Discs are supplied in cartridges containing 50 discs each, and there are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Gentamicin AST Discs can be dispensed using a Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which gentamicin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA tables, CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive:

- *Staphylococcus* spp.

Gram-Negative:

- Enterobacteriaceae
- *Pseudomonas aeruginosa*

Species with published breakpoints according to current CLSI M100:

Gram-Positive:

- *Staphylococcus* spp.

Gram-Negative:

- Enterobacteriaceae
- *Pseudomonas aeruginosa*
- *Acinetobacter* spp.

Species with published breakpoints according to current EUCAST breakpoint tables:

Gram-Positive:

- *Staphylococcus* spp. (*S. aureus*, Coagulase-negative *staphylococci*)

Gram-Negative:

- Enterobacteriaceae
- *Acinetobacter* spp.

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

Principles of Method

Gentamicin AST Discs are used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological traceability of calibrator and control material values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Gentamicin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Gentamicin AST Discs consist of 6mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Gentamicin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, please refer to table 2.

Table 1. Materials Provided with (CT0024B/CT0072B).	
Component description	Material description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge, 5 cartridges per pack.

Table 2. Description of Gentamicin AST Disc Reagents	
Reagent	Description of Function
Gentamicin sulphate	White or almost white hygroscopic powder. Antimicrobial agent. Mechanism of action involves inhibition of bacterial protein synthesis by binding to 30S ribosomes.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Gentamicin AST Discs can be dispensed using a Disc Dispenser, which is not included with the device.

Shelf-life & Storage Conditions

Unopened cartridges of Gentamicin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided (with an unsaturated (orange) desiccant), or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2 - 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Results summary for tested QC isolates Gentamicin (10 µg) CN10 Antimicrobial Susceptibility Test Discs (CT0024B) in accordance with EUCAST methodology

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3141063	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	19	25	22	20	20	20	2 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	19	26	22.5	20	20	20	2.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	17	23	20	18	18	18	2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3141064	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	19	25	22	20	20	20	2 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	19	26	22.5	20	20	20	2.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	17	23	20	19	19	19	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3138679	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	19	25	22	20	20	20	2 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	19	26	22.5	19	19	19	3.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	17	23	20	18	18	18	2 ± 0 (CV=0)

Table 4. Results summary for tested QC isolates Gentamicin (10 µg) CN10 Antimicrobial Susceptibility Test Discs (CT0024B), in accordance with CLSI methodology

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3141063	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	19	27	18	25	25	25	7 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	19	26	22.5	20	20	20	2.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	17	23	20	18	18	18	2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3141064	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	19	27	18	26	26	26	8 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	19	26	22.5	20	20	20	2.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	17	23	20	19	19	19	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3138679	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	19	27	18	24	24	24	6 ± 0 (CV=0)
	<i>Escherichia coli</i>	ATCC® 25922	MHA	19	26	22.5	19	19	19	3.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	17	23	20	18	18	18	2 ± 0 (CV=0)

Table 5. Results summary for tested QC isolates Gentamicin (30 µg) CN30 AST Discs (CT0024B), in accordance with EUCAST methodology

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3149024	<i>Enterococcus faecalis</i>	29212	MHA	12	18	15	16	16	16	1 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	51299	MHA	6	6	6	6	6	6	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
2994167	<i>Enterococcus faecalis</i>	29212	MHA	12	18	15	16	16	16	1 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	51299	MHA	6	6	6	6	6	6	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
2956088	<i>Enterococcus faecalis</i>	29212	MHA	12	18	15	15	15	15	0 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	51299	MHA	6	6	6	6	6	6	0 ± 0 (CV=0)

Clinical Performances Characteristics

Clinical performance of Gentamicin AST discs has been reported in numerous sources, all supporting the routine clinical use to determine antimicrobial susceptibility.

Azene and Bayeh⁴ investigated the antimicrobial susceptibility profile of pathogens cultured from wound infections. The study looked at swabs collected from 2003 to 2010 with a total of 559 samples being analysed of which 422 were culture positive. A total of 500 bacteria were isolated and identified using colony morphology, Gram-stain and biochemical tests as part of standard methods. The most frequently isolated pathogen was *S. aureus* (41.6%), followed by *Pseudomonas spp.* (18.4%), *E. coli* (16.4%), *Proteus spp.* (11.0%), *Enterobacter spp.* and *Citrobacter spp.* (4.2%), *Klebsiella spp.* (2.4%) and Coagulase-negative *Staphylococci* (CoNS) (1.8%). Antimicrobial susceptibility testing was performed using the Kirby-Bauer disk diffusion method and was interpreted in accordance with CLSI guidelines. Each isolate was tested for susceptibility against a panel of antibiotic discs (from Oxoid) including Gentamicin (10µg) CN10 Antimicrobial Susceptibility Test Discs. Results showed that, of the *S. aureus* isolates, resistance to gentamicin was 12.4%. For *Pseudomonas spp.*, resistance to gentamicin was 7.8%, in *E. coli* isolates, resistance was 14.4%, in *Proteus spp.* isolates, resistance was 9.6%. Furthermore, in *Citrobacter spp.*, *Enterobacter spp.* and CoNS isolates resistance was 16.7%, 4.8% and 14.3% respectively (Figure 7.). This article demonstrated that gentamicin discs were effective in the clinical analysis of Enterobacteriaceae from a clinical setting.

Haque and Abdus⁵ conducted a retrospective analysis of 443 consecutively collected midstream and/or catheter-catch urine samples from patients with a clinical suspicion of urinary tract infection. Of these, 189 produced a positive culture with the most predominant isolate being *E. coli* (59.3%) followed by *Staphylococcus saprophyticus* (19%), *Enterococcus spp.* (11.6%), *Klebsiella spp.* (5.5%), *Pseudomonas spp.* (2%), *Proteus spp.* (1.5%) and *Enterobacter spp.* (1%). Antimicrobial susceptibility testing was performed using the Kirby-Bauer disc diffusion method and was interpreted in accordance with CLSI guidelines. Isolates were tested for sensitivity against a panel of 11 antibiotics (from Oxoid) including Gentamicin (10µg) CN10 Antimicrobial Susceptibility Test Discs. Results showed resistance to gentamicin was 40.7% in *E. coli* isolates, 47.4% in *S. saprophyticus* isolates, 56.5% in *Enterococcus spp.*, 27.3% in *Klebsiella spp.*, 75% in *Pseudomonas spp.*, 33.3% in *Proteus spp.*, and 50% in *Enterobacter spp.* The study demonstrated that gentamicin discs were an effective tool in the determination and surveillance of antibiotic resistance in a clinical setting.

A study by Kot et al⁶ aimed to determine resistance profiles of methicillin-resistant *Staphylococcus aureus* (MRSA) clinical isolates from 2015 to 2017. A total of 112 isolates were collected and identifications confirmed with PCR. AST was carried out using the disc diffusion method in accordance with the Clinical Laboratory Standards Institute (CLSI) guidelines against gentamicin discs (10µl), as well as other antimicrobials. Results showed resistance to gentamicin was 8.0% in the *S. aureus* clinical isolates. The study as a whole concluded that isolates from 2017 were resistant to a higher number of antimicrobials when compared to 2016 isolates, indicating the rapid acquisition of resistance within the hospital environment. Results generated using gentamicin discs showed no inconsistencies in the context of the study.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established. In the event of malfunction do not use device.

References

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Azene, Mulugeta K., and Bayeh A. Beyene. 2011. "Bacteriology and Antibigram of Pathogens from Wound Infections at Dessie Laboratory, North-East Ethiopia." *Tanzania Journal of Health Research* 13 (4): 68–74.
<https://doi.org/10.4314/thrb.v13i4.64901>.
5. Haque, Rezwana, Most Laila Akter, and Md Abdus Salam. 2015. "Prevalence and Susceptibility of Uropathogens: A Recent Report from a Teaching Hospital in Bangladesh." *BMC Research Notes* 8 (September): 416.
<https://doi.org/10.1186/s13104-015-1408-1>.
6. Kot, Barbara, Kamila Wierchowska, Małgorzata Piechota, and Agata Gruzewska. 2020. "Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalised during 2015-2017 in Hospitals in Poland." *Medical Principles and Practice* 29 (1): 61–68.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of Symbols

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision information

Version	Date of issue and modifications introduced
1.0	2023-01-16

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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Referencije

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Azene, Mulugeta K., and Bayeh A. Beyene. 2011. "Bacteriology and Antibigram of Pathogens from Wound Infections at Dessie Laboratory, North-East Ethiopia." *Tanzania Journal of Health Research* 13 (4): 68–74.
<https://doi.org/10.4314/thrb.v13i4.64901>.
5. Haque, Rezwana, Most Laila Akter, and Md Abdus Salam. 2015. "Prevalence and Susceptibility of Uropathogens: A Recent Report from a Teaching Hospital in Bangladesh." *BMC Research Notes* 8 (September): 416.
<https://doi.org/10.1186/s13104-015-1408-1>.
6. Kot, Barbara, Kamila Wierzychowska, Małgorzata Piechota, and Agata Gruzewska. 2020. "Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalised during 2015-2017 in Hospitals in Poland." *Medical Principles and Practice* 29 (1): 61–68.

© 2023. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Reference

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Azene, Mulugeta K., and Bayeh A. Beyene. 2011. "Bacteriology and Antibigram of Pathogens from Wound Infections at Dessie Laboratory, North-East Ethiopia." *Tanzania Journal of Health Research* 13 (4): 68–74. <https://doi.org/10.4314/thrb.v13i4.64901>.
5. Haque, Rezwana, Most Laila Akter, and Md Abdus Salam. 2015. "Prevalence and Susceptibility of Uropathogens: A Recent Report from a Teaching Hospital in Bangladesh." *BMC Research Notes* 8 (September): 416. <https://doi.org/10.1186/s13104-015-1408-1>.
6. Kot, Barbara, Kamila Wierzchowska, Małgorzata Piechota, and Agata Gruzewska. 2020. "Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalised during 2015-2017 in Hospitals in Poland." *Medical Principles and Practice* 29 (1): 61–68.

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Litteratur

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Azene, Mulugeta K., and Bayeh A. Beyene. 2011. "Bacteriology and Antibigram of Pathogens from Wound Infections at Dessie Laboratory, North-East Ethiopia." *Tanzania Journal of Health Research* 13 (4): 68–74.
<https://doi.org/10.4314/thrb.v13i4.64901>.
5. Haque, Rezwana, Most Laila Akter, and Md Abdus Salam. 2015. "Prevalence and Susceptibility of Uropathogens: A Recent Report from a Teaching Hospital in Bangladesh." *BMC Research Notes* 8 (September): 416.
<https://doi.org/10.1186/s13104-015-1408-1>.
6. Kot, Barbara, Kamila Wierzbowska, Małgorzata Piechota, and Agata Gruzewska. 2020. "Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalised during 2015-2017 in Hospitals in Poland." *Medical Principles and Practice* 29 (1): 61–68.

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Viited

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Azene, Mulugeta K., and Bayeh A. Beyene. 2011. "Bacteriology and Antibigram of Pathogens from Wound Infections at Dessie Laboratory, North-East Ethiopia." Tanzania Journal of Health Research 13 (4): 68–74. <https://doi.org/10.4314/thrb.v13i4.64901>.
- Haque, Rezwana, Most Laila Akter, and Md Abdus Salam. 2015. "Prevalence and Susceptibility of Uropathogens: A Recent Report from a Teaching Hospital in Bangladesh." BMC Research Notes 8 (September): 416. <https://doi.org/10.1186/s13104-015-1408-1>.
- Kot, Barbara, Kamila Wierzchowska, Małgorzata Piechota, and Agata Gruzewska. 2020. "Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalised during 2015-2017 in Hospitals in Poland." Medical Principles and Practice 29 (1): 61–68.

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Viitteet

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Azene, Mulugeta K., and Bayeh A. Beyene. 2011. "Bacteriology and Antibigram of Pathogens from Wound Infections at Dessie Laboratory, North-East Ethiopia." *Tanzania Journal of Health Research* 13 (4): 68–74.
<https://doi.org/10.4314/thrb.v13i4.64901>.
5. Haque, Rezwana, Most Laila Akter, and Md Abdus Salam. 2015. "Prevalence and Susceptibility of Uropathogens: A Recent Report from a Teaching Hospital in Bangladesh." *BMC Research Notes* 8 (September): 416.
<https://doi.org/10.1186/s13104-015-1408-1>.
6. Kot, Barbara, Kamila Wierzchowska, Małgorzata Piechota, and Agata Gruzewska. 2020. "Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalised during 2015–2017 in Hospitals in Poland." *Medical Principles and Practice* 29 (1): 61–68.

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Références

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Azene, Mulugeta K., and Bayeh A. Beyene. 2011. "Bacteriology and Antibigram of Pathogens from Wound Infections at Dessie Laboratory, North-East Ethiopia." *Tanzania Journal of Health Research* 13 (4): 68–74. <https://doi.org/10.4314/thrb.v13i4.64901>.
5. Haque, Rezwana, Most Laila Akter, and Md Abdus Salam. 2015. "Prevalence and Susceptibility of Uropathogens: A Recent Report from a Teaching Hospital in Bangladesh." *BMC Research Notes* 8 (September): 416. <https://doi.org/10.1186/s13104-015-1408-1>.
6. Kot, Barbara, Kamila Wierchowska, Małgorzata Piechota, and Agata Gruzewska. 2020. "Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalised during 2015-2017 in Hospitals in Poland." *Medical Principles and Practice* 29 (1): 61–68.

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Referenzen

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Azene, Mulugeta K., and Bayeh A. Beyene. 2011. "Bacteriology and Antibigram of Pathogens from Wound Infections at Dessie Laboratory, North-East Ethiopia." *Tanzania Journal of Health Research* 13 (4): 68–74.
<https://doi.org/10.4314/thrb.v13i4.64901>.
5. Haque, Rezwana, Most Laila Akter, and Md Abdus Salam. 2015. "Prevalence and Susceptibility of Uropathogens: A Recent Report from a Teaching Hospital in Bangladesh." *BMC Research Notes* 8 (September): 416.
<https://doi.org/10.1186/s13104-015-1408-1>.
6. Kot, Barbara, Kamila Wierzchowska, Małgorzata Piechota, and Agata Gruzewska. 2020. "Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalised during 2015-2017 in Hospitals in Poland." *Medical Principles and Practice* 29 (1): 61–68.

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1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Azene, Mulugeta K., and Bayeh A. Beyene. 2011. "Bacteriology and Antibigram of Pathogens from Wound Infections at Dessie Laboratory, North-East Ethiopia." Tanzania Journal of Health Research 13 (4): 68–74.
<https://doi.org/10.4314/thrb.v13i4.64901>.
5. Haque, Rezwana, Most Laila Akter, and Md Abdus Salam. 2015. "Prevalence and Susceptibility of Uropathogens: A Recent Report from a Teaching Hospital in Bangladesh." BMC Research Notes 8 (September): 416.
<https://doi.org/10.1186/s13104-015-1408-1>.
6. Kot, Barbara, Kamila Wierzchowska, Małgorzata Piechota, and Agata Gruzewska. 2020. "Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalised during 2015-2017 in Hospitals in Poland." Medical Principles and Practice 29 (1): 61–68.

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Referințe

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Azene, Mulugeta K., and Bayeh A. Beyene. 2011. "Bacteriology and Antibigram of Pathogens from Wound Infections at Dessie Laboratory, North-East Ethiopia." *Tanzania Journal of Health Research* 13 (4): 68–74. <https://doi.org/10.4314/thrb.v13i4.64901>.
- Haque, Rezwana, Most Laila Akter, and Md Abdus Salam. 2015. "Prevalence and Susceptibility of Uropathogens: A Recent Report from a Teaching Hospital in Bangladesh." *BMC Research Notes* 8 (September): 416. <https://doi.org/10.1186/s13104-015-1408-1>.
- Kot, Barbara, Kamila Wierchowska, Małgorzata Piechota, and Agata Gruzewska. 2020. "Antimicrobial Resistance Patterns in Methicillin-Resistant *Staphylococcus aureus* from Patients Hospitalised during 2015-2017 in Hospitals in Poland." *Medical Principles and Practice* 29 (1): 61–68.

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