EN

Oxoid™ Levofloxacin Discs 5 µg (LEV5)

REF CT1587B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Levofloxacin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent levofloxacin. The discs are labelled on both sides with details of the antimicrobial (LEV) and amount present (µg): LEV5 (5µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Levofloxacin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which levofloxacin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Enterococcus* spp.
- *Streptococcus pneumoniae*
- *Streptococcus* spp. β-Hemolytic Group
- Methicillin-susceptible *Staphylococcus aureus*

Gram-Negative

- *Enterobacteriaceae*
- *Pseudomonas aeruginosa*
- *Haemophilus influenzae* and *parainfluenza*

Species with published breakpoints according to current CLSI M100 literature:

Gram-Positive

- *Staphylococcus* spp.
- *Enterococcus* spp.
- *Streptococcus pneumoniae*
- *Streptococcus* spp. β-Hemolytic Group
- *Streptococcus* spp. *Viridans* Group

Gram-Negative

- *Enterobacterales*
- *Pseudomonas aeruginosa*
- *Acinetobacter* spp.
- *Stenotrophomonas maltophilia*
- *Haemophilus influenzae* and *parainfluenzae*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Staphylococcus* spp. (*S. aureus*, coagulase-negative *Staphylococcus* spp.)
- *Enterococcus* spp. (uncomplicated UTI only)
- *Streptococcus pneumoniae*
- *Streptococcus* groups A, B, C and G
- *Bacillus* spp.

Gram-Negative

- *Enterobacterales*
- *Pseudomonas* spp.
- *Acinetobacter* spp.
- *Haemophilus influenzae*
- *Moraxella catarrhalis*
- *Pasteurella multocida*
- *Kingella kingae*
- *Aeromonas* spp.

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Levofloxacin AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Levofloxacin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Levofloxacin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Levofloxacin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT1587B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table2. Description of Levofloxacin AST Disc Reagents	
Reagent	Description of Function
Levofloxacin	Crystalline solid; soluble in PBS, ethanol, DMSO and DMF. Levofloxacin is used to treat a wide variety of bacterial infections. This fluoroquinolone antibacterial agent acts on the DNA-DNA-gyrase complex and topoisomerase IV.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Levofloxacin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Levofloxacin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2 °to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 1. Raw Data of Zone Sizes in Accordance with CLSI Methodology Taken from the 3 Most Recent Batches of Levofloxacin (5 µg) LEV5 AST Discs.

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3274023	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	37	33	31	31	31	-2 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	25	30	27.5	25	25	25	-2.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA + 5% SB	19	26	22.5	22	22	22	-0.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	32	40	36	35	35	35	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA + 5% SB	20	25	22.5	25	25	25	2.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
325	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	37	33	37	37	36	3.67 ± 0.58 (CV=0.16)

	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	25	30	27.5	29	29	29	1.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA + 5% SB	19	26	22.5	24	24	24	1.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	32	40	36	38	38	38	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA + 5% SB	20	25	22.5	22	22	22	-0.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3235885	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	37	33	30	30	30	-3 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 25923	MHA	25	30	27.5	25	25	25	-2.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA + 5% SB	19	26	22.5	22	22	22	-0.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	HTM	32	40	36	39	40	40	3.67 ± 0.58 (CV=0.16)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619	MHA + 5% SB	20	25	22.5	22	22	22	-0.5 ± 0 (CV=0)

SD: Standard Deviation

Table 2. Raw Data of Zone Sizes in Accordance with EUCAST Methodology Taken from the 3 Most Recent Batches of Levofloxacin (5 µg) LEV5 AST Discs.

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3274023	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	37	33	31	31	31	-2 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	23	29	26	25	25	25	-1 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	19	26	22.5	22	22	22	-0.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3258274	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	37	33	37	37	37	3.67 ± 0.58 (CV=0.15)

	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	23	29	26	27	27	27	1 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	19	26	22.5	24	24	24	1.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3235885	<i>Escherichia coli</i>	ATCC® 25922	MHA	29	37	33	30	30	30	-3 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	ATCC® 29213	MHA	23	29	26	24	24	24	-2 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	ATCC® 27853	MHA	19	26	22.5	22	22	22	-0.5 ± 0 (CV=0)

SD: Standard Deviation

Clinical Performance Characteristics

144 clinical isolates were included in AST susceptibility testing of Oxoid Levofloxacin discs via disk diffusion according to the EUCAST breakpoints guidelines⁴. Results showed 100% susceptibility to levofloxacin for the *S. pneumoniae*, *H. influenzae* and *M. catarrhalis* isolates.

A study shows that bacterial isolates were resistant to commonly prescribed drugs for urinary tract infections (UTIs), and the study suggested that levofloxacin is a promising drug in the management of urinary tract infections in the study centre.⁵ The article shows the Levofloxacin (5 µg) LEV5 AST Discs are effective with clinical analysis of *Salmonella* spp., *P. aeruginosa*, *Escherichia coli* and *S. aureus*, given the susceptibility rates and no out of specification or inconsistent results were reported.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005). Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. *PLoS ONE* 11(7): e0159183.
- Oli *et al.* (2017). "Bacteriology and Antiblogram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols


Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device

	Temperature limitation (storage temp.)
LOT	Lot / Batch Number
REF	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period-after-opening – the number of days the product can be used after opening
R_x only	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged
EC REP	Authorized representative in the European Community/ European Union
UDI	Unique device identifier
CE	European Conformity Mark
UK CA	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

Revision Information

Version	Date of issue and modifications introduced
2.0	2024-07-23

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005). Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. *PLoS ONE* 11(7): e0159183.
- Oli *et al.* (2017). "Bacteriology and AntibioGram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>

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Reference

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005). Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. *PLoS ONE* 11(7): e0159183.
- Oli *et al.* (2017). "Bacteriology and Antibigram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>

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Referenzen

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005). Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. *PLoS ONE* 11(7): e0159183.
- Oli *et al* (2017). "Bacteriology and Antiblogram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005). Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. *PLoS ONE* 11(7): e0159183.
- Oli *et al* (2017). "Bacteriology and Antibigram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>

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Nuorodos

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005). Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. *PLoS ONE* 11(7): e0159183.
- Oli *et al.* (2017). “Bacteriology and Antibigram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria.” *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>

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Referanser

- Clinical and Laboratory Standards Institute (CLSI). Se siste utgave:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Se siste versjon:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005). Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. *PLoS ONE* 11(7): e0159183.
- Oli *et al* (2017). "Bacteriology and AntibioGram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005). Idelevich et al. (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
- Oli et al (2017). "Bacteriology and Antibigram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." The Open Microbiology Journal 11 (1): 292–300. <https://doi.org/10.2174/1874285801711010292>

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Referencie

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. *PLoS ONE* 11(7): e0159183.
4. Oli *et al.* (2017). "Bacteriology and Antibigram of Urinary Tract Infection Among Female Patients in a Tertiary Health Facility in South Eastern Nigeria." *The Open Microbiology Journal* 11 (1): 292–300.
<https://doi.org/10.2174/1874285801711010292>

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