
	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Linezolid Discs 10μ (LZD10)

REF CT1649B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Linezolid Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent linezolid. The discs are labelled on both sides with details of the antimicrobial (VA) and amount present (μg): LZD10 (10μg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Linezolid AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Linezolid Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which linezolid has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Staphylococcus* spp.
- *Enterococcus* spp.
- *Streptococcus pneumoniae*
- *Streptococcus* spp. B-Hemolytic group
- *Streptococcus* spp. Viridans Group

Gram-Negative

- N/A

Species with published breakpoints according to current CLSI M100:

Gram-Positive

- N/A

Gram-Negative

- N/A

Species with published breakpoints according to current EUCAST breakpoint tables:

Gram-Positive

- *Staphylococcus* spp.
- *Enterococcus* spp.
- *Streptococcus* groups A, B, C, & G
- *Streptococcus pneumoniae*
- *Corynebacterium* spp.
- *Bacillus* spp.

Gram-Negative

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Linezolid AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology, refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Linezolid AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Linezolid AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Linezolid AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT1649B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Linezolid AST Disc Reagent	
Reagent	Description of Function
Linezolid	White or almost white powder. Antimicrobial agent. The drug mechanism of action involves the binding to 23S ribosomal RNA preventing the formation of the 70S initiation complex required for bacterial reproduction.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Linezolid AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Linezolid AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 1. Raw Data of Zone Sizes in Accordance with CLSI Methodology Taken from the Three Most Recent Batches of Linezolid AST Discs.

CLSI Methodology										
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance)
							1	2	3	
3266019	<i>Staphylococcus aureus</i>	ATCC 29213	MHA	21	27	24	23	23	23	-1 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	ATCC 29212	MHA	19	25	22	21	21	21	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA + 5% HB	23	29	26	27	27	27	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3164386	<i>Staphylococcus aureus</i>	ATCC 29213	MHA	21	27	24	21	21	21	-3 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	ATCC 29212	MHA	19	25	22	20	20	20	-2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA + 5% HB	23	29	26	27	27	27	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Average (mm) + SD
							1	2	3	
3207312	<i>Staphylococcus aureus</i>	ATCC 29213	MHA	21	27	24	24	24	24	0 ± 0 (CV=0)
	<i>Enterococcus faecalis</i>	ATCC 29212	MHA	19	25	22	20	20	20	-2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA + 5% HB	23	29	26	25	25	25	-1 ± 0 (CV=0)

Clinical Performance Characteristics

A study by Ferreira et al., 2012 aimed to evaluate the antimicrobial susceptibility in *Staphylococci* species isolated from UTI confirmed urine samples (a CFU count $\geq 100,000$ CFUs/mL) obtained from hospitals and health centres in Botucatu between 10th March and 14th November 2008. The susceptibility of these *Staphylococci* species was determined by Oxoid Linezolid (30 µg) LZD30 Antimicrobial Susceptibility Test Discs (CT1649B) according to CLSI criteria. Of the 101 *Staphylococci* spp. isolates tested all strains were susceptible to linezolid; this is in 100% accordance to the results obtained by E test and Vitek I. The study did not report any out of specification or inconsistent results associated with the Oxoid Linezolid AST Discs demonstrating that the discs are effective in the clinical analysis of *Staphylococci* spp. species, given no inconsistencies in results and the agreement with alternative susceptibility testing.

Hafeez et al. (2016) isolated 526 clinical samples from intensive care unit (ICU) patients (between January 2013 and January 2014) in Rawalpindi, Pakistan. Isolates were identified using colony morphology, Gram stain, catalase test and oxidase test. The most common Gram-positive pathogens isolated were (1) Coagulase Negative *Staphylococcus* (2) *Enterococcus* spp. (3) *Enterobacter* spp. (4) *Staphylococcus aureus* (5) *Streptococcus* spp. All isolated gram positive cocci were 100% susceptible to linezolid when tested with Oxoid Linezolid AST Discs (however, the concentration of linezolid present on the disc was not specified) when interpreted according to CLSI guidelines. The study did not report any out-of-specification or inconsistent results associated with the Oxoid Linezolid AST Discs, therefore these discs are effective in the clinical analysis of multiple Gram positive

pathogens ((1) Coagulase Negative Staphylococcus (2) *Enterococcus* spp. (3) *Enterobacter* spp. (4) *Staphylococcus aureus* (5) *Streptococcus* spp).

A study by Iorio et al., 2012 investigated the characteristics of 35 methicillin-resistant *S. epidermidis* (MRSE) isolates (14 blood, 7 nares, 4 surgical site were the predominant sites) in Rio de Janeiro hospitals confirmed by phenotype and PCR. Methicillin-resistance was analysed by AST discs included Oxoid Linezolid AST Discs (however, the concentration of linezolid present on the disc was not specified) in accordance with CLSI guidelines. The results demonstrated that all MRSE strains were 100% susceptible to linezolid (irrespective of the resistance profile) therefore, demonstrating Oxoid Linezolid AST Discs are effective in the clinical analysis of *S. epidermidis* given no out-of-specification or inconsistent results were reported in the study.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References


1. Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Ferreira, Adriano M., Mariana F. Bonesso, Alessandro L. Mondelli, Carlos H. Camargo, and Maria De Lourdes R.S. Cunha. 2012. "Oxacillin Resistance and Antimicrobial Susceptibility Profile of *Staphylococcus saprophyticus* and Other *Staphylococci* Isolated from Patients with Urinary Tract Infection." *Chemotherapy* 58: 482–91. <https://doi.org/10.1159/000346529>.
5. Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
6. Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Gonçalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant *Staphylococcus Epidermidis* ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests

	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-16

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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1. Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Ferreira, Adriano M., Mariana F. Bonesso, Alessandro L. Mondelli, Carlos H. Camargo, and Maria De Lourdes R.S. Cunha. 2012. "Oxacillin Resistance and Antimicrobial Susceptibility Profile of *Staphylococcus saprophyticus* and Other *Staphylococci* Isolated from Patients with Urinary Tract Infection." *Chemotherapy* 58: 482–91. <https://doi.org/10.1159/000346529>.
5. Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
6. Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Goncalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant *Staphylococcus Epidermidis* ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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Referencije

1. Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Ferreira, Adriano M., Mariana F. Bonesso, Alessandro L. Mondelli, Carlos H. Camargo, and Maria De Lourdes R.S. Cunha. 2012. "Oxacillin Resistance and Antimicrobial Susceptibility Profile of Staphylococcus epidermidis and Other Staphylococci Isolated from Patients with Urinary Tract Infection." *Chemotherapy* 58: 482–91. <https://doi.org/10.1159/000346529>.
5. Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
6. Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Goncalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant Staphylococcus Epidermidis ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

© 2023 Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Reference

1. Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Ferreira, Adriano M., Mariana F. Bonesso, Alessandro L. Mondelli, Carlos H. Camargo, and Maria De Lourdes R.S. Cunha. 2012. "Oxacillin Resistance and Antimicrobial Susceptibility Profile of *Staphylococcus saprophyticus* and Other *Staphylococci* Isolated from Patients with Urinary Tract Infection." *Chemotherapy* 58: 482–91. <https://doi.org/10.1159/000346529>.
5. Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/10.2374>.
6. Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Gonçalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant *Staphylococcus Epidermidis* ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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Litteratur

1. Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Ferreira, Adriano M., Mariana F. Bonesso, Alessandro L. Mondelli, Carlos H. Camargo, and Maria De Lourdes R.S. Cunha. 2012. "Oxacillin Resistance and Antimicrobial Susceptibility Profile of Staphylococcus saprophyticus and Other Staphylococci Isolated from Patients with Urinary Tract Infection." *Chemotherapy* 58: 482–91. <https://doi.org/10.1159/000346529>.
5. Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
6. Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Goncalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant Staphylococcus Epidermidis ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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Viited

- Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ferreira, Adriano M., Mariana F. Bonesso, Alessandro L. Mondelli, Carlos H. Camargo, and Maria De Lourdes R.S. Cunha. 2012. "Oxacillin Resistance and Antimicrobial Susceptibility Profile of Staphylococcus saprophyticus and Other Staphylococci Isolated from Patients with Urinary Tract Infection." *Chemotherapy* 58: 482–91. <https://doi.org/10.1159/000346529>.
- Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
- Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Goncalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant Staphylococcus Epidermidis ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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Bibliografia

1. Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Ferreira, Adriano M., Mariana F. Bonesso, Alessandro L. Mondelli, Carlos H. Camargo, and Maria De Lourdes R.S. Cunha. 2012. "Oxacillin Resistance and Antimicrobial Susceptibility Profile of Staphylococcus saprophyticus and Other Staphylococci Isolated from Patients with Urinary Tract Infection." *Chemotherapy* 58: 482–91. <https://doi.org/10.1159/000346529>.
5. Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
6. Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Gonçalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant Staphylococcus Epidermidis ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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Nuorodos

- Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ferreira, Adriano M., Mariana F. Bonesso, Alessandro L. Mondelli, Carlos H. Camargo, and Maria De Lourdes R.S. Cunha. 2012. "Oxacillin Resistance and Antimicrobial Susceptibility Profile of Staphylococcus saprophyticus and Other Staphylococci Isolated from Patients with Urinary Tract Infection." *Chemotherapy* 58: 482–91. <https://doi.org/10.1159/000346529>.
- Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
- Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Gonalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant Staphylococcus Epidermidis ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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1. Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Ferreira, Adriano M., Mariana F. Bonesso, Alessandro L. Mondelli, Carlos H. Camargo, and Maria De Lourdes R.S. Cunha. 2012. "Oxacillin Resistance and Antimicrobial Susceptibility Profile of *Staphylococcus epidermidis* and Other *Staphylococci* Isolated from Patients with Urinary Tract Infection." *Chemotherapy* 58: 482–91. <https://doi.org/10.1159/000346529>.
5. Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
6. Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Goncalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant *Staphylococcus epidermidis* ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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Referencie

- Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Ferreira, Adriano M., Mariana F. Bonesso, Alessandro L. Mondelli, Carlos H. Camargo, and Maria De Lourdes R.S. Cunha. 2012. "Oxacillin Resistance and Antimicrobial Susceptibility Profile of *Staphylococcus saprophyticus* and Other *Staphylococci* Isolated from Patients with Urinary Tract Infection." *Chemotherapy* 58: 482–91. <https://doi.org/10.1159/000346529>.
- Hafeez, Amira, Tehmina Munir, Sara Najeeb, Sabahat Rehman, Mehreen Gilani, Maliha Ansari, and Nadia Saad. 2016. "ICU Pathogens: A Continuous Challenge." *Journal of the College of Physicians and Surgeons Pakistan* 26 (7): 577–80. <https://doi.org/2374>.
- Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Gonçalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant *Staphylococcus Epidermidis* ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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