M.R.S. Agar

REF CM0361B, CM0361R, CM0361T

EN

Intended Use^{IVD}

M.R.S. Agar (CM0361B, CM0361R and CM0361T) is a medium for the growth of *Lactobacillus species* from oral, vaginal and faecal samples. M.R.S. Agar (Dehydrated) (CM0361B, CM0361R and CM0361T) are intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

These devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

The genus *Lactobacillus* was first described by Beijerinck in 1901. It belongs to class *Bacilli* and family *Lactobacillaceae*. *Lactobacilli* belong to lactic acid bacteria (LAB) group, which comprise gram positive, non-spore forming rods or coccobacilli. Their main characteristic is production of lactic acid and they are catalase negative, fermentative, anaerobic or aerotolerant, aciduric or acidophilic. These bacteria are found in food, or mucosal membrane of humans and animals (oral cavity, intestine and vagina). Some strains are considered probiotic, which confer health benefits. *Lactobacillus* are avirulent in general – normal resident of gastrointestinal tract of human. However, some have been associated with infections in immunocompromised patients with health complications, such as bacteraemia, endocarditis, and vascular graft infection. *Lactobacilli* have also been involved in development of dental caries due to their aciduric property¹.

Principle of Method

The medium contains peptone, Lab Lemco powder and yeast extract to supply nutrients. Glucose is the carbohydrate source. The medium is supplemented with sorbitan mono-oleate, sodium acetate, manganese sulphate and magnesium sulphate as growth factors.

Typical Formula

	<u>grams per litre</u>
Peptone	10.0
'Lab-Lemco' powder	8.0
Yeast extract	4.0
Glucose	20.0
Sorbitan mono-oleate	1.0ml
Dipotassium hydrogen phosphate	2.0
Sodium acetate 3H ₂ O	5.0
Triammonium citrate	2.0
Magnesium sulphate 7H ₂ O	0.2
Manganese sulphate 4H ₂ O	0.05
Agar	10.0

Materials Provided

CM0361B: 500g of dehydrated M.R.S. Agar CM0361R: 2.5kg of dehydrated M.R.S. Agar CM0361T: 5kg of dehydrated M.R.S. Agar

500g of dehydrated M.R.S. Agar yields approximately 8.1L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.

- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 04 and B 28.

Procedure

Suspend 62g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Mix well and pour into sterile containers.

Interpretation

Once the media is reconstituted, *Lactobacillus spp.* will show as pale straw colonies.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 72 h @ 30°C, microaerophilic

Positive Controls	
Inoculum level 10 - 100 cfu Colony count is 70 - 130% of the control medium count	
<i>Lactobacillus sakei</i> ATCC® 15521™	0.5-2 mm pale straw colonies
<i>Lactobacillus gasseri</i> ATCC® 19992™	0.5-2 mm pale straw colonies
<i>Lactococcus lactis</i> ATCC® 19435™	0.5-2 mm pale straw colonies
<i>Pediococcus pentosaceus</i> ATCC® 33316™	Pinpoint – 3 mm pale straw colonies

Limitations

Identifications are presumptive and should be confirmed.

The medium is not highly selective and other organisms may grow. Identification is presumptive and isolates should be confirmed using appropriate methods. Due to varying nutritional requirements, some lactic acid bacteria may be encountered which exhibit weak growth or fail to grow on this medium. It is important that adequate moisture vapour is present in the atmosphere above the agar because drying of the plates during incubation will concentrate the selective factors on the surface and make the medium inhibitory.

Performance Characteristics


Accuracy has been demonstrated through review of the QC data. Correct detection of *Lactobacillus* species is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of M.R.S. Agar (CM0361B/R/T) was demonstrated by an overall pass rate of 100% over six months of testing (November 2021 – May 2022). This shows that the performance is reproducible.

MRS Agar (CM0361B/R/T) device is tested in-house as part of the QC process since the products were launched in 1995. For target organisms, when using 10-100 cfu inoculum of *Lactobacillus sakei*, *Lactobacillus gasseri*, *Lactococcus lactis* and *Pediococcus pentosaceus* and incubating the device at 30°C for 72 hours, under microaerophilic conditions, the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

1. Liu, Dongyou. 2011. "Molecular Detection of Human Bacterial Pathogens ." In *Molecular Detection of Human Bacterial Pathogens* , by Dongyou Liu, 258.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Modifications Introduced
2.0	2023-12-15

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Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien