<https://www.thermofisher.com>

EN

Oxoid™ Norfloxacin Discs 10 µg (NOR10)

REF CT0434B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Norfloxacin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs that contain specific amounts of Norfloxacin. The discs are labelled on both sides with details of the antimicrobial (NOR10) and amount present (µg): NOR10 (10 µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Norfloxacin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. They are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients. These discs are used to determine susceptibility against bacteria for which norfloxacin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current CLSI M100 breakpoint tables (all associated with **10 µg** Norfloxacin):

Gram-Positive

- *Enterococcus* spp.
- *Staphylococcus* spp.

Gram-Negative

- *Enterobacterales*
- *Pseudomonas aeruginosa*

Species with published breakpoints according to current EUCAST breakpoint tables (all associated with **10 µg** Norfloxacin):

Gram-Positive

- *Staphylococcus* spp. (screen only)
- *Enterococcus* spp. (screen only)
- *Streptococcus* groups A, B, C and G (screen only)
- *Streptococcus pneumoniae* (screen only)
- *Aerococcus sanguincola* and *urinae* (screen only)
- *Bacillus* spp. (screen only)

Gram-Negative

- *Enterobacterales* (uncomplicated UTI only)

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Norfloxacin AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology, refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the

discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Norfloxacin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Norfloxacin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Norfloxacin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

Table 1. Materials Provided with CT0434B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Pale beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Norfloxacin AST Disc Reagents	
Reagent	Description of Function
Norfloxacin	Norfloxacin belongs to the broad-spectrum class of antibiotics, fluoroquinolones. Norfloxacin's mechanism of action causes the inhibition of the bacterial DNA gyrase, an enzyme essential for successful DNA replication. This subsequently inhibits DNA replication, transcription and repair, ultimately leading to cell lysis.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Norfloxacin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Norfloxacin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 1. Raw data of zone sizes, in accordance with CLSI methodology, taken from the 3 most recent batches of Norfloxacin (10 µg) NOR10 Antimicrobial Susceptibility Test Discs (CT0434B).

CLSI Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3195895	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	28	28	28	-3.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	17	28	22.5	25	25	25	2.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	22	29	25	27	27	27	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5 % SB	15	21	18	19	19	19	1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3179511	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	29	29	29	-2.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	17	28	22.5	23	23	23	0.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	22	29	25	26	26	26	1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5 % SB	15	21	18	20	20	20	2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3154432	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	29	29	29	-2.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	17	28	22.5	23	23	23	0.5 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	22	29	25	27	27	27	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5 % SB	15	21	18	21	21	21	3 ± 0 (CV=0)

Table 2. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the 3 most recent batches of Norfloxacin (10 µg) NOR10 Antimicrobial Susceptibility Test Discs (CT0434B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3195895	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	28	28	28	-3.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	29213	MHA	18	24	21.5	24	24	24	2.5 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	18	24	21.5	20	20	20	-1.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
3179511	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	29	29	29	-2.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	18	24	21.5	21	21	21	-0.5 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	18	24	21.5	22	22	22	0.5 ± 0 (CV=0)
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV) (%)
							1	2	3	
3154432	<i>Escherichia coli</i>	25922	MHA	28	35	31.5	29	29	29	-2.5 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	18	24	21.5	24	24	24	2.5 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	18	24	21.5	22	22	22	0.5 ± 0 (CV=0)

Batch	Organism	ATCC Number	Test Media	Lower Limit	Upper Limit	Reading (mm)			Average (mm) + SD
						1	2	3	
3154432	<i>Escherichia coli</i>	ATCC 25922	MHA	28	35	29	29	29	29 ± 0
	<i>Staphylococcus aureus</i>	ATCC 29213	MHA	18	24	24	24	24	24 ± 0
	<i>Streptococcus pneumoniae</i>	ATCC 49619	MHA-F	18	24	22	22	22	22 ± 0

Clinical Performance Characteristics

A retrospective study analysed bacterial pathogens.⁴ Antimicrobial susceptibility testing was done on wound swabs using disc diffusion technique per the standard of Kirby-Bauer method. The diameter of growth inhibition was measured, and samples were classified as sensitive, intermediate, or resistant according to Clinical and Laboratory Standards Institute. All the antibiotic discs, including the norfloxacin discs, were purchased from Oxoid. The study showed that norfloxacin was effective against *S. aureus*, with resistance of only 8.0 %.

A cross-sectional study analysed conjunctival swabs from communities in West Gojjam, Ethiopia from November 2009 to February 2010.⁵ The samples were transported to the Bahir Dar Regional Health Research Laboratory within 24 hours in transport media. Antimicrobial susceptibility tests were done using disc diffusion technique as per the standard of Kirby-Bauer method using Clinical and Laboratory Standards Institute standards. All the antibiotic discs, including norfloxacin discs, were purchased from Oxoid. The study demonstrated 86-94.4 % sensitivity to ciprofloxacin and norfloxacin. Norfloxacin showed low resistance levels to *S. aureus* and was effective against *E. coli*, *Pseudomonas* spp., and *Enterobacter* spp. The authors concluded that fluoroquinolones, norfloxacin, among others, are important antimicrobial treatment options for human infections.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References


- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Azene, M. K. and Beyene, B. A. (2011) Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Tanzania Journal of Health Research.*, 13 (4): 1-11.
- Abera, B. and Kibret, M.(2014). Azithromycin, Fluoroquinolone and Chloramphenicol Resistance of Non-chlamydia Conjunctival Bacteria in Rural Community of Ethiopia, *Indian Journal of Ophthalmology*, 62 (2): 236-239.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number
	Catalog Number


	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period-after-opening – the number of days the product can be used after opening
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged
	Authorized representative in the European Community/ European Union
	Unique device identifier
	European Conformity Mark
	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

Revision Information

Version	Date of issue and modifications introduced
2.0	2024-07-06

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Azene, M. K. and Beyene, B. A. (2011) Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Tanzania Journal of Health Research*, 13 (4): 1-11.
- Abera, B. and Kibret, M. (2014). Azithromycin, Fluoroquinolone and Chloramphenicol Resistance of Non-chlamydia Conjunctival Bacteria in Rural Community of Ethiopia, *Indian Journal of Ophthalmology*, 62 (2): 236-239.

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Litteratur

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Azene, M. K. and Beyene, B. A. (2011) Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Tanzania Journal of Health Research.*, 13 (4): 1-11.
- Abera, B. and Kibret, M. (2014). Azithromycin, Fluoroquinolone and Chloramphenicol Resistance of Non-chlamydia Conjunctival Bacteria in Rural Community of Ethiopia, *Indian Journal of Ophthalmology*, 62 (2): 236-239.

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Referanser

- Clinical and Laboratory Standards Institute (CLSI). Se siste utgave:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Se siste versjon:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Azene, M. K. and Beyene, B. A. (2011) Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Tanzania Journal of Health Research*, 13 (4): 1-11.
- Abera, B. and Kibret, M. (2014). Azithromycin, Fluoroquinolone and Chloramphenicol Resistance of Non-chlamydia Conjunctival Bacteria in Rural Community of Ethiopia, *Indian Journal of Ophthalmology*, 62 (2): 236-239.

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Referenser

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Azene, M. K. and Beyene, B. A. (2011) Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, Tanzania Journal of Health Research., 13 (4): 1-11.
- Abera, B. and Kibret, M.(2014). Azithromycin, Fluoroquinolone and Chloramphenicol Resistance of Non-chlamydia Conjunctival Bacteria in Rural Community of Ethiopia, Indian Journal of Ophthalmology, 62 (2): 236-239.

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