Oxoid™ Penicillin G Discs (1 IU and 10 IU) P1& P10

REF CT0152B and CT0043B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Thermo Scientific™ Oxoid™ Penicillin G Antimicrobial Susceptibility Test (AST) Discs are 6mm discs prepared by impregnating high quality absorbent paper with specified concentrations of penicillin. The discs are clearly marked on both sides with the corresponding concentration code — P1 and P10. The code designates the antimicrobial agent (penicillin G) and drug content (1 IU or 10IU).

Discs are supplied in cartridges containing 50 discs each. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack.

Oxoid Penicillin G AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Oxoid Penicillin G AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which penicillin G has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate, or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA³, CLSI M100^{1ac}, or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Staphylococcus* spp.
- *Enterococcus* spp.
- *Streptococcus* spp. (β Haemolytic Group)

Gram-Negative

- *Neisseria gonorrhoeae*

Species with published breakpoints according to current CLSI M100 literature:

Gram-Positive

- N/A

Gram-Negative

- N/A

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Staphylococcus* spp. (*S. aureus*, *S. lugdunensis*)
- *Streptococcus* spp.
- *Streptococcus* groups A, B, C and G
- Viridans group *Streptococci*
- *Streptococcus pneumoniae*
- *Listeria monocytogenes*
- *Corynebacterium* spp.
- *Aerococcus sanguinicola* and *urinae*
- *Haemophilus influenzae*

Gram-Negative

- *Pasteurella multocida*
- *Kingella kingae*

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Oxoid Penicillin G AST Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology, refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognized zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Penicillin G AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Oxoid Penicillin G AST Discs consist of 6mm diameter paper discs impregnated with a specific antimicrobial concentration. The discs are marked on both sides with an alpha-numeric code identifying the antimicrobial agent and concentration. Oxoid Penicillin G AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0152B and CT0043B	
Component Description	Material Description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Penicillin G AST Disc Reagent	
Reagent	Description of Function
Penicillin G	Penicillin G sodium, a water soluble benzylpenicillin, is a white to almost white crystalline powder that is almost odourless. Penicillin G sodium is designated 4-Thia-1-azabicyclo(3.2.0)-heptane-2-carboxylic acid, 3,3-dimethyl-7-oxo-6-[(phenylacetyl)amino]-, [2S-(2 α , 5 α , 6 β)]-, monosodium salt.

The concentration of antibiotic on the AST Disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Penicillin G AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Penicillin G AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2°C to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with CLSI methodology, taken from 3 batches of Penicillin G (10 IU) P10 Antimicrobial Susceptibility Test Discs (CT0043B).

CLSI Methodology										
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) %
							1	2	3	
3639426	<i>Staphylococcus aureus</i>	ATCC® 25923™	MHA	26	37	32.5	37	37	37	5.5 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC® 49226™	GCV	26	34	30	34	34	34	4 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA+5 % SB	24	30	27	30	30	30	3 ± 0 (CV=0)
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) %
							1	2	3	
3661903	<i>Staphylococcus aureus</i>	ATCC® 25923™	MHA	26	37	32.5	34	34	34	2.5 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC® 49226™	GCV	26	34	30	32	32	32	2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA+5 % SB	24	30	27	29	29	29	2 ± 0 (CV=0)
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) %
							1	2	3	
3671143	<i>Staphylococcus aureus</i>	ATCC® 25923™	MHA	26	37	32.5	33	33	33	1.5 ± 0 (CV=0)
	<i>Neisseria gonorrhoeae</i>	ATCC® 49226™	GCV	26	34	30	30	30	30	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA+5 % SB	24	30	27	28	28	28	1 ± 0 (CV=0)

Table 4. Raw data of zone sizes, in accordance with EUCAST methodology, taken from 3 batches of Penicillin G (1 IU) P1 Antimicrobial Susceptibility Test Discs (CT0152B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) %
							1	2	3	
3304679	<i>Staphylococcus aureus</i>	ATCC® 29213™	MHA	12	18	15	14	14	14	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA-F	16	22	19	21	21	21	2 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247™	MHA-F	6	9	7.5	6	6	6	-1.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766™	MHA-F	15	21	18	18	18	18	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) %
							1	2	3	
3294976	<i>Staphylococcus aureus</i>	ATCC® 29213™	MHA	12	18	15	14	14	14	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA-F	16	22	19	18	18	18	-1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247	MHA-F	6	9	7.5	6	6	6	-1.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766™	MHA-F	15	21	18	17	17	17	-1 ± 0 (CV=0)
Product Batch	Organism	ATCC® number™	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) %
							1	2	3	
3265997	<i>Staphylococcus aureus</i>	ATCC® 29213™	MHA	12	18	15	12	12	12	-3 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	ATCC® 49619™	MHA-F	16	22	19	18	18	18	-1 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49247™	MHA-F	6	9	7.5	6	6	6	-1.5 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	ATCC® 49766™	MHA-F	15	21	18	18	18	18	0 ± 0 (CV=0)

Clinical Performance Characteristics

In a study by Vikram Singh et al. in 2012⁴, a comparative assessment of CDS, CLSI disc diffusion and Etest techniques for antimicrobial susceptibility testing of *Neisseria gonorrhoeae* was conducted.⁴ As part of the study, Thermo Scientific™ Oxoid™ Penicillin G Antimicrobial Susceptibility Test Discs (0.5µg) P0.5 CT1168S were used for testing *N. gonorrhoeae* isolates using the CDC disc diffusion technique. Comparison of CDS, CLSI, and Etest for susceptibility testing of *N. gonorrhoeae* agreement for both was found as 100%. Penicillin demonstrated high level of correlation coefficient (0.99) on comparison of both the methods. Antibiotic susceptibility by the Etest method revealed that 105 (35.6%) isolates were resistant to penicillin, out of which 95 (32.2%) were b-lactamase positive and 10 (3.4%) were having chromosomally mediated resistance. Out of the remaining 190 (64.4%) isolates, 156 (52.9%) had reduced susceptibility to penicillin and 34 (11.5%) were susceptible to penicillin. These studies provide evidence that Oxoid Penicillin G AST Discs can perform effectively in a clinical setting.

In 2019, Åhman, Matuschek, and Kahlmeter⁵ published a study that compared the quality of antimicrobial susceptibility discs from various manufacturers.⁵ Discs were obtained from nine manufacturers in 2014 and 2017 to evaluate disc quality. Discs were tested in triplicate on Mueller Hinton medium against EUCAST quality control strains. Mean values were compared with targets and ranges. In both years for Oxoid Penicillin G AST Discs, the mean value was within ± 1 mm of the target value, which demonstrated excellent and consistent disc quality for both years.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.


References

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Vikram Singh, Manju Bala, Monika Kakran, V Ramesh. Comparative assessment of CDS, CLSI disc diffusion and Etest techniques for antimicrobial susceptibility testing of *Neisseria gonorrhoeae*: a 6-year study. *BMJ Open* 2012; 2:e000969. doi:10.1136/bmjopen-2012-000969.
5. J. Åhman, E. Matuschek, and G. Kahlmeter. 2019. "The Quality of Antimicrobial Discs from Nine Manufacturers—EUCAST Evaluations in 2014 and 2017." *Clinical Microbiology and Infection* 25 (3): 346–52. doi.org/10.1016/j.cmi.2018.05.021

© 2025 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number
	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Use-by date
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged
	Authorized representative in the European Community/ European Union
	Unique device identifier
	European Conformity Mark
	UK Conformity Mark
	Period-after-opening – the number of days the product can be used after opening
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

Revision Information

Version	Date of issue
3.0	2025-04-14 Typographical error corrected. Test media in Table 3 & 4 updated

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1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Vikram Singh, Manju Bala, Monika Kakran, V Ramesh. Comparative assessment of CDS, CLSI disc diffusion and Etest techniques for antimicrobial susceptibility testing of *Neisseria gonorrhoeae*: a 6-year study. *BMJ Open* 2012; 2:e000969. doi:10.1136/bmjopen-2012-000969.
5. J. Åhman, E. Matuschek, and G. Kahlmeter. 2019. "The Quality of Antimicrobial Discs from Nine Manufacturers—EUCAST Evaluations in 2014 and 2017." *Clinical Microbiology and Infection* 25 (3): 346–52. doi.org/10.1016/j.cmi.2018.05.021.

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Références

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Vikram Singh, Manju Bala, Monika Kakran, V Ramesh. Comparative assessment of CDS, CLSI disc diffusion and Etest techniques for antimicrobial susceptibility testing of *Neisseria gonorrhoeae*: a 6-year study. *BMJ Open* 2012; 2:e000969. doi:10.1136/bmjopen-2012-000969.
5. J. Ahman, E. Matuschek, and G. Kahlmeter. 2019. "The Quality of Antimicrobial Discs from Nine Manufacturers—EUCAST Evaluations in 2014 and 2017." *Clinical Microbiology and Infection* 25 (3): 346–52. doi.org/10.1016/j.cmi.2018.05.021.

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Nuorodos

1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Vikram Singh, Manju Bala, Monika Kakran, V Ramesh. Comparative assessment of CDS, CLSI disc diffusion and Etest techniques for antimicrobial susceptibility testing of *Neisseria gonorrhoeae*: a 6-year study. *BMJ Open* 2012; 2:e000969. doi:10.1136/bmjopen-2012-000969.
5. J. Åhman, E. Matuschek, and G. Kahlmeter. 2019. “The Quality of Antimicrobial Discs from Nine Manufacturers—EUCAST Evaluations in 2014 and 2017.” *Clinical Microbiology and Infection* 25 (3): 346–52. doi.org/10.1016/j.cmi.2018.05.021.

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Referințe

- Clinical Laboratory Standards Institute (CLSI).
 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Vikram Singh, Manju Bala, Monika Kakran, V Ramesh. Comparative assessment of CDS, CLSI disc diffusion and Etest techniques for antimicrobial susceptibility testing of *Neisseria gonorrhoeae*: a 6-year study. *BMJ Open* 2012; 2:e000969. doi:10.1136/bmjopen-2012-000969.
- J. Åhman, E. Matuschek, and G. Kahlmeter. 2019. "The Quality of Antimicrobial Discs from Nine Manufacturers—EUCAST Evaluations in 2014 and 2017." *Clinical Microbiology and Infection* 25 (3): 346–52. doi.org/10.1016/j.cmi.2018.05.021.

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