www.thermofisher.com

Purified Water

EN

REF BO0209M, BO0209V, EB0209C, BO1336E

Intended Use

A medium for the dilution of microorganisms isolated from clinical and non-clinical samples. Can also be used as a rinse solution for sterility testing and a reconstitution fluid for freeze dried supplements.

The device is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Pathogenic aerobic and anaerobic microorganisms encompass a broad range of microorganisms, which can often be found as part of the normal skin and mucosal flora of humans¹. *Escherichia coli* (*E. coli*) are non-spore forming bacteria that are able to grow in aerobic and anaerobic conditions². Whilst *Staphylococcus aureus* is also a nosocomial pathogen, with invasive medical interventions such as catheterization or immunosuppressant treatment, contributing to the incidence of peripheral intravenous catheter (PIVC) failure and catheter-related bloodstream infections (CRBSI)¹.

Infections from fastidious aerobic and anaerobic microorganisms can have a range of manifestations ranging from superficial to life-threatening^{1,2}. Intestinal pathogens like *Escherichia coli* increase their ability to persist and cause disease by utilizing various mechanisms to outcompete the host intestinal microbiota³. Additionally, infection with *Staphylococcus aureus*, a common commensal of the human upper respiratory tract and skin microbiomes, can cause superficial skin infections such as cellulitis and abscesses to life-threatening illnesses such as pneumonia, toxic shock syndrome and sepsis¹.

Principle of Method

Purified water is intended to be used as a medium for the dilution of microorganisms isolated from clinical and non-clinical samples. It contains no added substances to contribute to the growth of organisms.

Typical Formula

	grams per litre
Deionised water	N/A

Physical Appearance

Colour:	Colourless
Clarity	Clear
Fill weight:	
BO0209M	99-101g.
BO0209V	500.-504g
EB0209C	5.0 – 5.3g
BO1336E	10.0 – 10.3g
pH	1.0 – 14.0

Materials Provided

BO0209M: 10 x 100ml sirop bottles of Purified Water BO0209V: 10 x 500ml sirop bottles of Purified Water EB0209C: 200 x 5ml bijou bottles of Purified Water BO1336E: 24 x 10ml universal bottles of Purified Water

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store product in its original packaging at 2–25°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For in vitro diagnostic use only.

- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or bottles.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) Q 5, B42 and B44.

Procedure

Refer to local recommended guidelines such as (UK SMI) Q 5, B42 and B44

Interpretation

Refer to local recommended guidelines such as (UK SMI) Q 5, B42 and B44.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: Counts are performed on Tryptone Soya Agar incubated at 37° ± 2°C for 18-24 hours

Positive Control

Inoculum level: 50 – 150 cfu
Recovery after holding at 23-27°C for 30 minutes
shall be between 50% and 150% of the initial count.

Escherichia coli ATCC® 8739™

Limitations

Purified Water (BO0209M, BO0209V, EB0209C, BO1336E) devices are only intended to be used as a diluent and do not provide a diagnostic function.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct maintenance of microorganisms is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices. The precision of the devices was demonstrated by an overall pass rate of 100% obtained for the product over one month of testing (01.Jul.2022 – 12.Aug.2022; 10 batches) for BO0209M; over six months of testing (11.Jan.2022 – 29.Jul.2022; 10 batches) for BO0209V; and over one month of testing (01.Sep.2022 – 03.Oct.2022; 10 batches) for EB0209E. The precision of the devices was demonstrated by an overall pass rate of 95% obtained for the product over eleven months of testing (05.Aug.2021 – 18.Jul.2022; 20 batches) for BO1336E. This shows that the performance is reproducible.

Purified Water (BO0209M, BO0209V, EB0209C, BO1336E) devices have been tested in-house as part of the QC process since they were first manufactured; BO0209M in 2001, EB0209C in 2002, BO0209V in 2005 and BO1336E in 2018.

The Purified water is inoculated to achieve 50-150 cfu of *Escherichia coli* (ATCC®8739™) in the required volume. This is plated onto the control medium (Tryptone soya medium) immediately and the Purified water is sampled again after holding at 23-27°C for 30 minutes. The control plates are incubated at 37± 2°C for 18-24 hours. The counts at time zero and 30 minutes are compared, and the user can recover of 50%-150% of the initial inoculum after holding.

Bibliography

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species'. ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>
2. Public Health England. 2015. 'SMI ID 16: identification of Enterobacteriaceae'. ID 17 (3). <https://www.gov.uk/government/publications/smi-id-16-identification-of-enterobacteriaceae>
3. Ormsby, M. J., Logan, M., Johnson, S. A., McIntosh, A., Fallata, G., Papadopolou, R., ... & Wall, D. M. (2019). Inflammation associated ethanolamine facilitates infection by Crohn's disease-linked adherent- invasive Escherichia coli. *EBioMedicine*, 43, 325-332.

Symbol Legend

Symbol/Label	Meaning
	Manufacturer
IVD	In Vitro Diagnostic Medical Device
	Temperature limit
LOT	Batch Code
REF	Catalogue Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
EC REP	Authorized representative in the European Community/ European Union
UDI	Unique device identifier
CE	European Conformity Mark
UK CA	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

© 2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

MBD_BT_IFU-0714

For technical assistance please contact your local distributor.

Revision Information

Version	Date of issue and modifications introduced
3.0	2023-10-24

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Royaume-Uni