<https://www.thermofisher.com>

EN

Oxoid™ Tetracycline Discs 30 µg (TE30)

REFCT0054B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the Oxoid AST generic instructions for use supplied with the product and available online.

Thermo Scientific™ Oxoid™ Tetracycline Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent Tetracycline. The discs are labelled on both sides with details of antimicrobial (TE) and amount present (µg): TE30 µg.

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Oxoid Tetracycline AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disk should only be used once.

Intended Use

Oxoid Tetracycline AST Discs are intended to determine susceptibility against microorganisms for which Tetracycline has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

Used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing, Oxoid Tetracycline AST Discs are manufactured as an in vitro diagnostic device intended to be used in a diagnostic workflow to aid clinicians in treatment options. The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

These discs are suitable for use with fresh and pure clinical isolates from culture media.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100,^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- *Staphylococcus* spp.
- *Streptococcus pneumoniae*
- Group A β -hemolytic *streptococci*

Gram-Negative

- *Enterobacteriaceae*
- *Acinetobacter* spp.
- *Haemophilus influenzae* and *parainfluenzae*
- *Neisseria gonorrhoeae*

Species with published breakpoints according to current CLSI literature:

Gram-Positive

- *Staphylococcus* spp.
- *Enterococcus* spp.
- *Streptococcus pneumoniae*
- *Streptococcus* spp. β -Hemolytic Group
- *Streptococcus* spp. Viridans Group

Gram-Negative

- *Enterobacterales*
- *Acinetobacter* spp.
- *Burkholderia cepacia* complex
- *Haemophilus influenzae* and *parainfluenzae*
- *Neisseria gonorrhoeae*

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Staphylococcus* spp.
- *Streptococcus* groups A, B, C and G
- *Streptococcus pneumoniae*
- *Corynebacterium* spp.

Gram-Negative

- *Haemophilus influenzae*
- *Moraxella catarrhalis*
- *Pasteurella multocida* (Screen only)
- *Campylobacter jejuni* and *coli*
- *Kingella kingae*
- *Burkholderia pseudomallei* (Screen only)

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Oxoid Tetracycline AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Tetracycline AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Oxoid Tetracycline AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. The discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1 . Materials Provided with CT0054B	
Component Description	Material Description
Cartridge	K-Resin plastic cartridge for discs.
Spring/Cap/Plunger	Assembly components (cartridge).
Desiccant	Molecular sieve tablets 500mg. Light beige to brown small lozenge shaped FBE tablets.
Foil	Hard temper foil. Unprinted.

Table2. Description of Tetracycline AST Disc Reagent	
Reagent	Description of Function
Tetracycline	Tetracycline is a broad-spectrum polyketide antibiotic produced by the <i>Streptomyces</i> genus of Actinobacteria. It exerts a bacteriostatic effect on bacteria by binding reversibly to the bacterial 30S ribosomal subunit and blocking incoming aminoacyl-tRNA from binding to the ribosome acceptor site. It also binds to some extent to the bacterial 50S ribosomal subunit and may alter the cytoplasmic membrane causing intracellular components to leak from bacterial cells.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Oxoid Tetracycline AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Tetracycline AST Discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be

stored within the container at 2 - 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Results Summary for Tested QC Isolates in Accordance with EUCAST Methodology

EUCAST Methodology										
Product CT Code	Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Reading (mm)			Mean difference between experimental values and reference mid-point values (mm) + SD [Coefficient of variance (CV)]
							1	2	3	
CT0054B	3158883	<i>Staphylococcus aureus</i>	29213	MHA	23	31	26	26	26	26±0%
		<i>Streptococcus pneumoniae</i>	49619	MHA-F	28	34	30	30	30	30±0%
		<i>Haemophilus influenzae</i>	49766	MHA-F	28	34	31	31	31	31±0%
		<i>Campylobacter jejuni</i>	33560	MHA-F	30	38	33	34	34	33.7±0.6%
	3138683	<i>Staphylococcus aureus</i>	29213	MHA	23	31	27	27	27	27±0%
		<i>Streptococcus pneumoniae</i>	49619	MHA-F	28	34	32	32	32	32±0%
		<i>Haemophilus influenzae</i>	49766	MHA-F	28	34	29	29	29	29±0%
		<i>Campylobacter jejuni</i>	33560	MHA-F	30	38	35	35	35	35±0%
	3145157	<i>Staphylococcus aureus</i>	29213	MHA	23	31	28	28	28	28±0%
		<i>Streptococcus pneumoniae</i>	49619	MHA-F	28	34	33	33	33	33±0%
		<i>Haemophilus influenzae</i>	49766	MHA-F	28	34	31	31	31	31±0%
		<i>Campylobacter jejuni</i>	33560	MHA-F	30	38	38	38	38	38±0%

Table 4. Results Summary for Tested QC Isolates in Accordance with CLSI Methodology

CLSI Methodology										
Product CT Code	Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Reading (mm)			Average (mm)+ SD(%)
							1	2	3	
CT0054B	3158883	<i>Escherichia coli</i>	25922	MHA	18	25	23	23	23	23±0%
		<i>Staphylococcus aureus</i>	25923	MHA	24	30	28	28	28	28±0%
		<i>Haemophilus influenzae</i>	49247	HTM	14	22	21	19	19	19.7±1.2%
		<i>Neisseria gonorrhoeae</i>	49226	GCV	30	42	39	39	39	39±0%

	3138683	<i>Streptococcus pneumoniae</i>	49619	MHA + 5%HB	27	31	30	30	30	30±0%
		<i>Escherichia coli</i>	25922	MHA	18	25	25	25	25	25±0%
		<i>Staphylococcus aureus</i>	25923	MHA	24	30	29	29	29	29±0%
		<i>Haemophilus influenzae</i>	49247	HTM	14	22	19	19	19	19±0%
		<i>Neisseria gonorrhoeae</i>	49226	GCV	30	42	42	42	42	42±0%
		<i>Streptococcus pneumoniae</i>	49619	MHA + 5%HB	27	31	31	31	31	31±0%
	3145157	<i>Escherichia coli</i>	25922	MHA	18	25	23	23	23	23±0%
		<i>Staphylococcus aureus</i>	25923	MHA	24	30	27	27	27	27±0%
		<i>Haemophilus influenzae</i>	49247	HTM	14	22	19	19	19	19±0%
		<i>Neisseria gonorrhoeae</i>	49226	GCV	30	42	38	38	38	38±0%
		<i>Streptococcus pneumoniae</i>	49619	MHA + 5%HB	27	31	31	31	31	31±0%

Clinical Performance Characteristics

A study compared the ADAGIO automated system (Bio-Rad) against manual reading using a ruler when interpreting disk diffusion AST results on a collection of different fastidious bacteria.⁴ The study used 12 different Oxoid antibiotics including tetracycline (10 µg) to test 144 clinical isolates by disk diffusion according to EUCAST guidelines. Up to 12 antibiotics were tested for each group of microorganisms depending on the number of antibiotics for which specific EUCAST breakpoints were available. In conclusion this study provided evidence of the good performance of Oxoid Tetracycline AST Discs against *S. pneumoniae* and *H. influenzae* clinical isolates.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Idelevich, Evgeny A., Karsten Becker, Janne Schmitz, Dennis Knaack, Georg Peters, and Robin Köck. 2016. "Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria." PLoS ONE 11 (7). <https://doi.org/10.1371/journal.pone.0159183>

Glossary of symbols


Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device

	Temperature limitation (storage temp.)
LOT	Lot / Batch Number
REF	Catalog Number
	Do not re-use
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Expiry Date
	Period-after-opening – the number of days the product can be used after opening
R_x only	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged
EC REP	Authorized representative in the European Community/ European Union
UDI	Unique device identifier
CE	European Conformity Mark
UK CA	UK Conformity Mark
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

Revision Information

Version	Date of issue and modifications introduced
2.0	2024-10-16 Minor formatting changes

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Referencije

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Idelevich, Evgeny A., Karsten Becker, Janne Schmitz, Dennis Knaack, Georg Peters, and Robin Köck. 2016. "Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria." PLoS ONE 11 (7). <https://doi.org/10.1371/journal.pone.0159183>

© 2023. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Nuorodos

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Idelevich, Evgeny A., Karsten Becker, Janne Schmitz, Dennis Knaack, Georg Peters, and Robin Köck. 2016. “Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria.” PLoS ONE 11 (7). <https://doi.org/10.1371/journal.pone.0159183>

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).

4. Idelevich, Evgeny A., Karsten Becker, Janne Schmitz, Dennis Knaack, Georg Peters, and Robin Köck. 2016. "Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria." PLoS ONE 11 (7). <https://doi.org/10.1371/journal.pone.0159183>

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Referências

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.

2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Idelevich, Evgeny A., Karsten Becker, Janne Schmitz, Dennis Knaack, Georg Peters, and Robin Köck. 2016. "Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria." PLoS ONE 11 (7). <https://doi.org/10.1371/journal.pone.0159183>

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