	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Vancomycin Discs 5 µg/30 µg (VA5/VA30)

REF CT0188B & CT0058B

Antimicrobial Susceptibility Test Discs

Note: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Vancomycin Antimicrobial Susceptibility Test (AST) Discs are 6mm paper discs which contain specific amounts of the antimicrobial agent vancomycin. The discs are labelled on both sides with details of the antimicrobial (VA) and amount present (µg): VA5 (5µg) or VA30 (30µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Vancomycin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which Vancomycin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA tables, CLSI M100^{1ac}, or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive:

- *Staphylococcus spp.* (*Staphylococcus aureus* and coagulase-negative *Staphylococcus*)
- *Enterococcus spp.*
- *Streptococcus spp.* β-Hemolytic Group
- *Streptococcus spp.* Viridans group

* Exception to the recognised standard of CLSI M100, for *Streptococcus pneumoniae*, susceptibility test interpretive criteria are not recognised at this time.

Gram-Negative

- N/A

Species with published breakpoints according to current CLSI M100 (all associated with **30µg** vancomycin):

Gram-Positive:

- *Enterococcus spp.*
- *Streptococcus pneumoniae*
- *Streptococcus spp.* (β-Hemolytic Group)
- *Streptococcus spp.* Viridans Group

Gram-Negative

- N/A

Species with published breakpoints according to current EUCAST breakpoint tables (all associated with **5µg** vancomycin):

Gram-Positive:

- *Enterococcus spp.*

- *Streptococcus* groups A, B, C and G
- *Streptococcus pneumoniae*
- Viridans group streptococci.
- *Corynebacterium spp.*
- *Aerococcus sanguinicola* and *urinae*

Gram-Negative

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principles of Method

Vancomycin AST Discs may be used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological traceability of calibrator and control material values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Vancomycin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Vancomycin AST Discs consist of 6mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Vancomycin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, please refer to table 2.

Table 1. Materials provided with CT0188B & CT0058B	
Component description	Material description
Cartridge with spring, cap and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge, 5 cartridges per pack.

Table2. Description of Vancomycin AST Disc reagents	
Reagent	Description of Function
Vancomycin-Nonsterile (5µg – CT0188B)	Glycopeptide antimicrobial agent active against MRSA and enterococcal infections. Vancomycin inhibits cell wall formation by interrupting peptidoglycan synthesis.
Vancomycin-Nonsterile (30µg – CT0058B)	

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Vancomycin AST Discs can be dispensed using a Disc Dispenser, which is not included with the device.

Shelf-life & Storage Conditions

Unopened cartridges of Vancomycin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided (with an unsaturated (orange) desiccant), or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2 - 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Limitations

Non-wild type isolates were not available when developing the disk diffusion method for *Streptococcus* groups A, B, C and G, *Streptococcus pneumoniae*, *Corynebacterium spp.*, Viridans group streptococci, *Aerococcus sanguinicola* and *Aerococcus urinae*.

Vancomycin susceptible *Enterococci* exhibit sharp zone edges and do not exhibit colonies in the inhibition zone. Examine zone edges with transmitted light (plate held up to light). If the zone edge is fuzzy, colonies grow within the zone or if you are uncertain,

then perform confirmatory testing with PCR or report resistant even if the zone diameter is ≥ 12 mm. Isolates must not be reported susceptible before 24 h incubation.

Non-susceptible isolates of the Viridans group *Streptococci* are rare or not yet reported. The identification and antimicrobial susceptibility test result on any such isolate must be confirmed and the isolate sent to a reference laboratory.

Analytical Performance Characteristics

Table 3. Results summary for tested QC isolates, in accordance with EUCAST methodology

EUCAST Methodology										
Product Batch	Organism	ATCC number	Test Media	Reference Lower Limit	Reference Upper Limit	Reference mid-point value	Readings (mm)			Mean difference between experimental values and reference mid-point values (mm) + SD [Coefficient of variance (CV)]
							1	2	3	
3108134	<i>Enterococcus faecalis</i>	29212	MHA	10	16	13	13	13	13	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	17	23	20	20	20	20	0 ± 0 (CV=0)
Product Batch	Organism	ATCC number	Test Media	Reference Lower Limit	Reference Upper Limit	Reference mid-point value	Readings (mm)			Mean difference between experimental values and reference mid-point values (mm) + SD (CV)
							1	2	3	
2984645	<i>Enterococcus faecalis</i>	29212	MHA	10	16	13	13	13	13	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	17	23	20	18	18	18	-2 ± 0 (CV=0)
Product Batch	Organism	ATCC number	Test Media	Reference Lower Limit	Reference Upper Limit	Reference mid-point value	Readings (mm)			Mean difference between experimental values and reference mid-point values (mm) + SD (CV)
							1	2	3	
2972605	<i>Enterococcus faecalis</i>	29212	MHA	10	16	13	13	13	13	0 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA-F	17	23	20	20	20	20	0 ± 0 (CV=0)

Table 4. Results summary for tested QC isolates, in accordance with CLSI methodology

CLSI Methodology										
Product Batch	Organism	ATCC number	Test Media	Reference Lower Limit	Reference Upper Limit	Reference mid-point value	Readings (mm)			Mean difference between experimental values and reference mid-point values (mm) + SD [Coefficient of variance (CV)]
							1	2	3	
3111869	<i>Staphylococcus aureus</i>	25923	MHA	17	21	19	17	17	17	-2 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5 % SB	20	27	23.5	22	22	22	-1.5 ± 0 (CV=0)
Product Batch	Organism	ATCC number	Test Media	Reference Lower Limit	Reference Upper Limit	Reference mid-point value	Readings (mm)			Mean difference between experimental values and reference mid-point values (mm) + SD (CV)
							1	2	3	
3111868	<i>Staphylococcus aureus</i>	25923	MHA	17	21	19	18	18	18	-1 ± 0 (CV=0)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5 % SB	20	27	23.5	22	22	22	-1.5 ± 0 (CV=0)
Product Batch	Organism	ATCC number	Test Media	Reference Lower Limit	Reference Upper Limit	Reference mid-point value	Readings (mm)			Mean difference between experimental values and reference mid-point values (mm) + SD (CV)
							1	2	3	
2947426	<i>Staphylococcus aureus</i>	25923	MHA	17	21	19	20	20	19	0.67 ± 0.58 (CV=0.87)
	<i>Streptococcus pneumoniae</i>	49619	MHA + 5 % SB	20	27	23.5	21	21	21	-2.5 ± 0 (CV=0)

Clinical performance characteristics

Clinical performance of Vancomycin AST discs has been reported in numerous sources, all supporting the routine clinical use to determine antimicrobial susceptibility.

A double-centre retrospective study investigated blood cultures positive for *S. aureus* obtained from 1st January 2001 to 31 December 2014⁶. Antimicrobial susceptibility testing used Kirby–Bauer disk diffusion methods and used various antibiotic discs from Oxoid, including vancomycin. The study found that *S. aureus* had relatively low drug resistance to most antibiotics investigated with Vancomycin discs being an effective method of analysing resistance in *S. aureus*.

A study⁴ was run to evaluate automated reading systems against Oxoid antimicrobial discs, including vancomycin. 144 clinical isolates were included in AST susceptibility testing of vancomycin discs via disk diffusion according to the EUCAST breakpoints guidelines. Vancomycin discs were shown to be effective with clinical analysis of the *S. pneumoniae* and had good categorical agreement with the automated reading system tested. Additionally, a study⁵ consisting of 345 consecutive isolates of *S. aureus* collected from Başkent University Hospital showed that the most effective antibiotics against MRSA are glycopeptide antibiotics, and of 187 MRSA strains, no resistance was acquired against vancomycin.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established. In the event of malfunction do not use device.

References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests (–)
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. E A Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
5. L Savaş *et al.* (2005). Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Turk. J. Med. Sci.*, 35: 323-327.
6. J Zhou *et al.* (2015). *Staphylococcus aureus* bacteremias following liver transplantation: a clinical analysis of 20 cases. *Ther. Clin. Risk Management*, 11: 933-937.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device

	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-18

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. E A Idevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
5. L Savaş *et al.* (2005). Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Turk. J. Med. Sci.*, 35: 323-327.
6. J Zhou *et al.* (2015). *Staphylococcus aureus* bacteremias following liver transplantation: a clinical analysis of 20 cases. *Ther. Clin. Risk Management*, 11: 933-937.

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Referencije

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- E A Ilevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
- L Savaş *et al.* (2005). Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Turk. J. Med. Sci.*, 35: 323-327.
- J Zhou *et al.* (2015). *Staphylococcus aureus* bacteremias following liver transplantation: a clinical analysis of 20 cases. *Ther. Clin. Risk Management*, 11: 933-937.

© 2023. Thermo Fisher Scientific Inc. Sva prava pridržana. ATCC® je zaštitni znak društva ATCC. Svi ostali zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific Inc. i njegovih podružnica. Ove informacije nisu namijenjene za poticanje na uporabu tih proizvoda na način kojim se mogu povrijediti prava intelektualnog vlasništva drugih društava.

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Reference

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- E A Ilevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
- L Savaş *et al.* (2005). Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, Turk. J. Med. Sci., 35: 323-327.
- J Zhou *et al.* (2015). *Staphylococcus aureus* bacteremias following liver transplantation: a clinical analysis of 20 cases. Ther. Clin. Risk Management, 11: 933-937.

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Litteratur

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- E A Idevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
- L Savaş *et al.* (2005). Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Turk. J. Med. Sci.*, 35: 323-327.
- J Zhou *et al.* (2015). *Staphylococcus aureus* bacteremias following liver transplantation: a clinical analysis of 20 cases. *Ther. Clin. Risk Management*, 11: 933-937.

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Viited

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. E A Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
5. L Savaş *et al.* (2005). Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Turk. J. Med. Sci.*, 35: 323-327.
6. J Zhou *et al.* (2015). *Staphylococcus aureus* bacteremias following liver transplantation: a clinical analysis of 20 cases. *Ther. Clin. Risk Management*, 11: 933-937.

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Viitteet

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. E A Iidelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
5. L Savaş *et al.* (2005). Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Turk. J. Med. Sci.*, 35: 323-327.
6. J Zhou *et al.* (2015). *Staphylococcus aureus* bacteremias following liver transplantation: a clinical analysis of 20 cases. *Ther. Clin. Risk Management*, 11: 933-937.

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- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- E A Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
- L Savas *et al.* (2005). Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Turk. J. Med. Sci.*, 35: 323-327.
- J Zhou *et al.* (2015). *Staphylococcus aureus* bacteremias following liver transplantation: a clinical analysis of 20 cases. *Ther. Clin. Risk Management*, 11: 933-937.

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Atsauces

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- E A Idelevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
- L Savaş *et al.* (2005). Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Turk. J. Med. Sci.*, 35: 323-327.
- J Zhou *et al.* (2015). *Staphylococcus aureus* bacteremias following liver transplantation: a clinical analysis of 20 cases. *Ther. Clin. Risk Management*, 11: 933-937.

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Referencias

- Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- E A Idevlevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
- L Savaş *et al.* (2005). Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Turk. J. Med. Sci.*, 35: 323-327.
- J Zhou *et al.* (2015). *Staphylococcus aureus* bacteremias following liver transplantation: a clinical analysis of 20 cases. *Ther. Clin. Risk Management*, 11: 933-937.

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Referenser

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. E A Ilevich *et al.* (2016). Evaluation of an Automated System for Reading and Interpreting Disk Diffusion Antimicrobial Susceptibility Testing of Fastidious Bacteria. PLoS ONE 11(7): e0159183.
5. L Savaş *et al.* (2005). Prospective Analysis of Antibiotic Susceptibility Patterns of MRSA in a Turkish University Hospital, *Turk. J. Med. Sci.*, 35: 323-327.
6. J Zhou *et al.* (2015). *Staphylococcus aureus* bacteremias following liver transplantation: a clinical analysis of 20 cases. *Ther. Clin. Risk Management*, 11: 933-937.

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