Vitox Supplement

REF

SR0090A & SR0090H

EN

* This instructions for use (IFU) document is intended to be read in conjunction with the IFU for GC Agar Base (CM0367B),

Intended Use

Vitox Supplements (SR0090A and SR0090H) are enrichment supplements intended to be used in the preparation of Thayer Martin medium using Thermo Scientific™ GC Agar Base (CM0367) for the isolation of *Neisseria* and *Haemophilus* species from clinical samples (e.g., wounds, throat, genital, nose, etc).

Vitox Supplements (SR0090A and SR0090H) are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Neisseria* and *Haemophilus* infections.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Neisseria species are Gram-negative, non-motile, obligate human pathogens with no other natural host¹ (Public Health England 2015). The two main pathogens of this group are

N. gonorrhoeae and *N. meningitidis*, although some other species of *Neisseria* have been implicated as causes of infection in immunocompromised patients¹ (Public Health England 2015).

Haemophilus species are Gram-negative, non-motile bacteria² (Public Health England 2021). Of the 25 species of *Haemophilus* genus, 10 species, including *H. influenzae*, are known to be associated with humans² (Public Health England 2021). *H. influenzae* is well known for causing invasive diseases in younger population and people with underlying medical condition³ (ECDC 2014). *H. influenzae* has particularly severe manifestations in immunocompromised children and represents one of the leading causes of meningitis in children, a disease associated with severe sequelae and high mortality⁴ (National Institute for Health and Care Excellence 2010). *H. influenzae* tends to colonise the upper respiratory tract, causing a range of diseases from superficial skin infections such as cellulitis, to meningitis and bacteremia serotype² (Public Health England 2021).

Consequently, it is important to isolate *Neisseria* and *Haemophilus* species in clinical samples as some of these organisms, such as *N. gonorrhoeae*, *N. meningitidis* and *H. influenzae*, are considerable human pathogens associated with considerable morbidity and mortality^{1,2,5} (Public Health England 2014; 2015; 2021).

Principle of Method

Vitox Supplements are intended to be used in the preparation of Thayer Martin medium using GC Agar Base (CM0367B) for the isolation of *Neisseria* and *Haemophilus* species from clinical samples. The device contains glucose as a carbohydrate source supplying energy for bacterial growth. Essential vitamins, minerals, and amino acids are added as necessary factors to support the growth of *Neisseria* and *Haemophilus* strains.

Typical Formula

SR0090A Vitox Supplement:

mg/500ml at 2%v/v

Distilled Water	10.0ml
Glucose	1.0g
Vitamin B ₁₂	0.1
Adenine	10.0
L-Glutamine	100.0
Guanine	0.3
p-Aminobenzoic acid	0.13
L-Cysteine	11.0
NAD (Coenzyme 1)	2.5
Coccarboxylase	1.0
Iron (III) nitrate	0.2
Thiamine hydrochloride	0.03
Cysteine hydrochloride	259.0

SR0090H Vitox Supplement:

mg/2L at 2% v/v

Distilled Water	40.0ml
Glucose	4.0g
Vitamin B ₁₂	0.4
Adenine	40.0
L-Glutamine	400.0
Guanine	1.2
p-Aminobenzoic acid	0.52
L-Cysteine	44.0
NAD (Coenzyme 1)	10.0
Coccarboxylase	4.0
Iron (III) nitrate	0.8
Thiamine hydrochloride	0.12

Cysteine hydrochloride

1036.0

Materials Provided

SR0090A: 5 x rehydration Fluid (SR0090B) and 5 vitox vials (SR0090C)
SR0090H: 5 x 40ml rehydration fluid (SR0090J) and 5 vitox vials (SR0090K)

Materials Required but Not Supplied

- GC Agar Base (dehydrated) product code CM0376B
- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store product in its original packaging between 2°C and 8°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Warnings and Precautions

- For in vitro diagnostic use only.
- Wear protective gloves/protective clothing/eye protection/face protection
- Each vial is for single use. Do not re-use.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or vials.
- Do not use the product beyond the stated expiry date. Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 6, ID 12 and Q 5

Procedure

SR0090A:

- Allow supplement to equilibrate to room temperature before use.
- Add one vial of hydration fluid(SR0090B) to 1 vial of vitox (SR0090C)
- Aseptically add the vial contents to 1 litre (1% v/v) or 500ml (2% v/v) of sterile GC Agar Base (CM0367B) prepared as directed and cooled to 50°C.
- Mix well and aseptically dispense into suitable containers.


SR0090H:

- Allow supplement to equilibrate to room temperature before use.
- Add one vial of hydration fluid(SR0090J) to 1 vial of vitox (SR0090K)
- Aseptically add the vial contents to 4 litres (1% v/v) or 2 litres (2% v/v) of sterile GC Agar Base (CM0367B) prepared as directed and cooled to 50°C.
- Mix well and aseptically dispense into suitable containers.

Bibliography

1. Public Health England. 2015. 'Identification of Neisseria species' UK Standards for Microbiology Investigations. ID 6 Issue 7. <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>
2. Public Health England. 2021. "Identification of Haemophilus Species and the HACEK Group of Organisms.' UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>
3. European Centre for Disease Prevention and Control. 2014. 'Invasive Haemophilus influenzae disease' <https://www.ecdc.europa.eu/en/invasive-haemophilus-influenzae-disease/facts>
4. National Institute for Health and Care Excellence. 2010. 'Meningitis (Bacterial) and Meningococcal Septicaemia in under 16s: Recognition, Diagnosis and Management.'
5. Public Health England. 2014. 'Screening for *Neisseria Meningitidis*.' UK Standards for Microbiology Investigations. B 51 Issue 2. <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2024-05-24

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1. Public Health England. 2015. 'Identification of Neisseria species' UK Standards for Microbiology Investigations. ID 6 Issue 7. <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>
2. Public Health England. 2021. 'Identification of Haemophilus Species and the HACEK Group of Organisms.' UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>
3. European Centre for Disease Prevention and Control. 2014. 'Invasive Haemophilus influenzae disease' <https://www.ecdc.europa.eu/en/invasive-haemophilus-influenzae-disease/facts>
4. National Institute for Health and Care Excellence. 2010. 'Meningitis (Bacterial) and Meningococcal Septicaemia in under 16s: Recognition, Diagnosis and Management.'
5. Public Health England. 2014. 'Screening for Neisseria Meningitidis.' UK Standards for Microbiology Investigations. B 51 Issue 2. <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>

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
Litteratur

1. Public Health England. 2015. 'Identification of *Neisseria* species' UK Standards for Microbiology Investigations. ID 6 Issue 7. <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>
2. Public Health England. 2021. 'Identification of *Haemophilus* Species and the HACEK Group of Organisms.' UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>
3. European Centre for Disease Prevention and Control. 2014. 'Invasive *Haemophilus influenzae* disease' <https://www.ecdc.europa.eu/en/invasive-haemophilus-influenzae-disease/facts>
4. National Institute for Health and Care Excellence. 2010. 'Meningitis (Bacterial) and Meningococcal Septicaemia in under 16s: Recognition, Diagnosis and Management.'
5. Public Health England. 2014. 'Screening for *Neisseria Meningitidis*.' UK Standards for Microbiology Investigations. B 51 Issue 2. <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>

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Kirjandus

1. Public Health England. 2015. 'Identification of Neisseria species' UK Standards for Microbiology Investigations. ID 6 Issue 7. <https://www.gov.uk/government/publications/smi-id-6-identification-of-neisseria-species>
2. Public Health England. 2021. 'Identification of Haemophilus Species and the HACEK Group of Organisms.' UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>
3. European Centre for Disease Prevention and Control. 2014. 'Invasive Haemophilus influenzae disease' <https://www.ecdc.europa.eu/en/invasive-haemophilus-influenzae-disease/facts>
4. National Institute for Health and Care Excellence. 2010. 'Meningitis (Bacterial) and Meningococcal Septicaemia in under 16s: Recognition, Diagnosis and Management.'
5. Public Health England. 2014. 'Screening for *Neisseria Meningitidis*.' UK Standards for Microbiology Investigations. B 51 Issue 2. <https://www.gov.uk/government/publications/smi-b-51-screening-for-neisseria-meningitidis>

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien