remel PYR DISK w/ REAGENT

REF R30854301.....	Σ 50
REF R30854401.....	Σ 100

1. INTENDED USE

The PYR Disk is a rapid, colorimetric test for use in qualitative procedures for the presumptive identification of enterococci, group A streptococci and *Escherichia coli* by detection of pyrrolidonyl arylamidase (PYRase). The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

2. SUMMARY AND EXPLANATION

In 1981, Godsey *et al.* described a test to differentiate Group A streptococci and enterococci from other streptococci based on the ability to cleave L-pyrrolidonyl-β-naphthylamide (PYR).¹ In 1982, investigators used the PYR test in conjunction with the CAMP and bile esculin tests to presumptively identify streptococci.² Facklam *et al.* incorporated PYR substrate into agar and tested for enzyme activity after overnight incubation. Bosley *et al.* incorporated PYR substrate into a broth and tested for enzyme hydrolysis after 4 hours incubation.³ In 1985, Ellner *et al.* described a colorimetric method using filter paper strips containing the PYR substrate and found it easy to perform and cost-effective for presumptive identification of Group A streptococci and enterococci.^{4,5} More recently, York *et al.* reported the PYR test can be used for differentiating *E. coli* from other enteric Gram-negative bacilli.⁶ This method is recommended by the Clinical and Laboratory Standards Institute (CLSI) for rapid, abbreviated identification of *E. coli*.⁷

3. PRINCIPLE

Filter paper disks are impregnated with the substrate, L-pyrrolidonyl-β-naphthylamide. In the presence of the enzyme, pyrrolidonyl arylamidase (PYRase), the substrate is hydrolyzed resulting in the formation of β-naphthylamine. This compound produces a red color upon the addition of *p*-dimethylaminocinnamaldehyde (Color Developer). The test is routinely performed in less than four minutes.

4. REAGENTS AND MATERIALS SUPPLIED

1. PYR Disks:

Reactive Ingredient: L-pyrrolidonyl-β-naphthylamide
(1 vial of 50 or 100 disks)

2. Color Developer:

Reactive Ingredient: 0.01%
p-dimethylaminocinnamaldehyde (1 bottle, 4.5 ml)

3. Instructions for Use (IFU)

5. CONTENTS SYMBOLS

PYR Disks	PYR Disks
PYR Color Developer	PYR Color Developer

6. PRECAUTIONS

- This product is for *in vitro* diagnostic use and should be used by properly trained individuals.
- Standard precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and test materials after use.
- Directions should be read and followed carefully.
- Individual test disks are for single use and must be disposed of after use.
- Caution!** The Color Developer is toxic and may cause harm to the environment. Harmful by inhalation, contact with skin or eyes, or if swallowed. May impair fertility or cause harm to unborn child.

- Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.
- Color Developer is provided at the necessary working strength and should be dispensed directly from the dropper bottle.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
- In the event of malfunction do not use device.

Composition / information on ingredients

2-Methoxyethanol 109-86-4

Acetic acid 64-19-7

DANGER

US ONLY

US & EU

H315	Causes skin irritation
H319	Causes serious eye irritation
H335	May cause respiratory irritation
H336	May cause drowsiness or dizziness
H360	May damage fertility. May damage the unborn child
H373	May cause damage to organs through prolonged or repeated exposure
P201	Obtain special instructions before use
P202	Do not handle until all safety precautions have been read and understood
P281	Use personal protective equipment as required
P264	Wash face, hands and any exposed skin thoroughly after handling
P280	Wear eye/face protection
P260	Do not breathe dust/fume/gas/mist/ vapors/spray
P271	Use only outdoors or in a well-ventilated area
P308+313	If exposed or concerned: Get medical attention/advice
P304+P340	If INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing
P302+P352	If ON SKIN: Wash with plenty of soap and water
P332+P313	If skin irritation occurs: Get medical advice/attention
P362	Take off contaminated clothing and wash before reuse
P305+P351 +P338	If IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
P337+P313	If eye irritation persists: Get medical advice/attention
P405	Store locked up
P403+P233	Store in a well-ventilated place. Keep container tightly closed
P501	Dispose of contents/container to an approved waste disposal plant

Hazards not otherwise classified (HNOC)

None identified

WARNING! This product contains a chemical known in the State of California to cause birth defects or other reproductive harm.

Emergency Telephone Number

INFOTRAC - 24 Hour Number: 1-800-535-5053

Outside of the United States, call 24 Hour Number: 001-352-323-3500 (Call Collect)

7. STORAGE

This product is ready for use and no further preparation is necessary. Store product in the dark in its original container at 2-8°C until used. Do not freeze or overheat. Allow product to equilibrate to room temperature before use. Do not incubate prior to use.

8. PRODUCT DETERIORATION

This product should not be used if (1) the color of the disks has changed from white, (2) the expiration date has passed, or (3) there are other signs of deterioration. Protect disks from moisture by removing from the vial only those disks necessary for testing. Promptly replace the cap and return the vial to 2-8°C.

Note: Precipitate formation or color variation from straw to pink in the Color Developer is common and does not affect test performance.

9. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{8,9}

10. MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Forceps, micropipette, calibrated loop, (7) Demineralized water.

Optional Materials:

(1) Microscope slide, (2) Micropipette, 5-10 µl, (3) Applicator sticks.

11. PROCEDURE

Test isolates should be 18-24 hours old and growing on nonselective culture media. Ensure that the culture is pure or that selected colonies are well-isolated from surrounding organisms.

Streptococci: Verify that the test isolate is catalase-negative, Gram-positive cocci with microscopic and colonial morphology consistent with enterococci or Group A streptococci.

Enteric Gram-negative bacilli: When screening for *E. coli*, test only lactose-fermenting, Gram-negative bacilli which are indole-positive, oxidase-negative, and nonhemolytic.⁷

- Using forceps, place a disk on a clean microscope slide or in the lid of an agar plate, free from excess moisture.
- Moisten the disk slightly with 5-10 µl of demineralized water using an adjustable micropipette or a 10 µl inoculating loop. Do not over-saturate the disk. Alternatively, the disk may be placed on the surface of the agar medium for rehydration.
- With an applicator stick or a sterile inoculating loop, remove a heavy visible paste of the test isolate and rub it into the disk.
- Incubate at room temperature (RT) for 2 minutes.
- Add one (1) drop of Color Developer to the disk.
- Allow up to one minute for a pink to red color to develop.

12. INTERPRETATION

Positive Test - Pink to red color development within one minute of adding Color Developer

Negative Test - Cream, yellow, or no color change within one minute of adding Color Developer

Notes:

- If the PYR disk has been pre-moistened on a blood agar plate, the background of the disk may be pale to straw yellow. This will not alter the appearance of a positive reaction.
- Weak results may occur if insufficient cells are tested or if the test isolate is poor in the expression of PYRase activity. Confirm a weak result by repeating the procedure with a larger inoculum and incubating the disk for 5 minutes, rather than 2 minutes, prior to the addition of the Color Developer.

13. EXPECTED VALUES⁶⁻¹⁰

<i>Enterococcus</i> spp.	+
<i>Streptococcus</i> spp.	
Group A (<i>S. pyogenes</i>)*	+
Groups B, C, F, G	-
Group D	-
Viridans group	-
<i>Escherichia coli</i>	-
<i>Citrobacter</i> and <i>Klebsiella</i> spp.	+

*Group A streptococci in the *anginosus* group are PYR negative.⁹

14. QUALITY CONTROL

All lot numbers of PYR Disk w/ Reagent have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL	INCUBATION	RESULTS
<i>Citrobacter freundii</i> ATCC® 8090	Ambient, 2 min. @ RT	Positive
<i>Enterococcus faecalis</i> ATCC® 29212	Ambient, 2 min. @ RT	Positive
<i>Escherichia coli</i> ATCC® 25922	Ambient, 2 min. @ RT	Negative
<i>Streptococcus agalactiae</i> ATCC® 12386	Ambient, 2 min. @ RT	Negative
<i>Streptococcus pyogenes</i> ATCC® 19615	Ambient, 2 min. @ RT	Positive

15. LIMITATIONS

- PYR Disk w/ Reagent provides a high probability for presumptive separation of enterococci and *S. pyogenes* (Group A) from other streptococci; however, it is only part of the overall scheme for identification. Further biochemical and/or serological testing may be required for definitive identification. Consult appropriate references for further instructions.^{7,9}
- Aerococcus viridans*, *Lactococcus garvieae*, certain *staphylococci*, and most *Corynebacterium haemolyticum* strains are PYR-positive, as are some *Enterobacteriaceae* and other Gram-negative bacilli. Consult appropriate references for differential tests when necessary.^{7,9}
- Some strains of enterococci are beta-hemolytic; therefore, colony morphology must be carefully evaluated to differentiate between enterococci and *S. pyogenes*.^{7,9}
- False-negative reactions may occur if inadequate inoculum is used or if inoculum is removed from selective media (e.g., Columbia CNA, PEA, MacConkey agars, etc.).

16. PERFORMANCE CHARACTERISTICS

PYR Disk w/ Reagent was evaluated in three separate studies using 284 strains of *Streptococcus* and *Enterococcus* spp.¹¹ The results are shown below.

Organism Tested	PYR Disk w/ Reagent	
	Positive*	Negative
<i>Enterococcus</i> spp.	49	0
<i>E. faecalis</i>	61	0
<i>E. faecium</i>	19	0
<i>E. durans</i>	16	0
<i>E. avium</i>	2	0
<i>S. pyogenes</i> (Group A)	34	1
<i>S. agalactiae</i> (Group B)	0	7
<i>S. bovis</i> (Group D)	0	34
<i>Streptococcus</i> spp. Group C	0	8
<i>Streptococcus</i> spp. Group F	0	9
<i>Streptococcus</i> spp. Group G	0	12
<i>Streptococcus</i> spp. viridans group	0	25
<i>S. pneumoniae</i>	0	7

* Five strains (3 *E. durans* and 2 *E. faecium*) required 5-minute incubations.

In a study by York *et al.* the PYR Disk w/ Reagent was used to validate the accuracy of rapid spot tests for identification of *E. coli*.⁶ Participating laboratories obtained 1,064 clinically significant strains of Enterobacteriaceae with core criteria of indole-positive, oxidase-negative, nonswarming colonies. Organisms were identified using commercial kit systems and an algorithm which included lactose fermentation, beta-hemolysis, PYR, and 4-methyl-umbelliferyl-β-D-glucuronide (MUG) results. The study demonstrated that nonhemolytic lactose-fermenters which are indole-positive, oxidase-negative, and PYR-negative, can be reliably identified as *E. coli* with greater than 99% accuracy, significantly reducing reagent costs, technologist time, and time to reporting.

17. BIBLIOGRAPHY

- Godsey, J., R. Schulman, and L.A. Eriquez. 1981. Abstract C-84. Abstracts of the 81st General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

- Facklam, R.R., L.G. Thacker, B. Fox, and L. Eriquez. 1982. J. Clin. Microbiol. 15:987-990.
- Bosley, G.S., R.R. Facklam, and D. Grossman. 1983. J. Clin. Microbiol. 18:1275-1277.
- Ellner, P.D., D.A. Williams, M.E. Hosmer, and M.A. Cohenford. 1985. J. Clin. Microbiol. 22:880-881.
- Morgan, J.W. 1987. Lab. Med. 18:682-683.
- York, M.K., E.J. Baron, J.E. Clarridge, R.B. Thomson, and M.P. Weinstein. 2000. J. Clin. Microbiol. 38:3394-3398.
- Clinical and Laboratory Standards Institute (CLSI). 2002. Abbreviated Identification of Bacteria and Yeast. Approved Guideline, M35-A. CLSI, Wayne, PA.
- Carroll, K C, M A Pfaller et al., 2019. Manual of Clinical Microbiology: Twelfth Edition (ASM Books).
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott’s Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
- Inoue, K., K. Miki, K. Tamura, and R. Sakazaki. 1996. J. Clin. Microbiol. 34:1811-1812.
- Data on File. Remel Inc., Lenexa, KS.

18. PACKAGING

REF R30854301 PYR Disk w/ Reagent.....50 Tests/Kit
REF R30854401 PYR Disk w/ Reagent.....100 Tests/Kit

19. SYMBOL LEGEND

REF	Catalogue Number
IVD	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Not for near patient testing
	Do not re-use
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
	Keep away from sunlight
	Importer
EC REP	Authorized representative in the European Community
UK CA	UK Conformity Assessed
	Manufacturer

 Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA

For technical information contact your local distributor.
 ATCC is a registered trademark of American Type Culture Collection.

www.thermofisher.com/microbiology
 Tel: (800) 255-6730 • International: (913) 888-0939
 www.oxid.com/IFU
 Europe +800 135 79 135 • US 1 855 2360 190
 CA 1 855 805 8539 • ROW +31 20 794 7071

Version	Date of modifications introduced
IFU30854301	March 2023 Addition of Importer symbol to Section 19

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PYR Color Developer	PYR Color Developer
----------------------------	---------------------

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- Údaje v záznamoch. Remel Inc., Lenexa, KS.