INTEGRAL SYSTEM GARDNERELLA

REF. 71724

- ***Sistema colorimetrico per l'identificazione e l'antibiogramma di Gardnerella vaginalis***
- ***Colorimetric system for identification and susceptibility testing of Gardnerella vaginalis***
- ***Système colorimétrique pour l'identification et l'antibiogramme de Gardnerella vaginalis***
- ***Sistema colorimétrico para la identificación y el antibiograma de Gardnerella vaginalis***
- ***Sistema colorimétrico para a identificação e o antibiograma de Gardnerella vaginalis***
- ***Σύστημα για τον προσδιορισμό και το αντιβιογράμμα Gardnerella vaginalis (Κολπίτιδα από Gardnerella)***

Code F00213
Rev. 2 / 22.03.2004

INTEGRAL SYSTEM GARDNERELLA

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INTEGRAL SYSTEM GARDNERELLA

Colorimetric system for identification and susceptibility testing of Gardnerella vaginalis

DESCRIPTION

INTEGRAL SYSTEM GARDNERELLA is an 18-wells system containing culture media with desiccated biochemical substrates for identification and susceptibility testing of *Gardnerella vaginalis*. The system is inoculated with the bacterial suspension and incubated at $36\text{ }^{\circ}\text{C} \pm 1\text{ }^{\circ}\text{C}$ for 18-24 hours. Tests for identification and susceptibility testing of *Gardnerella vaginalis* are interpreted by assessing the change in colour of the various wells.

CONTENT OF THE PACKAGES

The packages contain:

- 20 Systems INTEGRAL SYSTEM GARDNERELLA
- Ninhydrin Reagent (1 Bottle of 5.0 ml)
- 1 Dropper Bottle for Ninidrine Reagent
- 20 Vials of Physiological Solution (7.0 ml/vial)
- H_2O_2 Reagent (1 Bottle of 5.0 ml)
- Sheet of instructions

ITEMS NECESSARY BUT NOT INCLUDED IN THE PACKAGES

- Sundry microbiology laboratory materials
- Vaseline oil for microbiological use (Vaseline oil 2 Bottles of 50 ml)..... cod. 80278

CONFIGURATION

The configuration of the system is shown in Table n°1.

Table n°1

Well	BIOCHEMICAL IDENTIFICATION
1-IPP	Hyppurate hydrolysis
2-AMD	Starch hydrolysis
3-ESC	Aesculin hydrolysis
4-ONPG	ONPG Test
5-CAT	Catalase Test
6-UR	Urea Hydrolysis
7-GLU	Glucose fermentation
8-MAL	Maltose fermentation
9-SOR	Sorbitol fermentation
10-RIB	Ribose fermentation
11-DES	Destrine fermentation
12-INU	Inuline fermentation
Well	SUSCEPTIBILITY TESTING
13-MZ	Metronidazole - 100 µg/ml
14-C	Chloramphenicol - 8 µg/ml
15-CXM	Cefuroxime - 16 µg/ml
16-AMS	Ampicillin/Sulbactam - 8/4 µg/ml
17-SXT	Co-Trimoxazole - 30 µg/ml
18-C	Growth control

PRINCIPLE OF THE METHOD

INTEGRAL SYSTEM GARDNERELLA permits identification and susceptibility testing of *Gardnerella vaginalis* coming from clinical material (vaginal swabs, urethral secreted, seminal liquid).

- Identification is based on biochemical reactions carried out using media with specific substrates in wells from **1-IPP** to **12-INU**.
- Susceptibility testing is assessed according to the growth or inhibition of micro-organisms in the media containing antibiotic and growth indicator in wells from **13-MZ** to **17-SXT**.
- Well **18-C** does not contain antibiotics, but only the culture medium and the indicator and it serves as microbe growth control.

COMPOSITION

Table n°2

Well	Content
1-IPP	Cultural medium for showing of hyppurate hydrolysis
2-AMD	Cultural medium for showing of starch hydrolysis
3-ESC	Cultural medium for showing aesculin hydrolysis
4-ONPG	Cultural medium for showing orto-nitrophenylgalactopyranoside hydrolysis
5-CAT	Cultural medium for showing catalase reaction
6-UR	Cultural medium for showing urea hydrolysis
7-GLU	Cultural medium for showing glucose fermentation
8-MAL	Cultural medium for showing maltose fermentation
9-SOR	Cultural medium for showing sorbitole fermentation
10-RIB	Cultural medium for showing rabinose fermentation
11-DES	Cultural medium for showing destrine fermentation
12-INU	Cultural medium for showing inuline fermentation
13-MZ	Medium with Metronidazole - 100 µg/ml
14-C	Medium with Chloramphenicol - 8 µg/ml
15-CXM	Medium with Cefuroxime - 16 µg/ml
16-AMS	Medium with Ampicillina/Sulbactam - 8/4 µg/ml
17-SXT	Medium with Co-Trimoxazole - 30 µg/ml
18-C	Medium without antibiotics

Physiological Solution (g/l) : Sodium chloride 9 g ; Distilled water 1000 ml ; pH 6,8 ± 0,2
Ninhydrin Reagent (g/l) : Ninhydrin 7 g ; Methoxyethanol 100 ml
H₂O₂ Reagent (g/l) : Hydrogen Peroxide 38 ml ; Distilled water 62 ml

COLLECTION AND CONSERVATION OF THE SAMPLES

The colonies to be subjected to the identification test and susceptibility testing with INTEGRAL SYSTEM GARDNERELLA must be collected from a cultural medium, selective and not selective, used for isolation of *Gardnerella vaginalis*.

TEST PROCEDURE

1. Take a system from the kit.
2. After checking with apposite researches that colonies developed on cultural media presumably belong to *Gardnerella* genus, take one or more colonies morphologically similar, well isolated from a culture media for *Gardnerella vaginalis* research and suspend into vial of Physiological Solution contained in the kit till to obtain a turbidity equivalent to 0,5-1 Mac Farland.
3. Transfer:
 - 0.2 ml of the bacterial suspension to each well of the system.
 - 0.1 ml of the bacterial suspension to well **1-IPP**.
4. Cover the well **6-UR** with 2 drops of sterile vaseline oil for microbiological use.
5. Cover the system with the apposite lid and incubate at $36\text{ }^{\circ}\text{C} \pm 1\text{ }^{\circ}\text{C}$ for 18-24 hours.

INTERPRETATION OF THE RESULTS

At the end of the incubation:

- Add 2 drops of Ninhydrin Reagent for Hippurate test to well **1-IPP**, after transferring the content of the amber dropper bottle which has to be stored in the dark freezing.
- Add 2 drops of H₂O₂ Reagent to well **5-CAT**.
- Observe the changes in colour of the wells from **1-IPP** to **12-INU** and interpret the results using Table n°3 (BIOCHEMICAL IDENTIFICATION).
Use Reactions scheme (Table n°4) to confirm identification of *Gardnerella vaginalis*.
- Observe the changes in colour of the wells from **13-MZ** to **17-SXT** and interpret the results as indicated in table n°3 (SUSCEPTIBILITY TESTING).
- Control well (**18-C**) must result positive (yellow). In case it results negative (red or orange), it is necessary verify inoculum vitality, its correct preparation and repeat the test using a new system.

Table n°3

Well	BIOCHEMICAL IDENTIFICATION	Well colour	
		Positive reaction	Negative reaction
1-IPP	Hyppurate hydrolysis	Blue-violet	Light blue
2-AMD	Starch hydrolysis	Yellow-green	Violet
3-ESC	Aesculin hydrolysis	Black	Yellow
4-ONPG	ONPG test	Yellow	Colorless
5-CAT	Catalase test	Presence of Bubbles	Absence of bubbles
6-UR	Urea hydrolysis	Red-fuchsia	Yellow
7-GLU	Glucose fermentation	Yellow	Red-orange
8-MAL	Maltose fermentation	Yellow	Red-orange
9-SOR	Sorbitol fermentation	Yellow	Red-orange
10-RIB	Ribose fermentation	Yellow	Red-orange
11-DES	Destrine fermentation	Yellow	Red-orange
12-INU	Inuline fermentation	Yellow	Red-orange

SUSCEPTIBILITY TESTING		
WELL COLOUR	BACTERIAL GROWTH	INTERPRETATION
Red	Inhibited	S = Sensitive
Orange	Intermediate	I = Intermediate sensitivity
Yellow	Good	R = Resistant

Table n°4

Microorganism	1-IPP	2-AMD	3-ESC	4-ONPG	5-CAT	6-UR	7-GLU	8-MAL	9-SOR	10-RIB	11-DES	12-INU
<i>Gardnerella vaginalis</i>	+	+	-	V	-	-	+	+	-	+	+	-

+ = Positive reaction

V = Variable reaction

- = Negative reaction

NB: Different biochemical reactions do not identify *Gardnerella vaginalis*

QUALITY CONTROL

Each lot of INTEGRAL SYSTEM GARDNERELLA is subjected to quality control using a suspension in Physiological Solution of strain *Gardnerella vaginalis* ATCC 14018 at concentrations equivalent to 0,5-1 Mac Farland. As negative control is used Physiological Solution.

FACTORS THAT MAY INVALIDATE THE TEST

Poor standardisation of the inoculum; application of the method to micro-organisms not in the genus *Gardnerella*; mixed or contaminated cultures; use of expired systems or expired reagents; incorrect application of the technique.

LIMITS AND WARNINGS

In some cases, for the final identification of *Gardnerella vaginalis* further confirmation tests are necessary.

PERFORMANCE

The bacterial identification results obtained with the INTEGRAL SYSTEM GARDNERELLA agree with those obtained using traditional test-tube methods⁽¹⁾.

The susceptibility testing results obtained with the INTEGRAL SYSTEM GARDNERELLA agree with those obtained using the radial diffusion methods according to Bauer et al.⁽²⁾ recommended by the Food and Drug Administration (FDA) ⁽¹⁾.

PRECAUTIONS

The product INTEGRAL SYSTEM GARDNERELLA cannot be classified as hazardous under current legislation, nor does it contain harmful substances in concentrations $\geq 1\%$. It therefore does not require a Safety Data Sheet to be available. INTEGRAL SYSTEM GARDNERELLA is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained operators, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes of temperature. In such conditions the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, INTEGRAL SYSTEM GARDNERELLA and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Code	Package
INTEGRAL SYSTEM GARDNERELLA	71724	20 tests

TABLE OF SYMBOLS

<u>SYMBOL</u>	<u>MEANING</u>	<u>SYMBOL</u>	<u>MEANING</u>
	In Vitro Diagnostic Medical Device		Do not reuse
	Manufacturer		Contains sufficient for <n> tests
	Catalogue number		Fragile, handle with care
	Use by		Caution, consult accompanying documents
	Temperature limitation		Batch code

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