AnaeroJar™ 2.5L, AnaeroBox™ 2.5L & AnaeroBox™ 3.5L

EN

[REF] AG0025A, AB0025A, & AB0035A

Consult the instructions for use of associated *Atmosphere Generation Systems (AGS)* products (AnaeroGen, CO₂Gen and CampyGen) for procedural details.

Intended Use

The AnaeroJar™ 2.5L, AnaeroBox™ 2.5L and AnaeroBox™ 3.5L (AG0025A, AB0025A, and AB0035A) are used to maintain optimal atmospheric conditions (e.g. microaerophilic, enriched CO₂ or anaerobic) for the culture of microorganisms. The accessories are used in a diagnostic workflow to aid clinicians in the growth, isolation and/or detection of bacteria grown on culture media devices.

The accessories are for professional use only, are not automated and nor are they companion diagnostics.

Summary and Explanation

These products are used to maintain the optimal atmospheric conditions for culture of microorganisms.

Principle of Method

AnaeroJar™ 2.5L (AG0025A), AnaeroBox™ 2.5L (AB0025A) and AnaeroBox™ 3.5L (AB0035A) are accessories which are intended to be used to maintain optimal atmospheric conditions for cultured microorganisms e.g., microaerophilic, enriched CO₂ and anaerobic conditions.

Description of Components

AnaeroJar 2.5L (AG0025A) is supplied as a 2.5L jar for holding up to 12 petri dishes, with 2.5L capacity polycarbonate jar base and lid secured to base by four self-venting clips. The Jar is designed to be used with AnaeroGen 2.5L Sachet (AN0025A), CO₂Gen™ 2.5L Sachet (CD0025A) and CampyGen™ 2.5L Sachet (CN0025A). The following products are individual components of AnaeroJar 2.5L (AG0025A), which are available separately.

- AnaeroJar Base (AG0026A)
- AnaeroJar Lid (AG0027A)
- AnaeroJar Handle (AG0028A)
- AnaeroJar Plate Carrier (AG0029A)
- AnaeroJar O-Ring (AG0030A)
- AnaeroJar Clips (AG0031A)

AnaeroBox™ 2.5L (AB0025A) is supplied as a 2.5L rectangular container for holding up to 12 petri dishes or 6 microplates, with a polycarbonate lid and a base with polyacetal latches. The Box is designed to be used with AnaeroGen 2.5L Sachet (AN0025A), CO₂Gen™ 2.5L Sachet (CD0025A) and CampyGen™ 2.5L Sachet (CN0025A). The following products are individual components of Anaerobox 2.5L (AB0025A), which are available separately.

- Anaerobox Lid (AB0025L)
- Anaerobox Seal (AB0025S)

AnaeroBox™ 3.5L (AB0035A) is supplied as a 3.5L rectangular container for holding up to 18 petri dishes, with a polycarbonate lid and a base with polyacetal latches. The Box is designed to be used with AnaeroGen 3.5L Sachet (AN0035A) and CampyGen™ 3.5L Sachet (CN0035A). The following products are individual components of AnaeroBox™ 3.5L (AB0035A), which are available separately.

- Anaerobox Lid (AB0035L)
- Anaerobox Seal (AB0035S)

Materials Required but Not Supplied

- Inoculated media
- Incubators
- Quality control organisms
- AnaeroGen 2.5L Sachet (AN0025A)
- AnaeroGen 3.5L Sachet (AN0035A)
- CO₂Gen™ 2.5L Sachet (CD0025A)
- CampyGen™ 2.5L Sachet (CN0025A)
- CampyGen™ 3.5L Sachet (CN0035A)
- Resazurin Anaerobic Indicator (BR0055B)

Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or product.

- Do not use the device if there are signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- As soon as the AGS paper sachet is exposed to the air, the reaction will start. It is therefore essential that the paper sachet is placed in the jar/box and the jar/box sealed within one minute.
- The reaction of the ascorbic acid with oxygen is exothermic. However, the temperature of the AGS paper sachet will not exceed 65°C.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Procedure

1. Place inoculated plates into the plate carrier. Disposable plastic Petri dishes should be of the vented variety to aid gas transfer between interior and exterior of the dishes.
2. When using the anaerobic system (i.e. AN0025A/AN0035A), prepare the Oxoid Anaerobic Indicator (BR0055B) by cutting and exposing 10mm of the fabric strip, insert into the smaller, upper clip on the dish carrier.
3. Lower the carrier into the polycarbonate base.
4. Tear open an AnaeroGen/CampyGen/CO₂Gen sachet at the tear-nick indicated and remove the paper sachet from within.
5. Immediately place the paper sachet in the appropriate clip in the plate carrier within the jar / box (see Technical insert).
6. Having inserted the sachet into the carrier immediately place the lid on the jar/ box, making sure the 'O' ring/seal is in place.
7. When using the anaerobic system the anaerobic indicator will change from pink to white giving a visual indication of anaerobiosis.
8. Remove jar/ box after the appropriate incubation period and open lid by carefully depressing the clips to release the jar/box lid from the base. Excessive downward pressure on the clips/seal should be avoided.
9. Occasionally, a slight vacuum may occur, producing a negative pressure and resulting in resistance to the removal of the lid (after release of the clips). This is overcome by using an appropriate object such as a small coin to turn the vacuum relief screw anticlockwise allowing an inflow of air. It is important, however, to ensure the valve is resealed, by turning clockwise, prior to further use.

Limitations

The jars/ boxes must be free from damage to achieve their intended purpose. See associated AGS product information for other limitations.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Consult instructions for use or consult electronic instructions for use
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier

	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in Japan	Made in Japan

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, England

For technical assistance please contact your local distributor.

Revision Information

Version	Date of issue and modifications introduced
3.0	2024-07-25

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Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.