## Sheep Blood Agar Base

REF **CM0854B**

EN

### Intended Use

Sheep Blood Agar Base (CM0854B) is a nonselective medium for the isolation and cultivation of fastidious pathogens and other microorganisms with clearly visible haemolysis from clinical samples (e.g., wounds, throat, genital, nose, groins etc.).

Sheep Blood Agar Base (CM0854B) is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The device is for professional use only, is not automated and is not a companion diagnostic.

### Summary and Explanation

In general *Streptococci* are grown in blood agar and the presumptive identification of the species is based on their various degrees of haemolysis in blood agar media. The addition of blood to culture media provides means of identification, as the ability to haemolyse blood can be used to classify certain microorganisms<sup>2</sup>. Any incomplete or partial lysis of red blood cells leading to a formation of green or brown colour around the colonies is known as alpha haemolysis. A beta haemolysis is characterised by a complete clearing of the medium around the colonies, and this is due to the complete lysis of red blood cells<sup>1,2</sup>. Consequently, blood agar are recommended for the primary culture of several human pathogens, including Gram-positive and Gram-negative microorganisms<sup>1,3</sup>.

The addition of blood to culture media provides means of identification, as the ability to haemolyse blood can be used to classify certain microorganisms<sup>1,3</sup>. Consequently, it is important to isolate *Haemophilus* species from clinical samples as early diagnosis is essential to prevent the morbidity and mortality associated with this bacterium in the vulnerable population<sup>4</sup>. Even after the introduction of pneumococcal vaccine as part of childhood and adult vaccination schedules there is considerable morbidity and mortality associated with *S. pneumoniae* infections<sup>5</sup>. Given the greater risk to patient outcomes and onwards transmission posed by these pathogenic Gram positive and Gram-negative species, it is a vital part of the diagnostic workflow to distinguish the clinically relevant strains like *Streptococcus*, *Staphylococcus* and *Haemophilus* species that can cause fatal infections especially in younger population.

### Principle of Method

The Sheep Blood Agar Base is a non-selective, general- purpose medium that facilitates the growth of a wide variety of microorganisms, including fastidious organisms, from clinical specimens in order to differentiate bacteria based on their haemolytic properties. Growth of the microorganisms is achieved through the inclusion of tryptone and peptone which provide nitrogen, carbon, and amino acids, and yeast extract which provide growth promoting B vitamins. Sodium chloride maintains the osmotic balance. Agar is the solidifying agent. The Sheep Blood Agar Base is formulated to be compatible with sheep blood. Defibrinated Sheep Blood is supplemented to the medium to provide additional growth factors and to visualise haemolytic reactions. On blood agar plates, beta haemolytic reactions can be visualised in the form a clear, colourless zone where the red blood cells have been lysed and haemoglobin destroyed to give a colourless zone. Alpha haemolytic reactions occur when bacteria reduce haemoglobin to methaemoglobin which produces a greenish zone around the colonies. In comparison to other blood agar bases supplemented with blood, Sheep Blood Agar Base showed considerable improvement and the expected beta haemolytic reactions with fastidious *Streptococci*.

### Typical Formula

	grams per litre
Tryptone	14.0
Peptone Neutralised	4.5
Yeast extract	4.5
Sodium chloride	5.0
Agar	12.5

### Materials Provided

CM0854B – 500g of dehydrated Sheep Blood Agar Base powder that yields approximately 12.5L after reconstitution.

### Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish
- 7% Defibrinated Sheep Blood

### Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.

- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 8°C.

### Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product ([www.thermofisher.com](http://www.thermofisher.com)).

### Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

### Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 7, ID 4, ID 12, and Q5.

### Procedure

Suspend 40g in 1 litre of distilled water and bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C and aseptically add 7% sterile sheep blood. Pour into sterile petri dishes and allow to set. Streak sample onto set plates and incubate for 18 - 24 h @ 37° ± 2°C aerobic.

### Interpretation

Any incomplete or partial lysis of red blood cells leading to a formation of green or brown colour around the colonies is known as alpha haemolysis. A beta haemolysis is characterised by a complete clearing of the medium around the colonies, and this is due to the complete lysis of red blood cells.

### Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 @ 37°C aerobic

Inoculum level: 10 – 100 cfu	
<i>Streptococcus pneumoniae</i> ATCC®6303™	Pinpoint – 1mm colourless colonies
<i>Streptococcus viridans</i> NCTC 1080	0.25 – 0.5 mm colourless colonies
<i>Streptococcus pyogenes</i> ATCC® 19615™	0.25 – 1 mm colourless/white colonies
<i>Staphylococcus aureus</i> ATCC® 25923™	1 – 2 mm white/straw colonies
<i>Pseudomonas aeruginosa</i> ATCC® 27853™	2 – 3 mm straw colonies, green pigmentation
Enriched with 7% v/v sheep blood	

<i>Streptococcus pneumoniae</i> ATCC®6303™	0.5 – 2 mm colourless colonies, α haemolysis
<i>Streptococcus viridans</i> NCTC 1080	0.5 – 1 mm grey/green colonies, α haemolysis
<i>Haemophilus influenzae</i> ATCC® 6303™	Pinpoint – 0.25 mm white/grey colonies

Zones of growth/no growth surrounding X, V and X+V factor discs (DD0003, DD0004 and DD0005) when plain plates are inoculated with the following organisms and incubated at 37°C for 18 hours:

	X	V	X+V
<i>Haemophilus influenzae</i> ATCC® 9344™	0	0	≥ 15 mm
<i>Haemophilus influenzae</i> ATCC® 19418™	0	0	≥ 15 mm
<i>Haemophilus influenzae</i> ATCC® 49247™	0	0	≥ 15 mm
<i>Haemophilus parainfluenzae</i> ATCC® 33392™	0	≥ 20 mm	≥ 20 mm

Zones of inhibition with Bacitracin discs (DD0002) shall be 10-20mm when 7% v/v horse blood plates are inoculated with *Streptococcus pyogenes* ATCC®19615 and incubated at 37°C for 18 hours.

<b>Testing performed in accordance with current CLSI M22 A</b>	
Reactions after incubation at 35°C for 18-24 hours Enriched with 7% v/v sheep blood Colony count is ≥ 70% of the control medium count	
<i>Streptococcus pneumoniae</i> ATCC® 6303™	0.5-2mm colourless colonies. A haemolysis
<i>Streptococcus pyogenes</i> ATCC® 19615™	0.5-1mm colourless / white colonies β haemolysis
<i>Staphylococcus aureus</i> ATCC® 25923™	1-2mm white, straw colonies
<i>Escherichia coli</i> * ATCC®29222™	1 – 2 mm cream colonies

### Limitations

As with all media, atypical organisms may give anomalous reactions. A small number of atypical strains may give a weak reaction and fail to grow, especially when low numbers are present in the sample. Further biochemical and serological tests must be carried out for further identification.

### Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of microorganism is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device, which must meet the defined acceptance criteria. The precision of Sheep Blood Agar Base (CM0854) was demonstrated by an overall pass rate of 100% over 12 months of testing (12.07.2021 – 09.08.2022).

Sheep Blood Agar Base (CM0854B) devices is tested in-house as part of the QC process since the products were launched in 1997 and 1998. For target organisms, when using optimum inoculum levels and incubating at the recommended time and temperature conditions, the user can recover organisms with colony size and morphology as listed in this document.

### Bibliography

- Public Health England. 2021a. 'Identification of *Streptococcus* Species, *Enterococcus* Species and Morphologically Similar Organisms'. UK Standards for Microbiology Investigations.ID 4 (3).  
<https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
- Public Health England. 2021b. "Identification of *Haemophilus* species and the HACEK group of organisms".

2021. ID 12 (4). <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>
- Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations.ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
  - European Centre for Disease Prevention and Control 2014. 'Invasive Haemophilus influenzae disease' <https://www.ecdc.europa.eu/en/invasive-haemophilus-influenzae-disease/facts>
  - Pascual, E. A., Robinot, D. C., Herrero, C. G., Irujo, M. N., Ponferrada, M. R., & Vilalta, M. P. (2021). Pediatric chest X-rays during the COVID- 19 pandemic. Radiología (English Edition), 63(2),106-114.

## Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Consult instructions for use
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

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### Literatura

- Public Health England. 2021a. 'Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms'. UK Standards for Microbiology Investigations.ID 4 (3). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
- Public Health England. 2021b. "Identification of Haemophilus species and the HACEK group of organisms". 2021. ID 12 (4). <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>
- Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations.ID 7 (4). <https://www.gov.uk/government/publications/smi->

- [id-7-identification-of-staphylococcus-species- micrococcus-species-and-rothia-species](#).
- European Centre for Disease Prevention and Control 2014. 'Invasive Haemophilus influenzae disease' <https://www.ecdc.europa.eu/en/invasive-haemophilus-influenzae-disease/facts>
  - Pascual, E. A., Robinot, D. C., Herrero, C. G., Irujo, M. N., Ponferrada, M. R., & Vilalta, M. P. (2021). Pediatric chest X-rays during the COVID- 19 pandemic. Radiología (English Edition), 63(2), 106-114.

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## Bibliographie

1. Public Health England. 2021a. 'Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms'. UK Standards for Microbiology Investigations.ID 4 (3). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
2. Public Health England. 2021b. "Identification of Haemophilus species and the HACEK group of organisms". 2021. ID 12 (4). <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>
3. Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations.ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. European Centre for Disease Prevention and Control 2014.' Invasive Haemophilus influenzae disease' <https://www.ecdc.europa.eu/en/invasive-haemophilus-influenzae-disease/facts>
5. Pascual, E. A., Robinot, D. C., Herrero, C. G., Irujo, M. N., Ponferrada, M. R., & Vilalta, M. P. (2021). Pediatric chest X-rays during the COVID- 19 pandemic. Radiología (English Edition), 63(2),106-114.

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[www.thermofisher.com](http://www.thermofisher.com)

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
## Bibliographie

1. Public Health England. 2021a. 'Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms'. UK Standards for Microbiology Investigations.ID 4 (3). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
2. Public Health England. 2021b. "Identification of Haemophilus species and the HACEK group of organisms". 2021. ID 12 (4). <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>
3. Public Health England 2020b. 'Identification of Staphylococcus Species, Micrococcus Species and Rothia Species'. UK Standards for Microbiology Investigations.ID 7 (4). <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
4. European Centre for Disease Prevention and Control 2014.' Invasive Haemophilus influenzae disease' <https://www.ecdc.europa.eu/en/invasive-haemophilus-influenzae-disease/facts>
5. Pascual, E. A., Robinot, D. C., Herrero, C. G., Irujo, M. N., Ponferrada, M. R., & Vilalta, M. P. (2021). Pediatric chest X-rays during the COVID- 19 pandemic. Radiologia (English Edition), 63(2),106-114.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Royaume-Uni

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Toate celelalte mărci comerciale aparțin Thermo Fisher Scientific Inc. și subsidiarelor acesteia.