www.thermofisher.com

Bile Aesculin Agar

REF CM0888K & CM0888B

EN

Intended Use **IVD**

Bile Aesculin Agar is a differential for the isolation and presumptive identification of enterococci and Group D streptococci from clinical samples. The device is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections. The device is for professional use only, is not automated and nor is it a companion diagnostic.

Summary and Explanation

The major use of Bile Aesculin Agar is to differentiate between enterococci / Group D streptococci and non-Group D streptococci. It may also be used for the presumptive identification of other groups of microorganisms.

The value of bile tolerance together with hydrolysis of aesculin^{1,2,3,4,5} as a means of presumptively identifying enterococci/Group D streptococci is widely recognised. The use of these parameters forms the basis of Bile Aesculin Agar which was described by Swan⁶ who concluded that the use of this medium is a valid alternative to Lancefield grouping for the recognition of enterococci/Group D streptococci. Facklam⁷ further confirmed its usefulness in differentiating enterococci/Group D streptococci from non-Group D streptococci while other workers^{8,9,10} have used the medium for presumptive identification of the *Klebsiella-Enterobacter- Serratia* group amongst the *Enterobacteriaceae*. Bile aesculin agar is also used to identify enterococci from water samples¹¹.

Principle of Method

Peptone supplies the nitrogen source and other nutrients to support the growth of the target bacteria. Agar is added as a solidifying agent. Enterococci / Group D streptococci hydrolyse aesculin to form aesculetin and dextrose. Aesculetin combines with the ferric citrate in the medium to form a dark brown or black complex which is indicative of a positive result. Bile salts inhibit Gram-positive bacteria other than enterococci / Group D streptococci.

Typical Formula

	Grams per litre
Peptone	8.0
Bile Salts	20.0
Ferric citrate	0.5
Aesculin	1.0
Agar	15.0

Materials Provided

CM0888K – 25Kg dehydrated culture media

CM0888B – 500g dehydrated culture media

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dishes

Storage

- Store product in its original packaging between 10°C and 30°C.
 - Keep container tightly closed.
 - The product may be used until the expiry date stated on the label.
 - Protect from moisture.
 - Store away from light.
 - Allow reconstituted product to equilibrate to room temperature before use.
- Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical

advice/attention.

- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 29.

Procedure

Refer to local protocols and guidelines.

Suspend 44.5g in 1 litre of distilled water and bring gently to the boil to dissolve completely. Sterilize by autoclaving at 121°C for 15 minutes. Allow to cool to 45°C, mix well and pour into sterile Petri dishes.

Using a sterile loop inoculate the medium with 4 to 5 colonies of the test strain and incubate at 37°C for 18 hours to 24 hours. The result is positive for bile salt tolerance and aesculin hydrolysis if blackening of the medium occurs.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 37°C for 18 - 24 hours

Positive Controls Inoculum of 10 - 100 colony forming units (cfu). Colony count is $\geq 50\%$ of the control medium count	
<i>Enterococcus faecalis</i> ATCC®19433™	0.25-1mm, brown colonies, aesculin hydrolysis
Inoculum level $10^4 - 10^6$	
<i>Enterococcus faecalis</i> ATCC®29212™	0.25-1mm, brown colonies, aesculin hydrolysis
<i>Enterococcus faecium</i> NCTC 7171	0.25-1mm, brown colonies, aesculin hydrolysis
<i>Enterococcus durans</i> ATCC®19432™	0.25-1mm, brown colonies, aesculin hydrolysis
<i>Streptococcus bovis</i> ATCC®27960™	0.25-1mm, brown colonies, aesculin hydrolysis
<i>Enterobacter aerogenes</i> ATCC®13048™	1-2mm, brown colonies, aesculin hydrolysis
Negative Control Inoculum level $10^4 - 10^6$	
<i>Streptococcus pyogenes</i> ATCC®19615™	No growth

Limitations

This test cannot be used in isolation to confirm the identification of enterococci / Group D streptococci, as other bacteria are also able to hydrolyse aesculin in the presence of bile but it can be used in combination with other tests such as salt tolerance.

Performance

Performance was evaluated using 24 bacterial strains including the following: *E. faecalis*, *E. faecium*, *Streptococcus pyogenes*,

E.aerogenes, *Klebsiella pneumoniae*, *Staphylococcus aureus* and *Escherichia coli*. All organisms gave expected growth characteristics according to the current product specification.

Bibliography

- (1) Facklam R. R. and Moody M. D. (1970). *Appl. Microbiol.* 20, 245-250.
- (2) Isenberg H. D., Goldberg D. and Sampson J. (1970). *Appl. Microbiol.* 20, 433-436.
- (3) Sabbaj J., Sutter V. L. and Finegold S. M. (1971). *Appl. Microbiol.* 22, 1008-1011.
- (4) Facklam R. (1972). *Appl. Microbiol.* 23, 1131-1139.
- (5) Facklam R. *et al* (1974). *Appl. Microbiol.* 27, 107-113.
- (6) Swan A. (1954). *J. Clin. Path.* 7, 160-163.
- (7) Facklam R. (1973). *Appl. Microbiol.* 26, 138-145.
- (8) Wasilaukas B. L. (1971). *Appl. Microbiol.* 21, 162-163.
- (9) Lindell S. S. and Quinn P. (1975). *J. Clin. Microbiol.* 1, 440-443.
- (10) Chan P. C. K. and Porschen R. K. (1977). *J. Clin. Microbiol.* 6, 528- 529.
- (11) Standing Committee of Analysts, *The Microbiology of Drinking Water (2012) - Part 5 - Isolation and enumeration of Enterococci by membrane filtration*, Methods for the Examination of Waters and Associated Materials, in this series, Environment Agency.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier

	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, England

For technical assistance please contact your local distributor.

Version	Modifications Introduced
3.0	2024-09-12

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.