www.thermofisher.com

Brilliance™ UTI Clarity

EN

REF CM1106B, CM1106W & CM1106T

Intended Use^{IVD}

Brilliance™ UTI Clarity (CM1106B, CM1106W and CM1106T) devices are differential media for the isolation and presumptive differentiation of the common microorganisms causing urinary tract infections (UTIs) and direct identification of *E. coli* from clinical samples, including urine samples.

Brilliance™ UTI Clarity devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having UTIs.

The devices are for professional use only, are not automated and nor are they companion diagnostics.

Summary and Explanation

UTIs occur due to the presence and growth of microorganisms in the urinary tract, giving rise to a range of clinical presentations¹. Prominent human pathogens are often associated with UTIs, with the most common of which being *Escherichia coli*, accounting for roughly 77% of isolates¹. Other microorganisms associated with acute, uncomplicated UTIs include *Proteus mirabilis*, *Staphylococcus saprophyticus*, other coagulase negative staphylococci, streptococci, and *Candida* species¹. Other microorganisms are associated with complicated UTIs, including *Klebsiella*, *Enterobacter* and *Proteus* spp., *Enterococcus* spp. and *Pseudomonas aeruginosa*¹. Although *S. aureus* rarely causes infection, it is associated with renal abnormality, secondary infections to bacteraemia, surgery, and catheterisation¹. Other coagulase negative Staphylococci have been associated with complicated infections in individuals with abnormalities of the urinary tract or underlying disease¹.

Principle of Method

Differentiation of the common microorganisms causing urinary tract infections (UTIs) is achieved through the inclusion of two chromogens that are targeted by specific enzymes: β-galactosidase and β-glucosidase. The action of these enzymes on the chromogens causes release of the coloured component inside the bacterial cell, resulting in coloured colonies. The colour produced depends on which enzymes the organisms produce. β-galactosidase is produced by *E. coli*, which grow as pink colonies. β-glucosidase is produced by enterococci which form blue or turquoise colonies. Cleavage of both chromogens by the KESC (*Klebsiella*, *Enterobacter*, *Serratia*, *Citrobacter*) group of coliforms gives dark blue or purple colonies.

Tryptophan acts as an indicator of tryptophan deaminase activity (TDA), resulting in colonies of *Proteus*, *Morganella* and *Providencia* spp. with brown halos. Peptone forms the nutritional source within the medium.

Typical Formula

	<u>grams per litre</u>
Peptone	9.0
Chromogenic mix	17.0
Tryptophan	1.0
Agar	10.0

Materials Provided

- CM1106B:** 500g of Brilliance™ UTI Clarity (Dehydrated)
- CM1106W:** 1.85kg of Brilliance™ UTI Clarity (Dehydrated)
- CM1106T:** 5kg of Brilliance™ UTI Clarity (Dehydrated)

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dishes

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the

relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 01, B 41 and Q 5.

Procedure

Suspend 37.0g in 1 litre of distilled water and mix well. Sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C. Mix well and dispense into sterile containers.

Interpretation

Once reconstituted; *Escherichia coli* will show as pink colonies, *Enterobacter aerogenes* will show as blue/purple colonies, *Proteus mirabilis* as straw/brown colonies, *Enterococcus faecalis* as blue/turquoise colonies, *Staphylococcus aureus* as cream colonies, *Citrobacter freundii* as 1-3mm blue/purple colonies and *Staphylococcus epidermidis* as white colonies

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18h @ 37°C

Positive Controls	
Inoculum level 10-100 cfu Colony count is ≥ 70% of the control medium count.	
<i>Escherichia coli</i> ATCC® 25922™	1-2mm, pink colonies.
<i>Enterobacter aerogenes</i> ATCC® 13048™	1-2mm, blue/ purple colonies
<i>Proteus mirabilis</i> NCTC 10975	1-5mm, straw/brown colonies
<i>Enterococcus faecalis</i> ATCC® 29212™	Pinpoint -1mm, blue/ turquoise colonies
<i>Staphylococcus aureus</i> ATCC® 25923™	0.5-2mm, cream colonies
<i>Citrobacter freundii</i> NCTC 8581	1-3mm, blue/purple colonies
<i>Staphylococcus epidermidis</i> ATCC® 14990™	0.25-1mm, white colonies

Indole may be detected by removing a few colonies, spreading onto filter paper and adding 1-2 drops of DMACA Indole reagent (dimethylamino cinnamaldehyde). *Escherichia coli* should be positive (blue/green) and *Enterobacter aerogenes* negative (colourless to pink).

Limitations

Organisms with atypical enzyme patterns may give anomalous reactions on Brilliance UTI Clarity.

S. saprophyticus gives similar colouration to *E. coli* but the colony morphology and a gram smear should give clear differentiation. Some streptococci can produce very pale blue colonies which should not be confused with enterococci. Some strains of *Shigella* can produce variable colony colours (i.e. on the same plate) and should be confirmed.

It is not mandatory, but an indole test can be performed using DMACA indole to confirm *E. coli* (do not use Kovac's reagent as the colour of the *E. coli* colonies may interfere with the red colour of a positive indole test). The reagent should not be applied directly to the plate, but the test should be performed on filter paper. This test will distinguish between *E. coli* and *Enterobacter*, and between *Proteus mirabilis* and other species.

Due to the colour of the colonies, care should be taken when using colonies directly to carry out identifications with kits that depend on a colour change. Customers performing susceptibility tests on colonies from this medium should be aware that a small number of antimicrobial/organism combinations may give results outside of the usual tolerances. It is preferable to sub-culture onto a non-selective medium, such as Tryptone Soya Agar or Columbia agar.

The medium should be always protected from light except during inoculation and after incubation.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of common microorganisms causing UTIs is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices, which must meet the defined acceptance criteria. The precision of Brilliance™ UTI Clarity Agar was demonstrated by an overall pass rate of 100% obtained for the product over 5 months of testing (November 2021 – April 2022; 10 batches).

Brilliance™ UTI Clarity are tested in-house as part of the QC process since the products were launched in 2006. For target organisms, when using 10-100 cfu inoculum of *Escherichia coli*, *Enterobacter aerogenes*, *Proteus mirabilis*, *Enterococcus faecalis*, *Staphylococcus aureus*, *Citrobacter freundii* and *Staphylococcus epidermidis* and incubating the device at 37°C for 18 hours, the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

- Public Health England. 2019. 'Investigation of urine'. UK Standards for Microbiology Investigations. B 41. Issue 8.7. <https://www.gov.uk/government/publications/smi-b-41-investigation-of-urine>

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date

	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2023-05-25

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc.