# remel

## BactiDrop™ Oxidase EN

REF R21540

### INTENDED USE

BactiDrop Oxidase is a reagent for use in qualitative procedures for the detection of the cytochrome oxidase enzyme activity in microorganisms. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

### SUMMARY AND EXPLANATION

The oxidase test was originally devised to identify *Neisseria* species, but later was used to separate the *Pseudomonadaceae* from the oxidase-negative *Enterobacteriaceae*.<sup>1,2</sup> Most gram-positive bacteria are oxidase-negative and many of the gram-negative bacilli, other than the *Enterobacteriaceae*, are variable.<sup>3</sup> This formulation is based on the oxidase reagent formula of Kovacs.<sup>1</sup>

### PRINCIPLE

The oxidase test is based on the bacterial production of an intracellular oxidase enzyme. The oxidase reaction is due to the presence of a cytochrome oxidase system, which activates the oxidation of reduced cytochrome by molecular oxygen, which in turn acts as an electron acceptor in the terminal stage of the electron transfer system. Microorganisms which produce oxidase enzyme, when in the presence of atmospheric oxygen, cytochrome c, and phenylenediamine oxidase reagent, oxidize to form a purple colored compound, indophenol. Ascorbic acid is incorporated in the reagent as a reducing agent to reduce auto-oxidation and to improve stability.<sup>3</sup>

### REAGENTS (CLASSICAL FORMULA)\*

N,N,N,N-tetramethyl-1,4-phenylenediamine  
(CAS 100-22-1) ..... 10.0 g  
Ascorbic Acid (CAS 50-81-7) ..... 2.0 g  
Demineralized Water (CAS 7732-18-5)..... 1000.0 ml

\*Adjusted as required to meet performance standards.

### PRECAUTIONS

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully. Refer to Safety Data Sheet on company website for detailed information on reagent chemicals.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction, do not use the device.

### STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 20-25°C until used. Do not freeze or overheat. Protect from light.

### PRODUCT DETERIORATION

This product should not be used if (1) there is evidence of dehydration, (2) the color has changed, (3) the expiration date has passed, or (4) there are other signs of deterioration. The expiration date applies to the product in its intact container when stored as directed.

### SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following recommended guidelines.<sup>4</sup>

### MATERIALS PROVIDED

50 BactiDrop Oxidase – glass ampule enclosed in a protective plastic sheath with a dropper dispenser tip (0.75/ampule)  
1 Ampule crusher (reusable)

### MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Whatman (No.1) filter paper or equivalent, (7) Platinum loop or needle, or wooden applicator stick.

### PROCEDURE

Place dropper in the assembled, reusable ampule crusher provided. Hold the dropper/crusher in an upright position and lightly tap the bottom to dislodge any bubbles that may have formed. Grasp the middle of the dropper/crusher with the thumb and forefinger, and with the tip pointing away, press gently to crush the ampule. Invert dropper and squeeze slightly to dispense in a drop-wise fashion. Discard remaining portion of partially used ampule at end of workday.

### Test Procedure:

1. Dispense 1-2 drops of reagent onto a strip of Whatman (No. 1) filter paper (or an equivalent).
2. Using a platinum loop or needle, or wooden applicator stick, smear the growth of an 18-24-hour old isolate onto the saturated filter paper.
3. Observe for color development within 10-30 seconds.

### INTERPRETATION

Positive Test - Development of a violet to purple color within 10-30 seconds  
  
Negative Test - Delayed reaction or no color development within 10-30 seconds

## QUALITY CONTROL

All lot numbers of BactiDrop Oxidase have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

| CONTROL                                      | INCUBATION                    | RESULTS  |
|----------------------------------------------|-------------------------------|----------|
| <i>Neisseria sicca</i><br>ATCC® 9913         | Aerobic, 30<br>seconds @ 25°C | Positive |
| <i>Pasteurella multocida</i><br>ATCC® 43137  | Aerobic, 30<br>seconds @ 25°C | Positive |
| <i>Pseudomonas aeruginosa</i><br>ATCC® 27853 | Aerobic, 30<br>seconds @ 25°C | Positive |
| <i>Escherichia coli</i><br>ATCC® 25922       | Aerobic, 30<br>seconds @ 25°C | Negative |

## LIMITATIONS<sup>5</sup>

1. Select colonies for testing that have not been grown on selective media or media containing glucose, as selective agents and fermentation may inhibit oxidase enzyme activity resulting in possible false-negative results.
2. Weak oxidase producers, such as *Pasteurella*, may appear negative within the time limits of the test. Results inconsistent with other biochemical reactions or with the organism should be repeated.
3. False negative results may occur if mixed cultures contain the two genera *Pseudomonas* and *Neisseria*. An inhibitory substance is produced by *Pseudomonas* species, which interferes with the production of oxidase by *Neisseria* species.
4. The oxidase reagent quickly auto-oxidizes and loses its sensitivity. Unused reagent should be discarded at the end of each workday.
5. The oxidase test is an oxidizing reaction; oxygen must reach the colonies for optimal results.
6. A platinum inoculating loop or needle, or a wooden applicator stick is recommended for use in oxidase testing.<sup>6</sup> The presence of any trace of iron (nichrome) can catalyze an oxidase reaction, resulting in a false-positive reaction.

## BIBLIOGRAPHY

1. Kovacs, N. 1956. Nature. 178:703.
2. Cowan, S.T. and K.J. Steel. 1966. Manual for the Identification of Medical Bacteria. Cambridge University Press, Cambridge, UK.
3. Steel, K.J. 1962. J. Appl. Bacteriol. 25:445-455.
4. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10<sup>th</sup> ed. ASM Press, Washington, D.C.
5. MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3<sup>rd</sup> ed. Lippincott Williams & Wilkins, Philadelphia, PA.
6. Steel, K.J. 1961. J. Gen. Microbiol. 25:297-306.

## Symbol Legend

|                                                                                   |                                        |
|-----------------------------------------------------------------------------------|----------------------------------------|
|  | Catalog Number                         |
|  | In Vitro Diagnostic Medical Device     |
|  | Batch Code (Lot Number)                |
|  | Consult Instructions for Use (IFU)     |
|  | Temperature Limitation (Storage Temp.) |
|  | Use By (Expiration Date)               |
|  | Contains sufficient for <n> tests      |
|  | Do not use if packaging damaged        |
|  | Authorized European Representative     |
|  | Manufacturer                           |
|  | UK Conformity Assessment               |
|  | European Conformity Assessment         |

Remel Inc.  
12076 Santa Fe Trail Drive  
Lenexa, KS 66215, USA  
[www.thermofisher.com/microbiology](http://www.thermofisher.com/microbiology)  
(800) 255-6730  
International: (913) 888-0939

For technical information contact your local distributor.

BactiDrop is a trademark of Thermo Fisher Scientific and its subsidiaries.

©2022 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue are a trademark of American Type Culture Collection.

CAS (Chemical Abstracts Service Registry No.)

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Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA
©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.