remel

BactiDrop™ Spot Indole EN

REF R21550

INTENDED USE

BactiDrop Spot Indole is a reagent for use in qualitative procedures to detect the ability of a microorganism to split indole from the tryptophan. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

SUMMARY AND EXPLANATION

In 1969, Lowrance, Reich and Traub found *p*-Dimethylaminocinnamaldehyde (DMACA) to be the most sensitive indole reagent, capable of detecting 3 µg of indole per milliliter of medium.¹ Vracko and Sherris, 1963, utilized indole-spot reagent for the presumptive separation of the *Proteus* sp. and *Escherichia coli*.² *p*-Dimethylamino-cinnamaldehyde is also used with anaerobes.

PRINCIPLE

Intracellular enzymes, collectively called "tryptophanase", mediate the production of indole by hydrolytic activity against the amino acid tryptophan. Indole combines with DMACA to form a blue-green compound. The reaction occurs by a condensation process formed by an acid splitting of the protein.³

REAGENTS (CLASSICAL FORMULA)*

<i>p</i> -Dimethylaminocinnamaldehyde	
(CAS 6203-18-5)	10.0 g
Hydrochloric Acid, Concentrate (CAS 7647-01-0)	100.0 ml
Demineralized Water (CAS 7732-18-5)	900.0 ml

*Adjusted as required to meet performance standards.

PRECAUTIONS

WARNING! Causes eye, skin, and respiratory tract irritation. Corrosive to metal.

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers and media after use. Directions should be read and followed carefully. Refer to Safety Data Sheet on company website for detailed information on reagent chemicals.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction, do not use the device

STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 20-25°C until used. Do not freeze or overheat. Protect from light.

PRODUCT DETERIORATION

This product should not be used if (1) there is evidence of dehydration, (2) the color has changed, (3) the expiration date has passed, or (4) there are other signs of deterioration. The expiration date applies to the product in its intact container when stored as directed.

SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{4,5}

MATERIALS PROVIDED

50 BactiDrop Spot Indole – glass ampule enclosed in a protective plastic sheath with a dropper dispenser tip (0.75/ampule)

1 Ampule crusher (reusable)

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Whatman filter paper (No. 1), or equivalent.

PROCEDURE

Place dropper in the assembled, reusable ampule crusher provided. Hold the dropper/crusher in an upright position and lightly tap the bottom to dislodge any bubbles that may have formed. Grasp the middle of the dropper/crusher with the thumb and forefinger, and with the tip pointing away, press gently to crush the ampule. Invert dropper and squeeze slightly to dispense in a dropwise fashion. Discard remaining portion of partially used ampule at end of workday.

Filter Paper Method:

1. Dispense 1-2 drops of BactiDrop Spot Indole onto a piece of No. 1 Whatman filter paper (or equivalent).
2. Smear the growth from an actively growing pure culture onto the saturated filter paper.
3. Observe for the development of a blue color within 1-3 minutes.

Swab Method:

1. Dispense 1-2 drops of BactiDrop Spot Indole onto the tip of a cotton swab.
2. Touch the tip of the swab to the top of the test colony growing on an agar plate.
3. Observe for the development of a blue color within 1-3 minutes.

INTERPRETATION

Positive Test -	Blue color development
Negative Test -	Pink color development

QUALITY CONTROL

All lot numbers of BactiDrop Spot Indole have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL

Bacteroides ovatus
ATCC® 8483

RESULTS

Positive

Escherichia coli
ATCC® 25922

Positive

Prevotella melaninogenica
ATCC® 25845

Negative

Proteus mirabilis
ATCC® 12453

Negative

LIMITATIONS

1. Test only colonies cultured on non-glucose containing media, as glucose inhibits indole production.³
2. Organisms from MacConkey Agar and EMB Agar cannot be tested; indicators in these media may cause a false positive reaction.³
3. Certain strains of *Proteus vulgaris*, *Providencia*, and *Aeromonas* will give a false negative reaction with the Spot Indole test.⁶
4. Media utilized in this test should be checked with known positive and negative control organisms to ensure adequate tryptophan content necessary for the indole reaction.
5. Because adjacent colonies are likely to take up diffused indole, positive tests are valid only if pure cultures are tested.³

BIBLIOGRAPHY

1. Lowrance, B.L., P. Reich, and W.H. Traub. 1969. Appl. Microbiol. 17:923-924.
2. Vracko, R. and J.C. Sherris. 1963. Am. J. Clin. Path. 39:429-432.
3. MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams and Wilkins, Philadelphia, PA.
4. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10th ed. ASM Press, Washington, D.C.
5. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2007. Bailey and Scott's Diagnostic Microbiology. 12th ed. Mosby Elsevier, St. Louis, MO.
6. Balzevic, D.J. and G.M. Ederer. 1975. Principles of Biochemical Tests in Diagnostic Microbiology. John Wiley & Sons, New York, NY.

Symbol Legend

	Catalog Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
	Use By (Expiration Date)
	Contains sufficient for <n> tests
	Do not use if packaging damaged
	Authorized European Representative
	Manufacturer
	UK Conformity Assessment
	European Conformity Assessment

Remel Inc
12076 Santa Fe Trail Drive
Lenexa, KS 66215, USA
www.thermofisher.com/microbiology
(800) 255-6730
International: (913) 888-0939

For technical information contact your local distributor.

BactiDrop is a trademark of Thermo Fisher Scientific and its subsidiaries.

©2022 Thermo Fisher Scientific Inc. All rights reserved

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

CAS (Chemical Abstracts Service Registry No.)

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Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA
©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.