Key Code TSMX7798C
www.oxid.com/ifu

Europe 1 855 236 09 10
CA 1 855 805 8539

US 1 855 236 0910
ROW +31 20 794 7071

Wellcolex Colour Shigella

REF R30858401 50

1. INTENDED USE

Wellcolex™ Colour Shigella is a qualitative latex slide agglutination test for the detection and species identification of *Shigella* isolates present on solid culture media. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

IVD For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Organisms of the genus *Shigella* cause bacillary dysentery, one of the most common diarrhoeal diseases throughout the world. The symptoms, which include acute abdominal pain and bloody diarrhoea, are associated with invasion of the gastric mucosa by the organisms, but further systemic penetration is unusual. The infectious dose may be as low as ten organisms, which assists the rapid dissemination of infection as a result of faecal contamination of the environment. Most cases of infection are relatively mild and require simple supportive therapy, but severe cases can be fatal if appropriate antibiotics are not provided. It is important that the causative organisms are identified accurately and as rapidly as possible to enable the judicious provision of antibiotic therapy, essential to minimise the prevalence of strains with multiple antibiotic resistance and to control outbreaks of disease⁵.

Shigella strains are classified into four species, *Shigella boydii*, *Shigella dysenteriae*, *Shigella flexneri* and *Shigella sonnei*. The classification may be achieved by serological tests on the basis of cell-associated O-antigens, according to which the species may be further subdivided into serological types. Wellcolex Colour Shigella provides a rapid simple serological procedure for species identification of the majority of *Shigella* serotypes encountered in clinical bacteriology.

3. PRINCIPLE OF THE PROCEDURE


In Wellcolex Colour Shigella, a suspension of bacteria is reacted with two test reagents, each consisting of a mixture of red and blue latex particles. Each coloured latex is coated with antibody specific for a different *Shigella* species. In the presence of homologous antigen one of the colours in a mixture will agglutinate, and the identity of the antigen is indicated by the colour of the aggregated particles with a contrasting change in the colour of the background. Each combination is easily distinguished from the others and from a negative result in which the particles remain in smooth purple suspension, and the occasional non-specific result in which all the particles agglutinate into purple aggregates against a cleared background.

4. REAGENTS

KIT CONTENTS	
WELLCOLEX COLOUR SHIGELLA	 50 (ZC51/R30858401)
• Latex Reagent 1	1 dropper bottle (white cap)
• Latex Reagent 2	1 dropper bottle (white cap)
• Red Positive Control	1 dropper bottle (red cap)
• Blue Positive Control	1 dropper bottle (blue cap)
• Disposable Sampling Sticks	2 bundles
• Disposable Reaction Cards	1 pack
• Disposable Sample Dispensers	2 packs
• Disposable Suspension Tubes	1 pack
• Instructions for Use	1
• Reading Guide	1

5. DESCRIPTION, PREPARATION FOR USE AND RECOMMENDED STORAGE CONDITIONS

See also **Warnings and Precautions**

Unless otherwise stated, all components should be stored at 2 to 8°C, under which condition they will retain activity until the expiry date of the kit.

LATEX REAGENT	1
Latex Reagent 1	
One dropper bottle containing a purple suspension of polystyrene latex particles in buffer containing 0.05% Bronidox® preservative. The latex particles are coated with rabbit antibody with the following specificity:	
Red Latex	<i>Shigella sonnei</i> (forms I & II)
Blue Latex	<i>Shigella flexneri</i> (types 1 to 6, X & Y)

LATEX REAGENT	2
Latex Reagent 2	
One dropper bottle containing a purple suspension of polystyrene latex particles in buffer containing 0.05% Bronidox® preservative. The latex particles are coated with rabbit antibody with the following specificity:	
Red Latex	<i>Shigella dysenteriae</i> (types 1 to 12)
Blue Latex	<i>Shigella boydii</i> (types 1 to 15)

RED CONTROL +	Red Positive Control
Contains killed bacterial suspension of organisms representative of <i>Shigella sonnei</i> and <i>Shigella dysenteriae</i> . Containing 0.05% Bronidox® and 0.5% formalin as preservative.	
BLUE CONTROL +	Blue Positive Control
Contains killed bacterial suspension of organisms representative of <i>Shigella flexneri</i> and <i>Shigella boydii</i> . Containing 0.05% Bronidox® and 0.5% formalin as preservative.	

Once opened, reaction cards should be stored at room temperature (18 to 30°C).

6. WARNINGS AND PRECAUTIONS

IVD For *in vitro* diagnostic use only. For professional use only.

Please refer to the Safety Data Sheet, available on the company website, and the product labelling for information on potentially hazardous components.

In the event of malfunction do not use device

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

6.1. HEALTH AND SAFETY INFORMATION

- 6.1.1

In accordance with the principles of Good Laboratory Practice, it is strongly recommended that all clinical specimens and materials should be handled as potentially infectious and used with all necessary precautions. Ideally, the procedure should be performed in a suitable microbiological safety cabinet¹.
- 6.1.2

The Red and Blue Positive Controls contain 0.5% formalin which is classified per applicable European Union and UK Regulations as a sensitizer. The following are the appropriate Hazard (H) and Precautionary (P) statements.

RED CONTROL + BLUE CONTROL + DANGER 	<table><tr><td>H350</td><td>May cause cancer</td></tr><tr><td>EUH208</td><td>Contains formaldehyde. May produce an allergic reaction</td></tr><tr><td>P201</td><td>Obtain special instructions before use</td></tr><tr><td>P281</td><td>Use personal protective equipment as required</td></tr><tr><td>P308+P313</td><td>IF exposed or concerned: Get medical advice/ attention Restricted to professional users</td></tr></table>	H350	May cause cancer	EUH208	Contains formaldehyde. May produce an allergic reaction	P201	Obtain special instructions before use	P281	Use personal protective equipment as required	P308+P313	IF exposed or concerned: Get medical advice/ attention Restricted to professional users
H350	May cause cancer										
EUH208	Contains formaldehyde. May produce an allergic reaction										
P201	Obtain special instructions before use										
P281	Use personal protective equipment as required										
P308+P313	IF exposed or concerned: Get medical advice/ attention Restricted to professional users										

- 6.1.3

Non-disposable apparatus should be sterilised after use, by autoclaving for 15 minutes at 121°C; disposables should be autoclaved or incinerated. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated areas swabbed with an effective disinfectant. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.
- 6.1.4

Do not pipette by mouth. Wear disposable gloves, gowns and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
- 6.1.5

If any of the reagents come into contact with the skin or eyes, wash the area extensively by immediately rinsing with plenty of water.
- 6.1.6

Avoid ingestion of the reagents.

6.2. ANALYTICAL PRECAUTIONS

- 6.2.1

Do not use the reagents beyond the stated expiry date.
- 6.2.2

Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
- 6.2.3

Allow all reagents and samples to come to room temperature (18 to 30°C) before use. Immediately after use ensure bottles are tightly sealed and return to fridge. Latex Reagents which show signs of aggregation when dispensed for the first time may have been frozen and should not be used.
- 6.2.4

Store Latex Reagents and Controls upright at 2 to 8°C. After prolonged storage some aggregation or drying around the top of the bottle may have occurred with the Latex Reagents. Under these circumstances the Latex Reagents should be shaken vigorously for a few seconds until resuspension is complete.
- 6.2.5

Personnel with **colour vision defects** should be able to see agglutination but may experience difficulty in differentiating the colour reactions; if a problem is encountered, reactions should be referred to someone with normal colour vision.
- 6.2.6

It is important to hold the dropper bottles vertically and that the drop forms at the tip of the nozzle. If the nozzle becomes wet, a drop of incorrect volume will form around the end and not at the tip; if this occurs, dry the nozzle before progressing.
- 6.2.7

Use of a flat-bed rotator is essential. For optimum performance ensure that the rotator is level and has been calibrated to run at 150 ± 5 rpm.
- 6.2.8

Do not touch the reaction areas on the cards. It is important to ensure that the Reaction Cards lie completely flat on the rotator otherwise clear agglutination patterns will not be seen.
- 6.2.9

Ensure that the Positive Controls are resuspended by shaking or vortexing vigorously for ten seconds. DO NOT vortex the Latex Reagents.

7. SPECIMEN COLLECTION AND PREPARATION OF CULTURES

For details of specimen collection and preparation of cultures, a standard textbook should be consulted^{2,3,4}. The test may be performed on non-lactose-fermenting colonies growing in primary culture on selective media (for example MacConkey Agar, Hektoen Enteric Agar or Xylose Lysine Deoxycholate Agar), in subculture from enrichment broth on these media or in pure culture (for example Nutrient Agar plates or slopes (slants)).

If it is necessary to confirm the results of the test, the bacterial suspension used for the colour test may be subcultured subsequently for further identification. The same suspension may also be used for testing with Wellcolex Colour Salmonella (ZC50/R30858301).

8. PROCEDURE

MATERIALS PROVIDED

Wellcolex Colour Shigella ZC51/R30858401 contains sufficient materials for 50 tests, see **Kit Contents**.

MATERIALS REQUIRED BUT NOT PROVIDED

1.

0.85% sterile saline.
2.

The test should be performed on a flat bed rotator with an orbital diameter of 1.5 inches (38 mm)/1.25 inches (32 mm) operating at approximately 150 rpm. **The speed and orbital diameter are critical to the performance of the assay.**
- No other materials or equipment are needed.

TEST PROCEDURE

CAUTION: Precautions appropriate to the handling of live cultures should be taken while performing the test.

Colony identification

- Step 1**

Dispense approximately 200 µl of saline into a Suspension Tube. The Disposable Sample Dispenser which is graduated at approximately 200 µl may be used.
- Step 2**

From an overnight culture pick one or two average-sized (1 to 2 mm) suspected *Shigella* colonies from the culture plate using the flat end of a Sampling Stick and carefully emulsify the bacteria in the saline. With small colonies more will need to be picked; the end of the Sampling Stick should be covered. Discard the Sampling Stick for safe disposal.
- Step 3**

Resuspend Latex Reagents 1 and 2 by shaking vigorously for a few seconds. Hold the bottle vertically and dispense one free-falling drop of each Latex Reagent into a separate circle on a flat Reaction Card. Burst any air bubbles with the end of a Sampling Stick.
- Step 4**

Using a Disposable Sample Dispenser held vertically, transfer one free-falling drop (40 µl) of bacterial suspension to each of the two circles. Take care not to dispense air bubbles. Discard the Dispenser for safe disposal.
- Step 5**

Using a Sampling Stick, mix the contents of each circle and spread to cover the whole area of the circle. The same stick may be used for both circles and then discarded for safe disposal.
- Step 6**

Place the card on a suitable flat-bed rotator and run at 150 ± 5 rpm for 2 minutes (see **Analytical Precautions**). Switch off and observe for agglutination without removing the card from the rotator. The card should be viewed from directly above at a normal reading distance (25 to 35 cm from the eyes). Do not use a magnifying lens. The patterns obtained are clear cut and can be recognised easily under normal lighting conditions. If there is any doubt about whether agglutination has occurred the test should be repeated using a 40 µl drop of saline. There should be no visible agglutination. This result should be used as a basis for comparison.
- Step 7**

Discard the used reaction circles for safe disposal. Ensure that the Latex Reagents are returned to the refrigerator (2 to 8°C).

9. RESULTS

READING OF RESULTS

Refer to Wellcolex Colour Shigella Reading Guide.

Positive

A colour change occurs in the reaction due to agglutination of one of the coloured latex suspensions in the mixture which causes a contrasting change in the colour of the background. (Figures 1 and 2). Usually a single colour will agglutinate in one Latex Reagent, but with a mixed *Shigella* culture, both colours may agglutinate in one Latex Reagent (in this event further isolation procedures should be undertaken), or a single colour in both Reagents. These possibilities are easily discriminated. **If the pattern of reaction differs from that shown in Figures 1 and 2, the speed of the rotator should be checked and adjusted if necessary.**

Negative

Neither Latex Reagent agglutinates and the smooth purple appearance remains substantially unchanged throughout the test (Figure 3). Note, however, that faint traces of granularity may be detected in negative patterns, depending on the visual acuity of the operator.

Non-specific

All the particles agglutinate giving rise to purple clumps against a cleared background (Figure 4). It is possible that some purple clumps may develop in the presence of a positive or negative reaction. If there has been a clear colour change in the test or if the background remains uniformly purple, these purple clumps should be ignored.

10. QUALITY CONTROL

Visual Inspection

The Latex Reagents should always be inspected for aggregation as they are dropped onto the Reaction Card and if there is evidence of clumping before addition of the bacterial suspension, that Reagent should not be used.

Control Procedures

Quality control testing should be run with each shipment and new kit lot number received. Each laboratory should follow their state and local requirements.

Positive Control Procedure

Confirm the performance of the Latex Reagents using the Positive Controls provided.

- Step 1**

Into 4 separate Reaction Circles, dispense 2 drops of Latex Reagent 1 and 2 drops of Latex Reagent 2 (one drop per circle).
- Step 2**

Add one drop of **each** Positive Control to 2 of the Reaction Circles (one containing Latex Reagent 1 and one containing Latex Reagent 2).
- Step 3**

Mix the reagents using a separate Disposable Sampling Stick for each Positive Control.
- Step 4**

Rotate the card at **150 ± 5 rpm for 2 minutes** on a mechanical rotator. After this time definitive agglutination should be visible in all 4 Reaction Circles.
- Step 5**

Discard the used Reaction Card for safe disposal.

The colour of the agglutinated latex in Reagent 1 and Reagent 2 should correspond to the colour of the Positive Control (Blue or Red).

Stock cultures of known *Shigella* species may be used in place of the Positive Control.

Negative Control Procedure

Repeat the test procedure using saline instead of the test sample. There should be no significant agglutination.

11. INTERPRETATION OF RESULTS

A positive result (coloured agglutination) indicates the presence and at the same time identifies the species of *Shigella* in the sample, as follows:

Colour of		Reagent	
Clumps	Background	1	2
Red	Blue	<i>Shigella sonnei</i>	<i>Shigella dysenteriae</i>
Blue	Pink	<i>Shigella flexneri</i>	<i>Shigella boydii</i>
None	Purple	Negative	
Purple	Cleared	Not interpretable (further isolation procedures required)	

The reagents will detect *Shigella sonnei* forms I and II, *Shigella flexneri* serotypes 1 to 6 and X and Y variants, *Shigella dysenteriae* serotypes 1 to 12 and *Shigella boydii* serotypes 1 to 15.

A negative result indicates that the organisms under the test do not belong to the serotypes of *Shigella* species covered by the reagents. If require they can be further characterised by alternative procedures^{3,4}. Mixed growths of non-lactose fermenters are frequently present in stool cultures, and if a negative result is obtained it may be necessary to repeat the test on other selected

lactose-negative colonies before discarding the culture as negative.

If a **non-specific reaction** is obtained the result is not interpretable, and alternative procedures for identification should be followed^{3,4}.

12. LIMITATIONS OF THE PROCEDURE

- 12.1.

Wellcolex Colour Shigella will identify *Shigella* isolates to the species level. If full serological identification of types within the species is required, this may be performed by alternative procedures; selection of the appropriate antisera may be guided by the results of the Wellcolex test.
- 12.2.

The genus *Shigella* is closely related to other members of the family *Enterobacteriaceae*, and many of the somatic antigens are shared with organisms from other species, particularly *Escherichia coli*, but including *Hafnia alvei*, *Plesiomonas shigelloides* and other species. These represent common antigens, not cross-reactions, and can be differentiated on the basis of biochemical tests. Many, but not all, ferment lactose.
- 12.3.

Some strains of *Shigella* carry envelope antigens reported to inhibit ‘O’ agglutination of suspensions. These have not been encountered in trials of the Reagents, but if necessary the envelope may be inactivated by heating the suspension at 100°C for two hours². If a biochemically well-characterised isolate still fails to react with the Reagents it may be sent to a reference laboratory for investigation if required.
- 12.4.

The quality of the agglutination pattern obtained is dependent on the speed and orbital diameter of the rotator and the use of a flat Reaction Card.
- 12.5.

Before implementation of WGS (Whole Genome Sequencing), *Shigella flexneri* belonging to ST1753 was misidentified as *Shigella boydii*. This type of *Shigella flexneri* has a mutation which means that the *S. flexneri* O-antigen is not expressed meaning the serology for *S. flexneri* in this instance is falsely negative⁷.

13. EXPECTED RESULTS

Strains belonging to the common serotypes of *S. sonnei*, *S. flexneri*, *S. dysenteriae* and *S. boydii* will give a red or blue agglutination with the corresponding component of Latex Reagent 1 or 2.

14. SPECIFIC PERFORMANCE CHARACTERISTICS

External Evaluation

A multicentre study, involving one Public Health Laboratory in the UK, three hospital laboratories in the USA, two in the Arabian Peninsula and three in South East Asia, was performed to evaluate Wellcolex Colour Shigella on routine cultures⁶.

Each laboratory performed tests on one or more of the following samples from each faecal specimen.

1.

Lactose negative colonies direct from primary selective-differential agar plates (MacConkey, XLD, DCA, SS and Hektoen).
2.

Lactose negative colonies direct from enrichment broth subcultures on selective-differential agar plates.
3.

Pure subcultures of lactose negative colonies on nutrient agar.

A rotator was used throughout. The results are shown in Table 1.

The performance of Wellcolex Colour Shigella was determined by comparison with results of traditional bacteriological analysis of the samples.

In this study the sensitivity of Wellcolex Colour Shigella from primary plate cultures was 98.9% (88/89), 100% (58/58) from enrichment broth subcultures and 100% (299/299) from pure cultures. (See Table 1).

The specificity of Wellcolex Colour Shigella in this study was 99.4% (521/524) from primary plate cultures, 100% (320/320) from enrichment broth subcultures and 95.2% (60/63) from pure cultures.

The predictive value of a positive result from plate cultures was 98.7% (445/451) and of a negative result 99.9% (901/902).

The incidence of *Shigella* in the samples studied was 33.0% (446/1353).

The occurrence of non-specific reactions with Wellcolex Colour Shigella was 3.0% (19/632) for primary plate cultures, 1.8% (7/385) for enrichment broth subcultures and 2.1% (8/385) for pure cultures. These samples have been excluded from Table 1.

Other organisms tested, and found to be negative, with Wellcolex Colour Shigella included *Escherichia coli*, *Proteus spp.*, *Salmonella spp.*, *Campylobacter jejuni*, *Enterobacter spp.* and *Citrobacter freundii*.

Additionally, 208 *Shigella* reference cultures, which included up to five strains from each of the serotypes of each species of *Shigella*, were tested in an independent study⁶. Wellcolex Colour Shigella correctly identified 195 cultures and gave an expected negative result with 13 isolates of *Shigella boydii* serotypes 16 to 18.

NOTE: Consult the scientific literature for further details of serological biochemical test procedures for the identification of *Shigella*.

Table 1 Identification of <i>Shigella</i> from Plate Cultures				
	Wellcolex Colour Shigella Result	Routine Result		Total
		Positive ^a	Negative	
Primary Culture	Positive	88	3 ^{b,c}	91
	Negative	1 ^d	521	522
Subculture From Broth	Positive	58	0	58
	Negative	0	320 ^e	320
Pure Subculture	Positive	299	3 ^{b,f}	302
	Negative	0	60	60
Total		446	907 ^g	1353

^a 33 *Shigella boydii*, 15 *Shigella dysenteriae*, 247 *Shigella flexneri* and 151 *Shigella sonnei*

^b 2 *Plesiomonas shigelloides* (1 *Shigella boydii* and 1 *Shigella sonnei* reaction)

^c *Hafnia alvei* (*Shigella boydii* reaction)

^d *Shigella sonnei*

^e includes 7 cultures where *Shigella* was isolated from primary culture plate only

^f *Escherichia coli* (*Shigella dysenteriae* reaction)

^g includes *Escherichia coli*, *Proteus spp.*, *Salmonella spp.*, *Campylobacter jejuni*, *Enterobacter spp.* and *Citrobacter freundii*

Please note: This study was carried out prior to implementation of WGS (Whole Genome Sequencing). Refer to limitations of the procedure (12.5)

BIBLIOGRAPHY

¹ **Advisory Committee on Dangerous Pathogens.** *The Approved List of biological agents* (July 2021) – issued by the Health and Safety Executive (HSE), UK .

² **Ewing, W.H.** (1986). *Edwards and Ewing’s Identification of Enterobacteriaceae, 4th Ed.*, Elsevier Science Publishing Co., Inc, New York. Pages 166-167.

³ **Isenberg, H.D., Washington, J.A., et al** (1985). *Manual of Clinical Microbiology, 4th Ed.*, Edited by Lennette, E.H., Balows, A., Hausler, W.J. and Shadomy, H.J. American Society for Microbiology, Washington D.C. Pages 73-98.

⁴ **Kelly, M.T., Brenner, D.J., et al** (1985). *Manual of Clinical Microbiology, 4th Ed.*, Edited by Lennette, E.H., Balows, A., Hausler, W.J. and Shadomy, H.J. American Society for Microbiology, Washington D.C. Pages 263-277.

⁵ **Dysentery (Shigellosis).** *Current WHO Guidelines and the WHO Essential Medicine List For Children.* DR Phoebe Williams & Prof James A Berkley. Nov 2016.

⁶ **Data on File.**

⁷ **Chattaway, M. A., Schaefer, U. et al** (2017). *Journal of Clinical Microbiology, Volume 55, Issue 2*, Feb 2017, Pages 616-623.

15. PACKAGING

REF	ZC51/R30858401.....	▽ 50
------------	---------------------	------

16. SYMBOL LEGEND

REF	Catalogue Number
IVD	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
▽ _N	Contains sufficient for <N> tests
	Not for near patient testing
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
UDI	Unique Device Identifier
EC REP	Authorized representative in the European Community
UK CA	UK Conformity Assessed
CE	European Conformity Assessment
	Manufacturer

Bronidox® is the registered trade name of Cognis UK Ltd.

Remel Europe Ltd.
Clipper Boulevard West, Crossways
Dartford, Kent, DA2 6PT, UK
www.thermofisher.com

For technical assistance please contact your local distributor.

Version	Date of modifications introduced
X7798C	June 2023 Updated to meet IVDR requirements

Printed in the United Kingdom

SEZNAM LITERATURY

¹ **Advisory Committee on Dangerous Pathogens.** *The Approved List of biological agents* (July 2021) – issued by the Health and Safety Executive (HSE), UK.

² **Ewing, W.H.** (1986). *Edwards and Ewing’s Identification of Enterobacteriaceae, 4th Ed.*, Elsevier Science Publishing Co., Inc, New York. Pages 166-167.

³ **Isenberg, H.D., Washington, J.A., et al** (1985). *Manual of Clinical Microbiology, 4th Ed.*, Edited by Lennette, E.H., Balows, A., Hausler, W.J. and Shadomy, H.J. American Society for Microbiology, Washington D.C. Pages 73-98.

⁴ **Kelly, M.T., Brenner, D.J., et al** (1985). *Manual of Clinical Microbiology, 4th Ed.*, Edited by Lennette, E.H., Balows, A., Hausler, W.J. and Shadomy, H.J. American Society for Microbiology, Washington D.C. Pages 263-277.

⁵ **Dysentery (Shigellosis).** *Current WHO Guidelines and the WHO Essential Medicine List For Children.* DR Phoebe Williams & Prof James A Berkley. Nov 2016.

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- ¹ **Advisory Committee on Dangerous Pathogens.** *The Approved List of biological agents* (July 2021) – issued by the Health and Safety Executive (HSE), UK.
- ² **Ewing, W.H.** (1986). *Edwards and Ewing’s Identification of Enterobacteriaceae, 4th Ed.*, Elsevier Science Publishing Co., Inc, New York. Pages 166-167.
- ³ **Isenberg, H.D., Washington, J.A., et al** (1985). *Manual of Clinical Microbiology, 4th Ed.*, Edited by Lennette, E.H., Balows, A., Hausler, W.J. and Shadomy, H.J. American Society for Microbiology, Washington D.C. Pages 73-98.
- ⁴ **Kelly, M.T., Brenner, D.J., et al** (1985). *Manual of Clinical Microbiology, 4th Ed.*, Edited by Lennette, E.H., Balows, A., Hausler, W.J. and Shadomy, H.J. American Society for Microbiology, Washington D.C. Pages 263-277.
- ⁵ **Dysentery (Shigellosis).** *Current WHO Guidelines and the WHO Essential Medicine List For Children.* DR Phoebe Williams & Prof James A Berkley. Nov 2016.
- ⁶ **Δεδομένα στο αρχείο.**
- ⁷ **Chattaway, M. A., Schaefer, U. et al** (2017). *Journal of Clinical Microbiology, Volume 55, Issue 2*, Feb 2017, Pages 616-623.

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Remel Europe Ltd.
Clipper Boulevard West, Crossways
Dartford, Kent, DA2 6PT, HB
www.thermofisher.com

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BIBLIOGRAFIA

- ¹ **Advisory Committee on Dangerous Pathogens.** *The Approved List of biological agents* (July 2021) – issued by the Health and Safety Executive (HSE), UK.
- ² **Ewing, W.H.** (1986). *Edwards and Ewing’s Identification of Enterobacteriaceae, 4th Ed.*, Elsevier Science Publishing Co., Inc, New York. Pagg. 166-167.
- ³ **Isenberg, H.D., Washington, J.A., et al** (1985). *Manual of Clinical Microbiology, 4th Ed.*, Edited by Lennette, E.H., Balows, A., Hausler, W.J. and Shadomy, H.J. American Society for Microbiology, Washington D.C. Pagg. 73-98.
- ⁴ **Kelly, M.T., Brenner, D.J., et al** (1985). *Manual of Clinical Microbiology, 4th Ed.*, Edited by Lennette, E.H., Balows, A., Hausler, W.J. and Shadomy, H.J. American Society for Microbiology, Washington D.C. Pagg. 263-277.
- ⁵ **Dysentery (Shigellosis).** *Current WHO Guidelines and the WHO Essential Medicine List For Children*. DR Phoebe Williams & Prof James A Berkley. Nov 2016.
- ⁶ **Dati in archivio.**
- ⁷ **Chattaway, M. A., Schaefer, U. et al** (2017). *Journal of Clinical Microbiology, Volume 55, Issue 2*, Feb 2017, Pagg. 616-623.

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REF	ZC51/R30858401.....	▽	50
------------	---------------------	---	----

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Remel Europe Ltd.
Clipper Boulevard West, Crossways
Dartford, Kent, DA2 6PT, Regno Unito
www.thermofisher.com

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LITERATŪRA

- ¹ **Advisory Committee on Dangerous Pathogens.** *The Approved List of biological agents* (July 2021) – issued by the Health and Safety Executive (HSE), UK.
- ² **Ewing, W.H.** (1986). *Edwards and Ewing’s Identification of Enterobacteriaceae, 4th Ed.*, Elsevier Science Publishing Co., Inc, New York. Pages 166-167.
- ³ **Isenberg, H.D., Washington, J.A., et al** (1985). *Manual of Clinical Microbiology, 4th Ed.*, Edited by Lennette, E.H., Balows, A., Hausler, W.J. and Shadomy, H.J. American Society for Microbiology, Washington D.C. Pages 73-98.
- ⁴ **Kelly, M.T., Brenner, D.J., et al** (1985). *Manual of Clinical Microbiology, 4th Ed.*, Edited by Lennette, E.H., Balows, A., Hausler, W.J. and Shadomy, H.J. American Society for Microbiology, Washington D.C. Pages 263-277.
- ⁵ **Dysentery (Shigellosis).** *Current WHO Guidelines and the WHO Essential Medicine List For Children.* DR Phoebe Williams & Prof James A Berkley. Nov 2016.
- ⁶ **Data on File.**
- ⁷ **Chattaway, M. A., Schaefer, U. et al** (2017). *Journal of Clinical Microbiology, Volume 55, Issue 2*, Feb 2017, Pages 616-623.

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 Remel Europe Ltd
Clipper Boulevard West, Crossways
Dartford, Kent, DA2 6PT, JK
www.thermofisher.com

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Bronidox® este denumirea comercială înregistrată a Cognis UK Ltd.