Bacillus cereus Selective Agar Base (PEMBA)

REF CM0617B, CM0617T, CM0617R

EN

Intended Use

Bacillus cereus Selective Agar Base (PEMBA)(CM0617B/R/T) is a selective medium intended for use in culture media. Bacillus cereus Selective Agar Base (PEMBA)(CM0617B/R/T) is intended to be used for the isolation of *Bacillus cereus* from faecal samples and food samples.

Bacillus cereus Selective Agar Base (PEMBA)(CM0617B/R/T) is used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having infections caused by *Bacillus cereus*.

The devices are for professional use only, are not automated, nor are they companion diagnostics

Summary and Explanation

Bacillus cereus is a Gram-positive bacterium that is widely distributed environmentally¹. It can be found in soil, water, and both fresh and decaying organic matter², with spores that germinate when in contact with organic matter or when inside a host.

B. cereus is a common cause of food-borne gastrointestinal illnesses³. However, *B. cereus* can also cause severe opportunistic extraintestinal infections, particularly in immunocompromised individuals⁴. This includes wound infection, abscess formation, bacteremia, endophthalmitis, pneumonia, meningitis, and septicemia^{2,4,5}. While gastrointestinal illness caused by *B. cereus* is often of low severity, mainly causing a diarrheal or emetic syndrome^{2,6}, extraintestinal infections caused by *B. cereus* can be severe or fatal⁵.

Moreover, *B. cereus* is capable of forming biofilms^{2,4} and shows resistance to most β -lactam antibiotics^{4,5} with resistance to penicillin being an identifying feature of *B. cereus*¹. *B. cereus* is also capable of resisting pasteurization and cleaning processes in hospital settings due to its capacity to sporulate⁶.

Principle of Method

A peptone level of 0.1% and the addition of sodium pyruvate improve egg yolk precipitation and enhance sporulation. Magnesium sulphate provides a source of energy for microbial growth. Disodium hydrogen phosphate and potassium dihydrogen phosphate are added to buffer pH. Bromothymol blue is added as a pH indicator to detect mannitol utilisation. Sodium chloride in the medium maintains osmotic balance. The medium is made selective by the addition of polymyxin B. The primary diagnostic features of the medium are the colonial appearance, precipitation of hydrolysed lecithin and the failure of *B. cereus* to utilise mannitol. Microscopic examination for presence of lipid globules in the vegetative cells is recommended as a rapid and confirmatory test for *B. cereus* and replaces the need for biochemical testing¹.

Typical Formula

	grams per litre
Peptone	1.0
Mannitol	10.0
Sodium chloride	2.0
Magnesium sulphate	0.1
Di-sodium hydrogen phosphate	2.5
Potassium Dihydrogen phosphate	0.25
Sodium pyruvate	10.0
Bromothymol blue	0.12
Agar	15.0

Materials Provided

CM0617B: 500g of Bacillus cereus Selective Agar powder that yields approximately 24.4L after reconstitution.

CM0617R: 2.5kg of Bacillus cereus Selective Agar powder that yields approximately 122L after reconstitution.

CM0617T: 5kg of Bacillus cereus Selective Agar powder that yields approximately 244L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Selective supplements (SR0047C, SR0099E)
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 8°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 29 (Public Health England, 2020).

Procedure

Suspend 20.5g in 475ml of distilled water. Bring to the boil to dissolve completely. Sterilize by autoclaving at 121°C for 15 minutes. Cool to 50°C and aseptically add 25ml of Egg Yolk Emulsion (SR0047) and the contents of 1 vial of Polymyxin B Supplement (SR0099E) reconstituted as directed. Mix well and pour into sterile Petri dishes.

Interpretation

Peacock blue colonies indicates *Bacillus cereus*.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 ± 2 h @ 37° ± 2°C aerobic

Positive Controls	
Inoculum level: 10-100 cfu	
Colony count is ≥ 70% of the control medium count	
<i>Bacillus cereus</i> ATCC®10876	3-5mm peacock blue colonies with halo
<i>Bacillus cereus</i> ATCC®14579	3-5mm peacock blue colonies with halo
Colony count is ≥ 50% of the control medium count	
<i>Staphylococcus aureus</i> * ATCC®25923	0.5-1mm yellow colonies, with or without halo
<i>Enterococcus faecalis</i> * ATCC®19433	0.5-1mm yellow colonies, no halo
Negative Controls	
Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Bacillus coagulans</i> ATCC®7050	No growth
<i>Pseudomonas aeruginosa</i> ATCC®27853	No growth
<i>Salmonella typhimurium</i> ATCC®14028	No growth

Testing performed in accordance with ISO11133:2014	
Incubation Conditions: 21 ± 3 h @ 37° ± 2°C aerobic	
Inoculum level: 10 ³ -10 ⁴ cfu	
<i>Bacillus cereus</i> ATCC®11778	3-5mm peacock blue colonies with halo

Limitations

Strains of *B. cereus* may occasionally show weak or negative reactions to egg yolk. *B. cereus* is indistinguishable from *Bacillus thuringiensis* on this medium. Some other *Bacillus* species, Gram-positive bacteria, and organisms resistant to polymyxin B may also be able to grow on this media.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of Gram-negative bacteria is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device, which must meet the defined acceptance criteria. The precision of *Bacillus cereus* Selective Agar Base (PEMBA)(CM0617B/R/T) was demonstrated by an overall pass rate of 100% obtained for the product over three years of testing (10.SEP.2019 – 25.JUL.2022; 10 batches). This shows that the performance is reproducible.

Bibliography

1. Public Health England. 2018. "Identification of Bacillus Species. UK Standards for Microbiology Investigations. ID 9 Issue 3.1."
2. Bottone, Edward J. 2010. "Bacillus Cereus, a Volatile Human Pathogen." Clinical Microbiology Reviews 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecasse, Daniel H. Bagilet, Laura G. Colombo, and Eduardo R. Gregorini. 2020. "Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound." Revista Argentina de Microbiología 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Hougenhouck-Tulleken, and Mohamed Said. 2020. "Post-Procedural Bacillus Cereus Septic Arthritis in a Patient with Systemic Lupus Erythematosus." African Journal of Laboratory Medicine 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, P C, and J M Kramer. 1983. "Non-Gastrointestinal Bacillus Cereus Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period." Journal of Clinical Pathology 36 (10): 1091–96. <https://doi.org/10.1136/jcp.36.10.1091>
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage, and Marc Germain. 2019. "Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review." Infection Control & Hospital Epidemiology 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests

	Do not use if packaging damaged and consult instructions for use
R _x only	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Manufacturer
EC REP	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2023-12-12

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1. Public Health England. 2018. "Identification of Bacillus Species. UK Standards for Microbiology Investigations. ID 9 Issue 3.1."
2. Bottone, Edward J. 2010. "Bacillus Cereus, a Volatile Human Pathogen." *Clinical Microbiology Reviews* 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo, and Eduardo R. Gregorini. 2020. "Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound." *Revista Argentina de Microbiología* 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Hougenhouck-Tulleken, and Mohamed Said. 2020. "Post-Procedural Bacillus Cereus Septic Arthritis in a Patient with Systemic Lupus Erythematosus." *African Journal of Laboratory Medicine* 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, P C, and J M Kramer. 1983. "Non-Gastrointestinal Bacillus Cereus Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period." *Journal of Clinical Pathology* 36 (10): 1091–96. <https://doi.org/10.1136/jcp.36.10.1091>.
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage, and Marc Germain. 2019. "Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review." *Infection Control & Hospital Epidemiology* 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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Litteratur

1. Public Health England. 2018. "Identification of Bacillus Species. UK Standards for Microbiology Investigations. ID 9 Issue 3.1."
2. Bottone, Edward J. 2010. "Bacillus Cereus, a Volatile Human Pathogen." Clinical Microbiology Reviews 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo og Eduardo R. Gregorini. 2020. "Bacillus Cereus Bakteriæmi hos en patient med et mave-tarmkanal." Revista Argentina de Microbiología 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Hougenhouck-Tulleken og Mohamed Said. 2020. "Post-Operational Bacillus Cereus Septic Arthritis in a Patient with Systemic Lupus Erythematosus." African Journal of Laboratory Medicine 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, PC og JM Kramer. 1983. "Non-Gastrointestinal Bacillus Cereus Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period." Journal of Clinical Pathology 36 (10): 1091-96. <https://doi.org/10.1136/jcp.36.10.1091>
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage og Marc Germain. 2019. "Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review." Infection Control & Hospital Epidemiology 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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Literaturverzeichnis

1. Public Health England. 2018. „Identification of Bacillus Species. UK Standards for Microbiology Investigations. ID 9 Ausgabe 3.1.“
2. Bottone, Edward J. 2010. „Bacillus Cereus, a Volatile Human Pathogen.“ Clinical Microbiology Reviews 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo, and Eduardo R. Gregorini. 2020. „Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound.“ Revista Argentina de Microbiología 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Hougenhouck-Tulleken, and Mohamed Said. 2020. „Post-Procedural Bacillus Cereus Septic Arthritis in a Patient with Systemic Lupus Erythematosus.“ African Journal of Laboratory Medicine 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, P C, and J M Kramer. 1983. „Non-Gastrointestinal Bacillus Cereus Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period.“ Journal of Clinical Pathology 36 (10): 1091–96. <https://doi.org/10.1136/jcp.36.10.1091>.
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage und Marc Germain. 2019. „Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review.“ Infection Control & Hospital Epidemiology 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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1. Public Health England. 2018. "Identification of Bacillus Species. UK Standards for Microbiology Investigations. ID 9 Issue 3.1."
2. Bottone, Edward J. 2010. "Bacillus Cereus, a Volatile Human Pathogen." Clinical Microbiology Reviews 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo, and Eduardo R. Gregorini. 2020. "Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound." Revista Argentina de Microbiología 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Houghenouck-Tulleken, and Mohamed Said. 2020. "Post-Procedural Bacillus Cereus Septic Arthritis in a Patient with Systemic Lupus Erythematosus." African Journal of Laboratory Medicine 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, P C, and J M Kramer. 1983. "Non-Gastrointestinal Bacillus Cereus Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period." Journal of Clinical Pathology 36 (10): 1091–96. <https://doi.org/10.1136/jcp.36.10.1091>
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage, and Marc Germain. 2019. "Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review." Infection Control & Hospital Epidemiology 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Reino Unido

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Kirjandus

1. Public Health England. 2018. "Identification of Bacillus Species. UK Standards for Microbiology Investigations. ID 9 Issue 3.1."
2. Bottone, Edward J. 2010. "Bacillus Cereus, a Volatile Human Pathogen." *Clinical Microbiology Reviews* 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo, and Eduardo R. Gregorini. 2020. "Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound." *Revista Argentina de Microbiología* 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Houghenhouck-Tulleken, and Mohamed Said. 2020. "Post-Procedural Bacillus Cereus Septic Arthritis in a Patient with Systemic Lupus Erythematosus." *African Journal of Laboratory Medicine* 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, P C, and J M Kramer. 1983. "Non-Gastrointestinal Bacillus Cereus Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period." *Journal of Clinical Pathology* 36 (10): 1091–96. <https://doi.org/10.1136/jcp.36.10.1091>.
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage, and Marc Germain. 2019. "Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review." *Infection Control & Hospital Epidemiology* 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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Bibliographie

1. Public Health England. 2018. "Identification of *Bacillus* Species. UK Standards for Microbiology Investigations. ID 9 Issue 3.1."
2. Bottone, Edward J. 2010. "Bacillus Cereus, a Volatile Human Pathogen". *Clinical Microbiology Reviews* 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo et Eduardo R. Gregorini. 2020. "Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound". *Revista Argentina de Microbiología* 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Houghenhouck-Tulleken et Mohamed Said. 2020. "Post-Procedural *Bacillus Cereus* Septic Arthritis in a Patient with Systemic Lupus Erythematosus". *African Journal of Laboratory Medicine* 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, PC, et JM Kramer. 1983. "Non-Gastrointestinal *Bacillus Cereus* Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period." *Journal of Clinical Pathology* 36 (10): 1091–96. <https://doi.org/10.1136/jcp.36.10.1091>
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage et Marc Germain. 2019. "Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review". *Infection Control & Hospital Epidemiology* 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Royaume-Uni

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
Bibliografija

1. Public Health England. 2018. "Identification of Bacillus Species. UK Standards for Microbiology Investigations. ID 9 Issue 3.1."
2. Bottone, Edward J. 2010. "Bacillus Cereus, a Volatile Human Pathogen." Clinical Microbiology Reviews 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo, and Eduardo R. Gregorini. 2020. "Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound." Revista Argentina de Microbiología 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Houghenhouck-Tulleken, and Mohamed Said. 2020. "Post-Procedural Bacillus Cereus Septic Arthritis in a Patient with Systemic Lupus Erythematosus." African Journal of Laboratory Medicine 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, P C, and J M Kramer. 1983. "Non-Gastrointestinal Bacillus Cereus Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period." Journal of Clinical Pathology 36 (10): 1091–96. <https://doi.org/10.1136/jcp.36.10.1091>
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage, and Marc Germain. 2019. "Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review." Infection Control & Hospital Epidemiology 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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Bibliografia

1. Public Health England. 2018. "Identification of *Bacillus* Species. UK Standards for Microbiology Investigations. ID 9 Issue 3.1."
2. Bottone, Edward J. 2010. "Bacillus Cereus, a Volatile Human Pathogen." *Clinical Microbiology Reviews* 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo ed Eduardo R. Gregorini. 2020. "Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound." *Revista Argentina de Microbiologia* 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Hougenhouck-Tulleken e Mohamed Said. 2020. "Post-Procedural Bacillus Cereus Septic Arthritis in a Patient with Systemic Lupus Erythematosus." *African Journal of Laboratory Medicine* 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, P C, and J M Kramer. 1983. "Non-Gastrointestinal Bacillus Cereus Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period." *Journal of Clinical Pathology* 36 (10): 1091-96. <https://doi.org/10.1136/jcp.36.10.1091>
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage e Marc Germain. 2019. "Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review." *Infection Control & Hospital Epidemiology* 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Regno Unito

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Bibliografi

1. Public Health England. 2018. "Identification of Bacillus Species. UK Standards for Microbiology Investigations. ID 9 Issue 3.1."
2. Bottone, Edward J. 2010. "Bacillus Cereus, a Volatile Human Pathogen." Clinical Microbiology Reviews 23 (2): 382-98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo, and Eduardo R. Gregorini. 2020. "Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound." Revista Argentina de Microbiología 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Hougenhouck-Tulleken, and Mohamed Said. 2020. "Post-Procedural Bacillus Cereus Septic Arthritis in a Patient with Systemic Lupus Erythematosus." African Journal of Laboratory Medicine 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, P C, and J M Kramer. 1983. "Non-Gastrointestinal Bacillus Cereus Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period." Journal of Clinical Pathology 36 (10): 1091–96. <https://doi.org/10.1136/jcp.36.10.1091>.
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage, and Marc Germain. 2019. "Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review." Infection Control & Hospital Epidemiology 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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1. Public Health England. 2018. „Identification of Bacillus Species. UK Standards for Microbiology Investigations. ID 9, wydanie 3.1.”
2. Botone, Edward J. 2010. „Bacillus Cereus, a Volatile Human Pathogen.” Clinical Microbiology Reviews 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo i Eduardo R. Gregorini. 2020. „Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound.” Revista Argentina de Microbiología 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Houghenhouck-Tulleken i Mohamed Said. 2020. „Post-Procedural Bacillus Cereus Septic Arthritis in a Patient with Systemic Lupus Erythematosus.” African Journal of Laboratory Medicine 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, P C, i J M Kramer. 1983. „Non-Gastrointestinal Bacillus Cereus Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period.” Journal of Clinical Pathology 36 (10): 1091–96. <https://doi.org/10.1136/jcp.36.10.1091>
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage i Marc Germain. 2019. „Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review.” Infection Control & Hospital Epidemiology 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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Bibliografie

1. Public Health England. 2018. „Identification of Bacillus Species. UK Standards for Microbiology Investigations. ID 9 Issue 3.1.”
2. Bottone, Edward J. 2010. „Bacillus Cereus, a Volatile Human Pathogen.” *Clinical Microbiology Reviews* 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo, and Eduardo R. Gregorini. 2020. „Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound.” *Revista Argentina de Microbiología* 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Hougenhouck-Tulken, and Mohamed Said. 2020. „Post-Procedural Bacillus Cereus Septic Arthritis in a Patient with Systemic Lupus Erythematosus.” *African Journal of Laboratory Medicine* 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, P C, and J M Kramer. 1983. „Non-Gastrointestinal Bacillus Cereus Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period.” *Journal of Clinical Pathology* 36 (10): 1091–96. <https://doi.org/10.1136/jcp.36.10.1091>.
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage, and Marc Germain. 2019. „Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review.” *Infection Control & Hospital Epidemiology* 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Regatul Unit

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Referenser

1. Public Health England. 2018. "Identification of *Bacillus* Species. UK Standards for Microbiology Investigations. ID 9 Issue 3.1."
2. Bottone, Edward J. 2010. "Bacillus Cereus, a Volatile Human Pathogen." *Clinical Microbiology Reviews* 23 (2): 382–98. <https://doi.org/10.1128/CMR.00073-09>.
3. Acosta Pedemonte, Noelia B., Nicolás S. Rocchetti, Juan Villalba, Damián Lerman Tenenbaum, Claudio J. Settecase, Daniel H. Bagilet, Laura G. Colombo och Eduardo R. Gregorini. 2020. "Bacillus Cereus Bacteremia in a Patient with an Abdominal Stab Wound." *Revista Argentina de Microbiología* 52 (2): 115–17. <https://doi.org/10.1016/j.ram.2019.07.003>.
4. Mitton, Barend, Roxanne Rule, Nontombi Mbelle, Wesley van Hougenhouck-Tulleken and Mohamed Said. 2020. "Post-Procedural *Bacillus Cereus* Septic Arthritis in a Patient with Systemic Lupus Erythematosus." *African Journal of Laboratory Medicine* 9 (1). <https://doi.org/10.4102/ajlm.v9i1.1119>.
5. Turnbull, P C, and J M Kramer. 1983. "Non-Gastrointestinal *Bacillus Cereus* Infections: An Analysis of Exotoxin Production by Strains Isolated over a Two-Year Period." *Journal of Clinical Pathology* 36 (10): 1091–96. <https://doi.org/10.1136/jcp.36.10.1091>.
6. Lewin, Antoine, Caroline Quach, Virginie Rigourd, Jean-Charles Picaud, Thérèse Perreault, Pierre Frange, Marc-Christian Domingo, Cindy Lalancette, Gilles Delage, and Marc Germain. 2019. "Bacillus Cereus Infection in Neonates and the Absence of Evidence for the Role of Banked Human Milk: Case Reports and Literature Review." *Infection Control & Hospital Epidemiology* 40 (7): 787–93. <https://doi.org/10.1017/ice.2019.110>.

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Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien