www.thermofisher.com

CCDA Selective Supplement*

REF SR0155E, SR0155H

EN

*This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Oxoid Limited's Campylobacter Blood-Free Selective Agar Base (product codes: CM0739B, CM0739R and CM0739T)

Intended Use

CCDA Selective Supplement (SR0155E, SR0155H) are selective supplements intended to be used with Campylobacter Blood-Free Selective Agar Base (Dehydrated) (CM0739B, CM0739R, or CM0739T) for the isolation of *Campylobacter* species from faecal samples.

CCDA Selective Supplements, when added to CCDA Campylobacter Blood-Free Selective Agar Base, are intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having campylobacteriosis.

The devices are for professional use only, are not automated, and nor are they companion diagnostics.

Summary and Explanation

Campylobacter species are Gram-negative rods in appearance. Two of the most medically relevant *Campylobacter* species isolated from human samples are *Campylobacter coli* and *Campylobacter jejuni*.¹ The symptoms of *Campylobacter* spp. infections are characterised by vomiting, nausea, body aches, abdominal pain and in some cases bloody diarrhoea.² In acute stages, the symptoms of *Campylobacter* infections mimic ulcerative colitis and appendicitis.³ Long term complications of *Campylobacter* infections include arthritis, Guillain Barr Syndrome (GBS) and Irritable Bowel Syndrome (IBS).

The implications of *Campylobacter* infections can vary from person to person.⁴ In most cases, pregnant ladies, patients with AIDS, other blood disorders like thalassemia, patients undergoing chemotherapy and those with a weakened immune system are at high risk.⁵ Moreover, the resistance shown by the *Campylobacter* species to fluoroquinolones like ciprofloxacin complicates treatment.⁶ Campylobacter Blood-Free Selective Agar Base (Dehydrated) (CM0739B/R/T) with CCDA Selective Supplements (SR0155E/H) are therefore a key tool for the clinician to support the isolation of *Campylobacter* species from faecal samples.

Principle of Method

CCDA Selective Supplements (SR0155E and SR0155H) added to DCM contains cefoperazone, a cephalosporin that acts as a broad-spectrum antimicrobial with activity against a wide range of microorganism including Enterobacteriaceae and *Pseudomonas* species. Amphotericin B is added to inhibit fungi, this is required when medium is incubated at 37°C to prevent overgrowth of yeasts.

Typical Formula

	Milligrams per 500ml
Cefoperazone	16.0 mg
Amphotericin B	5.0 mg

Materials Provided

SR0155E: 10 vials - each supplements 500 ml of medium SR0155H: 10 vials - each supplements 2 L of medium

Materials Required but Not Supplied

- Dehydrated culture media base CM0739B/R/T
- Inoculating loops
- Swabs
- Collection containers
- Incubators/ Anaerojar (AG0025A, AB0035A, AB0025A) / CampyGen (CN0025A, CN0035A)
- Quality control organisms

Storage

- Store product in its original packaging between 2°C and 8°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Warnings and Precautions

Signal Word: Danger:

Hazard Statements

H317 - May cause an allergic skin reaction

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled

Precautionary Statements

P261 - Avoid breathing dust

P280 - Wear protective gloves

P285 - In case of inadequate ventilation wear respiratory protection

P342 + P311 - If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician

P302 + P352 - IF ON SKIN: Wash with plenty of soap and water

P333 + P313 - If skin irritation or rash occurs: Get medical advice/attention

P304 + P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing

Each vial is for single use. Do not re-use.

For in vitro diagnostic use only.

For professional use only.

Inspect the product packaging before first use.

Do not use the product if there is any visible damage to the packaging and vial.

Do not use the product beyond the stated expiry date.

Do not use the device if signs of contamination are present.

It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local the recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) Q5.

Procedure


- Allow supplement to equilibrate to room temperature before use.
- Aseptically add 2ml of sterile distilled water to 1 vial and mix gently to dissolve completely. Avoid frothing.
- SR0155E: Add the vial contents to 500ml of sterile Campylobacter Blood-Free Selective Agar Base (CM0739) prepared as directed and cooled to 50°C.
- SR0155H: Add the vial contents to 2l of sterile Campylobacter Blood-Free Selective Agar Base (CM0739) prepared as directed and cooled to 50°C.
- Mix well and pour into sterile Petri dishes.

Bibliography

1. Public Health England. 2014 'Identification of Campylobacter species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
2. Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
3. Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between Campylobacter spp./Arcobacter spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
4. Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with Campylobacter jejuni infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
5. Centers for Disease Control and Prevention. 2019. 'Campylobacter (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
6. Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant Campylobacter jejuni infections. *Frontiers in Public Health*, 9, 820.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit

	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

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Literatura

1. Public Health England. 2014 'Identification of *Campylobacter* species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
2. Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
3. Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
4. Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
5. Centers for Disease Control and Prevention. 2019. '*Campylobacter* (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
6. Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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1. Public Health England. 2014 'Identification of Campylobacter species ID_23i3.1 l. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
2. Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
3. Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between Campylobacter spp./Arcobacter spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
4. Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with Campylobacter jejuni infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
5. Centers for Disease Control and Prevention. 2019. 'Campylobacter (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
6. Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant Campylobacter jejuni infections. *Frontiers in Public Health*, 9, 820.

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1. Public Health England. 2014 'Identification of *Campylobacter* species ID_23i3.1 l. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
2. Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
3. Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
4. Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter jejuni* infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
5. Centers for Disease Control and Prevention. 2019. '*Campylobacter* (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
6. Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter jejuni* infections. *Frontiers in Public Health*, 9, 820.

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Bibliografi

1. Public Health England. 2014 'Identification of *Campylobacter* species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
2. Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
3. Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between *Campylobacter* spp./*Arcobacter* spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
4. Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with *Campylobacter* jejuni infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
5. Centers for Disease Control and Prevention. 2019. '*Campylobacter* (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
6. Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant *Campylobacter* jejuni infections. *Frontiers in Public Health*, 9, 820.

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Literatúra

1. Public Health England.2014 'Identification of Campylobacter species ID_23i3.1 I. <https://www.gov.uk/government/publications/smi-id-23-identification-of-campylobacter-species>.
2. Skirrow M.B. (1977) *Brit. Med. Jour.* 2, 9-11.
3. Akar, M., Aydin, F., Yurci, M. A., Abay, S., Ateş, İ., & Deniz, K. (2018). The possible relationship between Campylobacter spp./Arcobacter spp. and patients with ulcerative colitis. *European Journal of Gastroenterology & Hepatology*, 30(5), 531–538.
4. Ramos, A. P., Leonhard, S. E., Halstead, S. K., Cuba, M. A., Castañeda, C. C., Dioses, J. A., Tipismana, M. A., Abanto, J. T., Llanos, A., & Gourlay, D. (2021). Guillain-Barré syndrome outbreak in Peru 2019 associated with Campylobacter jejuni infection. *Neurology- Neuroimmunology Neuroinflammation*, 8(2).
5. Centers for Disease Control and Prevention. 2019. 'Campylobacter (Campylobacteriosis). <https://www.cdc.gov/campylobacter/technical.html>
6. Rodrigues, J., Cha, W., Mosci, R. E., Mukherjee, S., Newton, D. W., Lephart, P., Salimnia, H., Khalife, W. T., Rudrik, J. T., & Manning, S. D. (2021). Epidemiologic associations vary between tetracycline and fluoroquinolone resistant Campylobacter jejuni infections. *Frontiers in Public Health*, 9, 820.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien