Nutrient Agar

EN

REF CM0003B, CM0003K, CM0003R & CM0003T

Intended Use

Nutrient Agar (CM0003B/K/R/T) devices are a general purpose medium for the cultivation of non-fastidious microorganisms which can be enriched for special purposes. Can be used to maintain microorganisms from clinical samples (e.g. urine and wound).

Nutrient Agar (CM0003B/K/R/T) devices are intended to be used in a diagnostic workflow to support the clinician for the maintenance of bacteria from clinical samples.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Nonfastidious microorganisms are organisms that grow and replicate without special nutritional supplements or conditions. They can synthesize all organic molecules required for their growth. Nonfastidious microorganisms can grow in nutrient sparse environments and can survive in a wider temperature range. This groups of bacteria are involved in causing illnesses and infections in humans and includes *E. coli* and *S. aureus*^{1,2,3,4,5}.

The first step in processing a clinical specimen in order to isolate organisms in the laboratory is the inoculation of culture medium. Cell culture is the most commonly used method for cultivating and isolating intracellular bacteria from clinical samples. There are two possible approaches. First, culture conditions can be standardised to allow for the efficient growth of a diverse range of bacteria. The second strategy consists isolating a previously known bacterium and exposing it to specially adapted culture conditions in order to optimise bacterial growth. Non-selective culture media contain no inhibitors and should allow the growth of the majority of microorganisms found in clinical samples⁶.

Principle of Method

Lab-Lemco beef extract is combined with peptone and sodium chloride to form the basic bouillon described by Loeffler and other early bacteriologists. Yeast extract is added to provide vitamins and minerals to help speed the growth of most organisms. .

Typical Formula

	grams per litre
'Lab-Lemco' powder	1.0
Yeast extract	2.0
Peptone	5.0
Sodium chloride	5.0
Agar	15.0

Materials Provided

CM0003B: 500g of Nutrient Agar
CM0003K: 25kg of Nutrient Agar
CM0003R: 2.5kg of Nutrient Agar
CM0003T: 5kg of Nutrient Agar

500g of Nutrient Agar yields approximately 17.9L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.

- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) Q 5.

Procedure

Suspend 28g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilize by autoclaving at 121°C for 15 minutes. Cool to 50°C. Mix well and pour into sterile containers.

Interpretation

Once the medium is reconstituted:

- White/ straw colonies indicates *Staphylococcus aureus*
- Straw colonies indicates *Pseudomonas aeruginosa*
- Colourless colonies indicates *Streptococcus pyogenes*
- Grey/green colonies indicates *Streptococcus pneumoniae*

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 24 ± 2 hours @ 37 ± 2°C

Positive Controls Inoculum level: 10 – 100 cfu Colony count is ≥ 70% of the control medium count	
<i>Staphylococcus aureus</i> ATCC® 25923™	1-2 mm white/ straw colonies
<i>Pseudomonas aeruginosa</i> ATCC® 27853™	1-3 straw colonies
Enriched with 7% v/v horse blood	
<i>Streptococcus pyogenes</i> ATCC® 19615™	0.25 – 1mm colourless colonies, β haemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6303™	1-2 mm grey/ green colonies, α haemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6305™	0.5-1 mm grey/green colonies, α haemolysis
Enriched with 7% v/v horse blood Incubation Conditions: 37 ± 2°C for 24 ± 2 hours in 5% CO ₂	
<i>Haemophilus influenzae</i> ATCC® 19418™	Pinpoint-0.25mm colourless colonies
Testing performed in accordance with ISO11133:2014 Inoculum level: 50-120 cfu Incubation Conditions: 30 ± 2°C for 24 ± 2 hours	
<i>Yersinia enterocolitica</i> ATCC® 23715™	1-3 mm straw colonies
<i>Yersinia enterocolitica</i> ATCC® 9160™	1-3mm straw colonies

Incubation Conditions: 37 ± 2°C for 24 ± 2 hours	
<i>Escherichia coli</i> ATCC® 25922™	1-3 mm straw colonies
<i>Escherichia coli</i> ATCC® 8739™	1-3 mm straw colonies
<i>Salmonella</i> Typhimurium ATCC® 14028™	1-3 mm straw colonies
<i>Salmonella</i> Enteritidis ATCC® 13076™	1-3 mm straw colonies

Limitations

Nutrient Agar (CM0003B/K/R/T) are non-selective media intended to be used for the maintenance of microorganisms isolated from other culture media products previously inoculated using clinical samples, and there are no expectations of specimens negative for the assay but maintain non-fastidious microorganisms in a condition that can mimic the test condition.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Maintenance of nonfastidious microorganisms is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Nutrient Agar (CM0003B/K/R/T) devices was demonstrated by an overall pass rate of 100% obtained for the product over two years of testing (27.NOV.2020 – 05.SEPT.2022; 10 batches). This shows that the performance is reproducible.

Nutrient Agar (CM0003B/K/R/T) devices are tested in-house as part of the QC process since the products were launched in 1996. For target organisms, when using 10-100 cfu inoculum of *Staphylococcus aureus* ATCC® 25923™ and *Pseudomonas aeruginosa* ATCC® 27853™ or when enriched with 7% v/v horse blood and using 10 – 100 cfu of using 10-100 cfu inoculum of *Streptococcus pyogenes* ATCC® 19615™, *Streptococcus pneumoniae* ATCC® 6303™ or *Streptococcus pneumoniae* ATCC® 6305™ and incubating the device at 37°C ± 2°C for 22 – 26 hours, the user can recover organisms with colony size and morphology as listed in this document. When using 10 – 100 cfu of *Haemophilus influenzae* ATCC® 19418™ on medium enriched with 7% v/v horse blood and incubating the device at 37°C ± 2°C for 22 – 26 hours in a 5% CO₂ atmosphere, the user can recover organisms with colony size and morphology as listed in this document.


When testing in accordance with ISO11133:2014 and using 50-120 cfu inoculum of *Yersinia enterocolitica* ATCC® 23715™ and *Yersinia enterocolitica* ATCC® 9160™ and incubating for 30 ± 2°C for 22 – 26 hours, or *Escherichia coli* ATCC® 25922™ and *Escherichia coli* ATCC® 8739™ at 37 ± 2°C for 22 – 26 hours the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

- Centers for Disease Control and Prevention. 2021. 'Urinary Tract Infection'. <https://www.cdc.gov/antibiotic-use/uti.html>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use

	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
2.0	2024-06-14

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Typisk formel

	<u>Gram pr. liter</u>
"Lab-Lemco"-pulver	1,0
Gærekstrakt	2,0
Pepton	5,0
Natriumchlorid	5,0
Agar	15,0

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© 2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.