

VWR® Pipet Tip compatibility chart

Rack Color ●			VWR Filtered, Volume / Cat. No.										VWR Non-Filtered, Volume / Cat. No.						
			2 µl ●	10 µl ●	10 µl XL ●	20 µl ●	100 µl ●	200 µl ●	300 µl ●	1000 µl ●	1250 µl ●		10 µl ●	10 µl XL ●	200 µl ●	250 µl ●	300 µl ●	1000 µl ●	1250 µl ●
Brand & Model	Volume Range	Type	732-3709	732-3630	732-3631	76322-134	732-3633	76322-150	732-3636	732-3637	732-3638		613-6462	613-6463	613-6465	613-6468	613-6469	613-6470	613-6471
Biohit mLINE®	0.1 - 3 µl	1 ch	X		X								X	X					
Biohit mLINE	0.5 - 10 µl	1 ch	X	X	X								X	X					
Biohit mLINE	2 - 20 µl	1 ch				X									X	X	X		
Biohit mLINE	5 - 100 µl	1 ch				X	X	X							X	X	X		
Biohit mLINE	20 - 200 µl	1 ch				X	X	X	X						X	X	X		
Biohit mLINE	100 - 1000 µl	1 ch								X	X							X	X
Biohit mLINE	5 - 100 µl	8 ch				X	X	X							X	X	X		
Biohit mLINE	30 - 300 µl	8/12 ch				X	X	X	X						X	X	X		
Biohit mLINE	5 - 100 µl	12 ch				X	X	X							X	X	X		
Biohit mLINE	0.5 - 10 µl	12 ch	X	X	X								X	X					
Biohit ProLINE®	0.1 - 3 µl	1 ch	X	X									X	X					
Biohit ProLINE	0.5 - 10 µl	1 ch	X	X	X								X	X				X	X
Biohit ProLINE	100 - 1000 µl	1 ch								X	X								
Biohit ProLINE	0.5 - 10 µl	8 ch	X	X															
Biohit ProLINE® Plus	0.1 - 3 µl	1 ch	X	X	X								X	X					
Biohit ProLINE Plus	0.5 - 10 µl	1 ch											X	X					
Biohit ProLINE Plus	2 - 20 µl	1 ch				X									X	X	X		
Biohit ProLINE Plus	100 - 1000 µl	1 ch								X	X							X	X
Biohit ProLINE Plus	20 - 200 µl	1 ch				X	X	X	X						X	X	X		
Biohit ProLINE Plus	200 µl	1 ch				X	X	X	X						X	X	X		
Biohit ProLINE Plus	0.5 - 10 µl	8 ch	X	X	X								X	X					
Biohit ProLINE Plus	10 - 100 µl	8 ch				X	X	X	X						X	X	X		
Biohit ProLINE Plus	30 - 300 µl	8/12 ch				X	X	X	X						X	X	X		
Biohit ProLINE Plus	10 - 100 µl	12 ch				X	X	X	X						X	X	X		
BRAND Transferpette® S	0.1 - 1 µl	1 ch	X		X								X	X					
BRAND Transferpette S	0.1 - 2.5 µl	1 ch	X		X								X	X					
BRAND Transferpette S	0.5 - 10 µl	1 ch	X	X	X								X	X					
BRAND Transferpette S	2 - 20 µl	1 ch				X	X								X	X			
BRAND Transferpette S	10 - 100 µl	1 ch				X	X								X	X			
BRAND Transferpette S	20 - 200 µl	1 ch				X	X	X	X						X	X	X		

	Rack Color ●		VWR Filtered, Volume / Cat. No.									VWR Non-Filtered, Volume / Cat. No.						
			2 µl ●	10 µl ●	10 µl XL ●	20 µl ●	100 µl ●	200 µl ●	300 µl ●	1000 µl ●	1250 µl ●	10 µl ●	10 µl XL ●	200 µl ●	250 µl ●	300 µl ●	1000 µl ●	1250 µl ●
Brand & Model	Volume Range	Type	732-3709	732-3630	732-3631	76322-134	732-3633	76322-150	732-3636	732-3637	732-3638	613-6462	613-6463	613-6465	613-6468	613-6469	613-6470	613-6471
BRAND Transferpette S	100 - 1000 µl	1 ch								X	X						X	X
BRAND Transferpette S	15 - 300 µl	12 ch				X	X	X	X					X	X	X		
BRAND Transferpette S	20 - 200 µl	12 ch				X	X	X	X					X	X	X		
CappBravo	0.1 - 2 µl	1 ch	X		X	X	X					X	X					
CappBravo	0.5 - 10 µl	1 ch	X	X	X							X	X					
CappBravo	2 - 20 µl	1 ch				X								X	X	X		
CappBravo	5 - 50 µl	1 ch				X	X	X	X					X	X	X		
CappBravo	10 - 100 µl	1 ch				X	X	X	X					X	X	X		
CappBravo	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
CappBravo	0.5 - 10 µl	12 ch	X	X	X							X	X					
CappBravo	5 - 50 µl	12 ch				X	X	X	X					X	X	X		
Eppendorf Reference®	0.1 - 2.5 µl	1 ch	X		X							X	X					
Eppendorf Reference	0.5 - 10 µl	1 ch	X	X	X							X	X					
Eppendorf Reference	2 - 20 µl	1 ch				X								X	X			
Eppendorf Reference	10 - 100 µl	1 ch				X	X							X	X			
Eppendorf Reference	100 - 1000 µl	1 ch								X	X						X	X
Eppendorf Reference® 2	0.1 - 2.5 µl	1 ch	X		X							X	X					
Eppendorf Reference 2	0.5 - 10 µl	1 ch	X	X	X							X	X					
Eppendorf Reference 2	2 - 20 µl	1 ch				X	X	X						X	X	X		
Eppendorf Reference 2	10 - 100 µl	1 ch				X	X							X	X	X		
Eppendorf Reference 2	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
Eppendorf Reference 2	30 - 300 µl	1 ch				X	X	X	X					X	X	X		
Eppendorf Reference 2	100 - 1000 µl	1 ch								X	X						X	X
Eppendorf Research®	0.5 - 10 µl	1 ch	X		X							X	X					
Eppendorf Research	0.1 - 2.5 µl	1 ch	X		X							X	X					
Eppendorf Research	2 - 20 µl	1 ch				X												
Eppendorf Research	10 - 100 µl	1 ch				X	X							X	X			
Eppendorf Research	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
Eppendorf Research	30 - 300 µl	1 ch				X	X	X	X					X	X	X		
Eppendorf Research	100 - 1000 µl	1 ch								X	X						X	X
Eppendorf Research	0.5 - 10 µl	8 ch	X	X								X	X					
Eppendorf Research	30 - 300 µl	8 ch				X	X	X	X					X	X	X		
Eppendorf Research	10 - 100 µl	12 ch				X	X							X	X			
Eppendorf Research plus	0.1 - 2.5 µl	1 ch	X		X							X	X					
Eppendorf Research plus	0.5 - 10 µl	1 ch	X	X	X							X	X					
Eppendorf Research plus	10 - 100 µl	1 ch				X	X							X	X			
Eppendorf Research plus	30 - 300 µl	1 ch				X	X	X	X					X	X	X		
Eppendorf Research plus	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
Eppendorf Research plus	100 - 1000 µl	1 ch								X	X						X	X
Eppendorf Research plus	0.5 - 10 µl	12 ch	X	X	X							X	X					
Eppendorf Research plus	20 - 300 µl	12 ch				X	X	X	X					X	X	X		
Eppendorf Research plus (Blue)	50 - 200 µl	1 ch	X		X							X	X					
Eppendorf Research plus (Light Grey)	2 - 20 µl	1 ch	X	X	X							X	X					
Eppendorf Research plus (Yellow)	2 - 20 µl	1 ch				X								X	X			
Eppendorf Research plus (Yellow)	50 - 200 µl	1 ch				X	X	X	X					X	X	X		
Eppendorf Xplorer®	0.5 - 10 µl	8 ch	X	X	X							X	X					
Eppendorf Xplorer® plus	50 - 1200 µl	8 ch								X	X						X	X
Finnpipette®	100 - 1200 µl	8 ch								X	X						X	X
Finnpipette, Electronic	1 - 10 µl	8 ch	X	X	X							X	X					
Finnpipette, Electronic	5 - 50 µl	8 ch				X	X							X	X	X		

	Rack Color ●		VWR Filtered, Volume / Cat. No.									VWR Non-Filtered, Volume / Cat. No.						
			2 µl ●	10 µl ●	10 µl XL ●	20 µl ●	100 µl ●	200 µl ●	300 µl ●	1000 µl ●	1250 µl ●	10 µl ●	10 µl XL ●	200 µl ●	250 µl ●	300 µl ●	1000 µl ●	1250 µl ●
Brand & Model	Volume Range	Type	732-3709	732-3630	732-3631	76322-134	732-3633	76322-150	732-3636	732-3637	732-3638	613-6462	613-6463	613-6465	613-6468	613-6469	613-6470	613-6471
Finnpipette, Electronic	30 - 300 µl	12 ch				X	X	X	X					X	X	X		
Finnpipette® II	0.5 - 10 µl	1 ch	X	X	X							X	X					
Finnpipette® 4500	0.5 - 10 µl	1 ch	X	X	X							X	X					
Finnpipette 4500	2 - 20 µl	1 ch				X								X	X	X		
Finnpipette 4500	5 - 50 µl	1 ch				X	X	X	X					X	X	X		
Finnpipette 4500 (Micro)	5 - 50 µl	1 ch	X	X	X							X	X					
Finnpipette 4500	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
Finnpipette 4500	100 - 1000 µl	1 ch								X	X						X	X
Finnpipette 4500	200 - 1000 µl	1 ch								X	X						X	X
Finnpipette 4500	5 - 50 µl	12 ch				X	X	X						X	X	X		
Finnpipette® F1	0.2 - 2 µl	1 ch	X	X	X							X	X					
Finnpipette F1	0.5 - 5 µl	1 ch	X		X							X	X					
Finnpipette F1 (Micro)	2 - 20 µl	1 ch	X	X	X							X	X					
Finnpipette F1	10 - 100 µl	1 ch				X	X	X	X					X	X	X		
Finnpipette F1	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
Finnpipette F1	30 - 300 µl	1 ch				X	X	X	X					X	X	X		
Finnpipette F1 (Pink Top)	1 - 10 µl	1 ch	X	X	X							X	X					
Finnpipette F1 (Yellow Top)	1 - 10 µl	1 ch												X	X	X		
Finnpipette® F2	0.2 - 2 µl	1 ch	X	X	X							X	X					
Finnpipette F2	0.5 - 5 µl	1 ch	X	X	X							X	X					
Finnpipette F2	1 - 10 µl	1 ch	X	X	X							X	X					
Finnpipette F2	10 - 100 µl	1 ch				X	X	X	X					X	X	X		
Finnpipette F2	20 - 200 µl	1 ch				X	X	X						X	X	X		
Finnpipette F2	100 - 1000 µl	1 ch								X	X						X	X
Finnpipette F2	10 - 100 µl	8 ch				X	X	X	X					X	X	X		
Finnpipette F2	30 - 300 µl	12 ch				X	X	X	X					X	X	X		
Finnpipette F2 (Yellow Cap)	2 - 20 µl	1 ch				X								X	X	X		
Finnpipette® Focus	5 - 50 µl	1 ch				X	X							X	X	X		
Finnpipette Focus	30 - 300 µl	1 ch				X	X	X	X					X	X	X		
Finnpipette Focus	100 - 1000 µl	1 ch								X	X						X	X
Finnpipette® Novus	30 - 300 µl	8 ch				X	X	X	X					X	X	X		
Gilson PIPETMAN®	20 - 300 µl	8 ch				X	X	X	X					X	X	X		
Gilson PIPETMAN® L	0.2 - 2 µl	1 ch	X	X	X							X	X					
Gilson PIPETMAN L	1 - 10 µl	1 ch	X	X	X							X	X					
Gilson PIPETMAN L	2 - 20µl	1 ch				X								X	X	X		
Gilson PIPETMAN L	10 - 100 µl	1 ch				X	X							X	X	X		
Gilson PIPETMAN L	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
Gilson PIPETMAN L	100 - 1000 µl	1 ch								X	X						X	X
Gilson PIPETMAN L	0.5 - 10 µl	8 ch	X	X	X							X	X					
Gilson PIPETMAN L	2 - 20 µl	8 ch				X								X	X	X		
Gilson PIPETMAN L	20 - 300 µl	8 ch				X	X	X	X					X	X	X		
Gilson PIPETMAN L	0.5 - 10 µl	12 ch	X	X	X							X	X					
Gilson PIPETMAN L	2 - 20 µl	12 ch				X								X	X	X		
Gilson PIPETMAN L	20 - 200 µl	12 ch				X	X	X	X					X	X	X		
Gilson PIPETMAN L	20 - 300 µl	12 ch				X	X	X	X					X	X	X		
Gilson PIPETMAN L P2L	0.2 - 2 µl	1 ch	X	X	X							X	X					
Gilson PIPETMAN L P10L	1 - 10 µl	1 ch	X	X	X							X	X					
Gilson PIPETMAN® Neo	20 - 200 µl	12 ch				X	X	X	X					X	X	X		
Gilson PIPETMAN® P10	1 - 10 µl	1 ch	X	X	X							X	X					
Gilson PIPETMAN P100	10 - 100 µl	1 ch				X	X							X	X			
Gilson PIPETMAN P1000	100 - 1000 µl	1 ch								X	X						X	X
Gilson PIPETMAN P2	0.2 - 2 µl	1 ch	X	X	X							X	X					
Gilson PIPETMAN P20	2 - 20 µl	1 ch				X								X				

	Rack Color ●		VWR Filtered, Volume / Cat. No.									VWR Non-Filtered, Volume / Cat. No.						
			2 µl ●	10 µl ●	10 µl XL ●	20 µl ●	100 µl ●	200 µl ●	300 µl ●	1000 µl ●	1250 µl ●	10 µl ●	10 µl XL ●	200 µl ●	250 µl ●	300 µl ●	1000 µl ●	1250 µl ●
Brand & Model	Volume Range	Type	732-3709	732-3630	732-3631	76322-134	732-3633	76322-150	732-3636	732-3637	732-3638	613-6462	613-6463	613-6465	613-6468	613-6469	613-6470	613-6471
Gilson PIPETMAN P200	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
Gilson PIPETMAN P200	20 - 200 µl	8 ch				X	X	X	X					X	X	X		
Gilson PIPETMAN P200	2 - 20 µl	12 ch				X	X	X	X					X	X	X		
Gilson PIPETMAN® Ultra	1 - 10 µl	1 ch			X								X					
Gilson PIPETMAN Ultra	2 - 20 µl	1 ch				X								X	X	X		
Gilson PIPETMAN Ultra	20 - 100 µl	1 ch				X	X							X	X	X		
Gilson PIPETMAN Ultra	20 - 300 µl	1 ch				X	X	X	X					X	X	X		
Gilson PIPETMAN Ultra	100 - 1000 µl	1 ch								X	X						X	X
Hamilton SoftGrip	0.2 - 2 µl	1 ch	X	X	X							X	X					
Hamilton SoftGrip	1 - 10 µl	1 ch			X								X					
Hamilton SoftGrip	2.5 - 25 µl	1 ch				X	X	X	X					X	X	X		
Hamilton SoftGrip	10 - 100 µl	1 ch				X	X	X	X					X	X	X		
Hamilton SoftGrip	30 - 300 µl	1 ch				X	X	X	X					X	X	X		
Hamilton SoftGrip	100 - 1000 µl	1 ch								X	X						X	X
Hamilton SoftGrip	30 - 300 µl	8 ch				X	X	X	X					X	X	X		
Hirschmann labopette®	0.1 - 3 µl	1 ch	X		X							X	X					
Hirschmann labopette	0.5 - 10 µl	1 ch	X	X	X							X	X					
Hirschmann labopette	2 - 20 µl	1 ch				X								X	X	X		
Hirschmann labopette	10 - 100 µl	1 ch				X	X							X	X	X		
Hirschmann labopette	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
Hirschmann labopette	100 - 1000 µl	1 ch								X	X						X	X
HTL Pipette DP2	0.2 - 2 µl	1 ch	X		X							X	X					
HTL Pipette DP10	0.5 - 10 µl	1 ch	X	X	X							X	X					
HTL Pipette DP20	2 - 20 µl	1 ch				X	X							X	X	X		
HTL Pipette DP100	10 - 100 µl	1 ch				X	X	X						X	X	X		
HTL Pipette DP200	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
HTL Pipette DP1000	100 - 1000 µl	1 ch								X	X						X	X
Labnet Excel™ LE1200-8	100 - 1200 µl	8 ch								X	X						X	X
Nichiryo Nichipet EX	0.5 - 10 µl	1 ch	X	X	X							X	X					
Nichiryo Nichipet EX	2 - 20 µl	1 ch				X	X							X	X	X		
Nichiryo Nichipet EX	10 - 100 µl	1 ch				X	X	X						X	X	X		
Nichiryo Nichipet EX	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
Nichiryo Nichipet EX	100 - 1000 µl	1 ch								X	X						X	X
Nichiryo Nichipet EX II	0.5 - 10 µl	1 ch	X	X	X							X	X					
Nichiryo Nichipet EX II	2 - 20 µl	1 ch				X	X							X	X	X		
Nichiryo Nichipet EX II	10 - 100 µl	1 ch				X	X	X						X	X	X		
Nichiryo Nichipet EX II	20 - 200 µl	1 ch				X	X	X	X					X	X	X		
Nichiryo Nichipet EX II	100 - 1000 µl	1 ch								X	X						X	X
Nichiryo Benchmate II	0.1 - 2 µl	1 ch	X	X	X							X	X					
Nichiryo Benchmate II	2 - 20 µl	1 ch				X								X	X	X		
Nichiryo MultiMate	0.5 - 10 µl	8 ch	X	X	X							X	X					
Nichiryo MultiMate	5 - 50 µl	8 ch				X								X	X	X		
Nichiryo MuLtiMate	0.5 - 10 µl	12 ch	X	X	X							X	X					
Nichiryo MultiMate	50 - 300 µl	12 ch						X	X						X	X		
Sartorius Picus®	0.2 - 2 µl	1 ch	X		X							X	X					
Sartorius Picus	0.2 - 10 µl	1 ch	X	X	X							X	X					
Sartorius Picus	5 - 120 µl	1 ch						X	X						X	X		
Sartorius Picus	10 - 300 µl	1 ch				X	X	X	X					X	X	X		
Sartorius Picus	50 - 1000 µl	1 ch								X	X						X	X
Sartorius Tacta®	0.1 - 3 µl	1 ch	X		X							X	X					
Sartorius Tacta	0.5 - 10 µl	1 ch	X	X	X							X	X					
Sartorius Tacta	2 - 20 µl	1 ch				X	X	X	X					X	X	X		

Rack Color ●			VWR Filtered, Volume / Cat. No.									VWR Non-Filtered, Volume / Cat. No.						
			2 μl ●	10 μl ●	10 μl XL ●	20 μl ●	100 μl ●	200 μl ●	300 μl ●	1000 μl ●	1250 μl ●	10 μl ●	10 μl XL ●	200 μl ●	250 μl ●	300 μl ●	1000 μl ●	1250 μl ●
Brand & Model	Volume Range	Type	732-3709	732-3630	732-3631	76322-134	732-3633	76322-150	732-3636	732-3637	732-3638	613-6462	613-6463	613-6465	613-6468	613-6469	613-6470	613-6471
Sartorius Tacta	10 - 100 μl	1 ch				X	X	X						X	X	X		
Sartorius Tacta	20 - 200 μl	1 ch				X	X	X	X					X	X	X		
Sartorius Tacta	100 - 1000 μl	1 ch								X	X						X	X
Socorex Calibra® 822	1 - 10 μl	1 ch	X	X	X							X	X					
Socorex Calibra 822	10 - 100 μl	1 ch				X	X							X	X	X		
Socorex Calibra 822	20 - 200 μl	1 ch				X	X	X	X					X	X	X		
Socorex Calibra 822	100 - 1000 μl	1 ch								X	X						X	X
VWR® EHP	2 - 20 μl	1 ch				X								X	X	X		
VWR EHP	50 - 250 μl	1 ch				X	X	X	X					X	X	X		
VWR EHP	100 - 1000 μl	1 ch								X	X						X	X
VWR® UHP	0.5 - 10 μl	1 ch	X	X	X							X	X					
VWR UHP	0.1 - 2 μl	1 ch	X	X	X							X	X					
VWR UHP	2 - 20 μl	1 ch				X								X	X	X		
VWR UHP	10 - 100 μl	1 ch				X	X	X	X					X	X	X		
VWR UHP	20 - 200 μl	1 ch				X	X	X	X					X	X	X		
VWR UHP	100 - 1000 μl	1 ch								X	X						X	X