

Mycoplasma Off™

Surface disinfectant for the complete elimination of mycoplasma from cell culture work area and microbiology labs.

Background

Mycoplasma contamination represent a constant threat in cell culture laboratories. Routine disinfection of laboratory work areas, safety cabinets and devices are essential measures against mycoplasma contamination. However, cleaning lab surfaces and equipment with untreated wipes or with ineffective or unspecific disinfectants can favour germ propagation rather preventing it. Mycoplasma Off™ is a ready-to-use solution for cleaning and disinfecting a variety of laboratory surfaces and devices.

Mycoplasma Off™ is highly efficient and rapidly effective against mycoplasma and also against common bacteria and fungi such as *Staphylococcus aureus*, *Enterococcus hirae*, *Pseudomonas aeruginosa*, *Candida albicans*. It is also effective against a wide range of viruses including enveloped viruses incl. coronavirus SARS-CoV-2 and adenovirus.

Protection and Precaution Information

Do not use for alcohol-sensitive surfaces. To avoid damage or discoloration, it is recommended to spot test the surfaces (e.g. coated or sensitive surfaces) to be treated prior to use. Wear disposable gloves and safety glasses when using the product. Ensure adequate ventilation after usage. For research use only.



Instructions for Use

For disinfecting and cleaning of surfaces and laboratory equipment. Spray Mycoplasma Off™ directly onto the surface and dry with paper towels. For decontamination of **electronic laboratory devices** we recommend Mycoplasma Off™ Wipes.

Active ingredients: Propan-1-ol, Ethanol, and Methoxypropoxypropanol

Cat. No.:
15-1000 (1 L Spray)
15-5000 (5 L Refill)
15-1001 (50 Wipes in dispenser box)
15-5001 (50x50 refill sachets)

Germany:
Minerva Biolabs GmbH
Schkopauer Ring 13 · 12681 Berlin · Germany
Telefon +49 30 200 04 37-0
order@minerva-biolabs.com
www.minerva-biolabs.com

USA / Canada:
Minerva Biolabs Inc.
400 Route 518 · Skillman · NJ 08558-2401 · USA
Phone +1 (609) 309-5003
order@minervabiolabs.us
www.minervabiolabs.us