

Plate Cells Quickly and Evenly


Rattler™ Plating Beads

- Fast & Easy-To-Use: Simply pour beads onto the plate, and shake to quickly spread cells evenly over the entire surface. Ideal for processing multiple plates simultaneously.
- No Flaming Required: 4.5 mm glass plating beads are provided sterile in polycarbonate bottles.
- **Reusable:** Just clean and autoclave after use.

Plate Cells Quickly and Evenly Over the Entire Surface



Shake Beads to Spread Cells



Pour off Rattler™ Beads, Autoclave and Reuse



Product	Cat. No.	Size	Specifications
Rattler™ Plating Beads - 230 g/bottle	S1001 S1001-5	1 bottle 5 bottles	Material: Solid, glass 4.5mm beads can be washed, autoclaved, and reused. Packaging: Polycarbonate, autoclavable wide mouth bottle. The bulk format is supplied non-sterile as a 25 kg bag.
Rattler™ Plating Beads - bulk format (non-sterile)	S1001-B	25 kg bag	

www.zymoresearch.com



