

TaqMan® RNase P Fast 384-Well Instrument Verification Plate

Pub. Part no. 4455286 Rev. D Rev. Date October 2011

Part Number	Part	Storage Conditions
4455280	TaqMan® RNase P Fast 384-Well Instrument Verification Plate	Store at -15 to -25°C

Note: For safety and biohazard guidelines, refer to the “Safety” section in the appropriate user guide for your instrument. For every chemical, read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Description

The TaqMan® RNase P Fast 384-Well Instrument Verification Plate is used to verify the performance of Applied Biosystems® Real-Time PCR instruments. Each well contains reaction mix consisting of human genomic DNA template at a defined copy number, 1× TaqMan® Fast Universal PCR Master Mix, and RNase P primers and probe (FAM™/MGB dye-labeled).

For information on plate use, refer to the appropriate user guide for your instrument.

For more information, visit our web site: <http://www.lifetechnologies.com>



For Research Use Only. Not intended for any animal or human therapeutic or diagnostic use.

NOTICE TO PURCHASER: PLEASE REFER TO THE APPROPRIATE INSTRUMENT USER GUIDE, PRODUCT INSERT AND PROTOCOL FOR LIMITED LABEL LICENSE OR DISCLAIMER INFORMATION.

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

TaqMan is a registered trademark of Roche Molecular Systems, Inc.

© 2011 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

For support visit www.appliedbiosystems.com/support
www.lifetechnologies.com

