

TISSUE PROCESSOR ; MODEL NO.: RD-1301

ABOUT THE PRODUCT

After the power failure has the function of memory and program program memory,the machine in the operation of power failure,can continue to run according to the original program for up to six hours.The operation of the machine and dehydration basket by gear and rack transmission system,it will effectively reduce the mechanical failure.

- 20 editable programs can be stored in the system.
- With the protection function to prevent accident stop of cylinder.
- 7inch color touching LCD screen ,Human-machine conversation ,easy operation.
- Fully intelligent design, enabling timely automatic determination and deal from an abnormal event.
- Turn-on:within a year that set arbitrary within years, months, days, hours, minutes, seconds at any time start-up.
- Internal dry heating mechanism and triple protection,Precise temperature control,escape from harmless water vapour.
- Ion clarification technique and decomposition of harmful gases through a pipe.Also with glass cover, environmentally friendly and safe.





DETAIL PRESENTATION



INNOVATION • QUALITY • RELIABLE



- Internal dry heating mechanism and triple protection.
- With automatic reset and automatic positioning function.
- Human-machine conversation interface provides clear instructions for each step.
- High-quality LCD screen above the machine,offers clear display and simple operation.
- Green inner-cycling air purification system to efficiently adsorb, remove poisonous gas.
- 20 editable programs can be stored in the system,the user can choose different system according to the specimen.
- Turn-on: within a year that set arbitrary within years, months, days, hours, minutes, seconds at any time start-up.
- Fully intelligent design, enabling timely automatic determination and recovery from an abnormal event,to ensure the specimen from damage.
- The operation of the machine and dehydration basket by gear and rack transmission system, it will effectively reduce the mechanical failure.
- After the power failure has the function of memory and program memory, the machine in the operation of power failure ,can continue to run according to the original program for up to six hours.

Model No. RD-1301 Specification List		
Photo		
Capacity of single cup	1000ML	1500ML
Number of cups	12 (last 3 are paraffin cups)	12 (last 3 are paraffin cups)
Capacity of basket	60 specimens (Stander cassette)	90 specimens (Stander cassette)
Single cylinder working time	Adjustable within 0-99 hours	Adjustable within 0-99 hours
Temperature range of paraffin cup	Adjustable within RT-99C	Adjustable within RT-99°C
Frequency of agitation	1 times/20 seconds	1 times/20 seconds
Tissue protection station	7th station	7th station
Liquid shake times	1-15 times adjustable	1-15 times adjustable
Voltage and frequency	AC220V±10% 50Hz AC110V±10% 60Hz	AC220V±10% 50Hz AC110V±10% 60Hz
Dimensions	1200*430*490mm	1200*650*620mm