



Liofilchem® Giotto 2



- ☒ Automatic reading and interpretation of Liofilchem® ID and AST systems.
- ☒ Ideal for clinical and industrial microbiology lab.
- ☒ Color digital imaging technology.
- ☒ Electronic archive of test results and images.
- ☒ Automatic updates by built-in internet connection.
- ☒ Remote technical assistance.
- ☒ Compact format.



Convenient, Fast and Easy to Use




LIOFILCHEM® s.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY

Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

Liofilchem® Giotto 2



Analysis in 3 quick steps

1. Select the Liofilchem® ID/AST system.
2. Enter the sample data.
3. Read. The user reviews and confirms, or can adjust any results.

Save, print out and browse results

- ✓ Automatic save of images and results of performed analysis.
- ✓ Print out the test results or send by e-mail.
- ✓ Recall data by test code, patient name, microorganism or other criteria.



One-click Data Protection:
simply click on Backup Button,
the backup takes few seconds

Quick Calibration in one minute,
with the caliper to place into the tray



Upload Images for
Remote Assistance

Software Update Alert



LIOFILCHEM® s.r.l.

Via Sciozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY

Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

Rev.1 / 25.05.2016