



Lab plastic essentials

# Versi-Dry surface protectors Protect laboratory surfaces from damage and contamination

Thermo Scientific<sup>™</sup> Versi-Dry<sup>™</sup> Surface Protectors can help make your lab a safer place by quickly absorbing spills and cushioning breakable labware. Use Versi-Dry protectors as a bench protector, shelf liner, cabinet liner or for wiping up isolated spills. The skid-resistant polyethylene backing protects your benchtop and other lab surfaces from damage and contamination caused by liquid drips and spills. The cloth-like top layer is longlasting and will not fall apart, even when wet.

#### Versi-Dry surface protectors

- Thousands of Thermo Scientific<sup>™</sup> Thirsty-Cells<sup>™</sup> protectors quilted into the surface quickly absorb spills. Available in white regular absorbency (750ml/m2) or blue super absorbency (1,050 ml/m2)
- Cushioned surface minimizes breakage of fragile glassware when dropped or tipped over
- Low-pile sheets lie flat and are skid-resistant helping to stabilize labware and prevent tipping
- Remains fl at and won't curl up or warp with repeated wetting and drying
- Polyethylene backing prevents chemicals from breaking through to protected surfaces beneath

## thermo scientific

### Versi-Dry resist the following chemical spills:

- Acetone
- Ammonium Hydroxide, 30%\*
- Chlorine Bleach
- Chlorobenzene
- Chloroform
- Ethanolamine
- Hydrochloric Acid, 10%\*

- Hydrofluoric Acid, 10%\*
- Methylene Chloride
- Methyl Ethyl Ketone (MEK)
- Methyl Iso Butyl Ketone (MIBK)
- Mineral Spirits
- Nitric Acid, 10%\*
- Phosphoric Acid, 10%\*

- Potassium Hydroxide, 30%\*
- Sulfuric Acid, 10%
- Sodium Hydroxide, 10N\*
- Toluene
- Xylene
- \*Highest concentration tested

#### Ordering information

Description	Size	Color	No. per case	Cat. No
Regular Versi-Dry Surface Protectors-absorbs 750ml/m <sup>2</sup>				
Mat	46 x 51 cm (18 x 20 in)	White	350	62080-00
Dispenser Rolls	51 cm x 91 m (20 in x 300 ft)	White	1	62050-00
	51 cm x 30 m (20 in x 100 ft)	White	4	62070-00
Rolls	51 cm x 91 m (20 in x 300 ft)	White	2	62060-00
	51 cm x 46 m (20 in x 150 ft)	White	2	62065-00
Super Versi-Dry Surface Protectors – absorbs 1,050ml/m <sup>2</sup>				
Small Mat	46 x 51 cm (18 x 20 in)	Blue	350	74018-00
Large Mat	46 x 102 cm (18 x 40 in)	Blue	175	74040-00
Table Mat	51 x 110 cm (20 x 43 in)	Blue	150	74043-00
Fume Hood Mat	46 x 51 cm (18 x 20 in)	Blue	100	74218-00
Floor Mat	46 x 102 cm (18 x 40 in)	Blue	50	74240-00
Benchtop Roll	51 cm x 76 m (20 in x 250 ft)	Blue	1	74000-00
	51 cm x 30 m (20 in x 100 ft)	Blue	4	74050-00
Floor Mat Roll	102 cm x 61 m (40 in x 200 ft)	Blue	1	74200-00

Get a free sample of Super Versi-Dry protectors. Go to **thermofisher.com/versidry** for details.

Learn more at thermofisher.com/breaktheglasshabit

### thermo scientific

For general laboratory use. © 2022 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. COL25839 0122