

# FastDigest restriction enzymes

## Save more with the FastDigest Value Pack


Experience the power and convenience of Thermo Scientific™ FastDigest™ restriction enzymes by testing 13 of our most popular restriction enzymes, along with both colorless and green 10X FastDigest buffers, before adding to your library.

Enjoy the flexibility of having a small stock of restriction enzymes for your workflow along with the benefits of FastDigest enzymes, including:

- One universal buffer for all 176 FastDigest enzymes
- Designed to complete digestion in 5–15 minutes
- Direct loading on gels with FastDigest green buffers

### Thermo Scientific™ FastDigest™ Value Pack

Includes:	Size	Cat. No.
FastDigest BamHI	20 reactions	K1991
FastDigest BglII	20 reactions	
FastDigest EcoRI	20 reactions	
FastDigest Eco32I	20 reactions	
FastDigest HindIII	20 reactions	
FastDigest KpnI	20 reactions	
FastDigest NdeI	20 reactions	
FastDigest NotI	20 reactions	
FastDigest PstI	20 reactions	
FastDigest Sall	20 reactions	
FastDigest SmaI	20 reactions	
FastDigest XbaI	20 reactions	
FastDigest XhoI	20 reactions	
FastDigest Green Buffer	2 x 1 mL	
FastDigest Buffer	2 x 1 mL	



Order a FastDigest Value Pack with 13 of our most popular restriction enzymes, as well as both colorless and green 10X FastDigest buffers.

## Ordering information

Product	Reactions	Cat. No.
FastDigest Alul	100	FD0014
FastDigest Alw21I	100	FD0024
FastDigest Alw26I*	100	FD0034
FastDigest Alw44I	200	FD0044
FastDigest BamHI	800/2,500	FD0054/5
FastDigest BcnI	200	FD0064
FastDigest BglI	200	FD0074
FastDigest BglII	100/200	FD0083/4
FastDigest Bpu1102I	50	FD0094
FastDigest Bsp119I	200	FD0124
FastDigest Bsp120I	200	FD0134
FastDigest Bsu15I	50/100	FD0143/4
FastDigest BsuRI	400	FD0154
FastDigest Cfr10I	50	FD0184
FastDigest Cfr13I	100	FD0194
FastDigest Csp6I	100	FD0214
FastDigest DraI	200	FD0224
FastDigest Eam1104I*	50	FD0234
FastDigest Eam1105I	100	FD0244
FastDigest Ecl136II	100	FD0254
FastDigest EcoO109I	200	FD0264
FastDigest EcoRI	800/2,500	FD0274/5
FastDigest Eco31I*	50/100	FD0293/4
FastDigest Eco32I	200/400	FD0303/4
FastDigest Eco47I	200	FD0314
FastDigest Eco47III	20	FD0324
FastDigest Eco52I	50	FD0334
FastDigest Eco57I*	20	FD0344
FastDigest Eco72I	200	FD0364
FastDigest Eco81I	50	FD0374
FastDigest Eco88I	200	FD0384
FastDigest Eco91I	200	FD0394
FastDigest Eco105I	50	FD0404
FastDigest Eco130I	200	FD0414
FastDigest Eco147I	100	FD0424
FastDigest Scal	100	FD0434
FastDigest EheI	20/50	FD0443/4
FastDigest Esp3I*	20	FD0454
FastDigest GsuI*	20	FD0464
FastDigest Hin1I	100	FD0474
FastDigest Hin6I	200	FD0484
FastDigest HincII	100	FD0494
FastDigest HindIII	800/2,500	FD0504/5
FastDigest HpaII	200	FD0514
FastDigest KpnI	300	FD0524
FastDigest Kpn2I	100	FD0534
FastDigest MspI	400	FD0544
FastDigest MvaI	200	FD0554
FastDigest MluI	100	FD0564
FastDigest NcoI	20/100/300	FD0573/4/5
FastDigest NdeI	100/300/1,000	FD0583/4/5
FastDigest NotI	20/50/250	FD0593/4/6
FastDigest PaeI	20	FD0604
FastDigest PstI	800/2,500	FD0614/5
FastDigest PvuI	20	FD0624
FastDigest PvuII	400	FD0634
FastDigest Sall	200	FD0644
FastDigest SdaI	50	FD0654
FastDigest SmaI	100/200	FD0663/4

Product	Reactions	Cat. No.
FastDigest TaqI	400	FD0674
FastDigest XbaI	300/750	FD0684/5
FastDigest XhoI	400/1,200	FD0694/5
FastDigest Bst1107I	100	FD0704
FastDigest Van91I	100	FD0714
FastDigest BclI	300	FD0724
FastDigest Mph1103I	100	FD0734
FastDigest CpoI	20	FD0744
FastDigest MunI	20/50	FD0753
FastDigest Psp5I	50	FD0764
FastDigest SspI	100	FD0774
FastDigest Pfl23I	40	FD0784
FastDigest HinfI	400	FD0804
FastDigest MboI	50	FD0814
FastDigest MboII*	50	FD0824
FastDigest BspTI	150	FD0834
FastDigest Pfl23II	50	FD0854
FastDigest BseGI*	100	FD0874
FastDigest BseNI*	100	FD0884
FastDigest Bsh1285I	100	FD0894
FastDigest Acc65I	100	FD0904
FastDigest VspI	200	FD0914
FastDigest Bsh1236I	100	FD0924
FastDigest Bsp1407I	50/100	FD0933
FastDigest Psp1406I	20	FD0944
FastDigest Mva1269I*	50	FD0964
FastDigest NheI	50/100	FD0973/4
FastDigest TruI	50	FD0984
FastDigest AatII	50	FD0994
FastDigest BshNI	300	FD1004
FastDigest BpII*	20	FD1014
FastDigest BstXI	100	FD1024
FastDigest KspAI	50	FD1034
FastDigest MnlI*	50	FD1074
FastDigest BseDI	20	FD1084
FastDigest RsaI	100	FD1124
FastDigest SacI	100/200	FD1133
FastDigest Tail	50	FD1144
FastDigest BspLI	20	FD1154
FastDigest BfmI	20	FD1164
FastDigest Bpu10I	20	FD1184
FastDigest SdaI	20	FD1194
FastDigest BseLI	100	FD1204
FastDigest MlsI	50	FD1214
FastDigest NsbI	50	FD1224
FastDigest Adel	50	FD1234
FastDigest SmlI	200	FD1244
FastDigest BcuI	20/50	FD1253
FastDigest BseMI*	20	FD1264
FastDigest MbiI	100	FD1274
FastDigest PagI	50	FD1284
FastDigest TatI	20	FD1294
FastDigest XagI	100	FD1304
FastDigest BpII	20	FD1314
FastDigest PstII	100	FD1334
FastDigest MssI	100	FD1344
FastDigest TasI	100	FD1354
FastDigest TaaI	20	FD1364
FastDigest SchI*	100	FD1374

Product	Reactions	Cat. No.
FastDigest XapI	50/100	FD1383/4
FastDigest Cail	50	FD1394
FastDigest BseMII*	20	FD1404
FastDigest ApaI	300	FD1414
FastDigest Bme1390I	152	FD1424
FastDigest BoxI	100	FD1434
FastDigest BseSI	50	FD1444
FastDigest BseXI*	20	FD1454
FastDigest BshTI	20	FD1464
FastDigest XceI	50	FD1474
FastDigest XmiI	50	FD1484
FastDigest NmuCI	20	FD1514
FastDigest PdiI	50	FD1524
FastDigest PdmI	100	FD1534
FastDigest PstII	150	FD1554
FastDigest XmaJI	20	FD1564
FastDigest Hpy8I	20	FD1574
FastDigest OliI	20	FD1634
FastDigest SatI	20	FD1644
FastDigest TauI	20	FD1654
FastDigest FspAI	20	FD1664
FastDigest DpnI	50/100	FD1703/4
FastDigest BseJI	200	FD1714
FastDigest AasI	25	FD1724
FastDigest HpyF10VI	20	FD1734
FastDigest BveI*	50	FD1744
FastDigest PfoI	20	FD1754
FastDigest FspBI	50	FD1764
FastDigest PteI	50	FD1784
FastDigest SsiI	20	FD1794
FastDigest FaqI*	20	FD1814
FastDigest SfiI	150	FD1824
FastDigest HinIII	100	FD1834
FastDigest HhaI	200	FD1854
FastDigest HpyF3I	50	FD1884
FastDigest SgsI	100	FD1894
FastDigest CseI*	20	FD1904
FastDigest LglI*	20	FD1934
FastDigest AjuI*	20	FD1954
FastDigest Ppu21I	50	FD1974
FastDigest RseI	20	FD2004
FastDigest MreI	20	FD2024
FastDigest BspOI	20	FD2044
FastDigest AanI	20	FD2064
FastDigest Lsp1109I*	100	FD2074
FastDigest MauBI	20	FD2084
FastDigest SfaAI	100	FD2094
FastDigest TscAI	100	FD2104
FastDigest CsiI	50	FD2114
FastDigest BmsI*	20	FD2124
FastDigest PteI	20	FD2134
FastDigest FokI*	100	FD2144
FastDigest RruI	100	FD2154
FastDigest KfiI	20	FD2164
FastDigest SqaAI	50	FD2174
FastDigest BfoI	200	FD2184
FastDigest PacI	25	FD2204
FastDigest PspFI	20	FD2224

\* These products are also Type IIS restriction enzymes.

Find out more at [thermofisher.com/fastdigest](https://thermofisher.com/fastdigest)

**ThermoFisher**  
SCIENTIFIC